



PrimeStore® MTM

Technical Product Specification

Product LH-PS-MTM-12

INTENDED USE

PrimeStore® Molecular Transport Medium (PS-MTM) is intended for inactivation of infectious samples and stabilization and transportation of microbial RNA and DNA.

Product Name: PrimeStore® MTM

Product Number LH-PS-MTM-12

- FDA Cleared Class II Medical Device
- Molecular Transport Medium (MTM)
- 2mL fill in a 12mL tube

Tube Description

- Tube Center 0.616in O.D. x 3.850in Height (15.65mm x 97.79mm)
 - Upper End O.D. 15.85mm
- Conical bottom with skirt
- Polypropylene
- Freestanding

Cap Description

- Red
- Threaded with pre-assembled screw caps

Product Description

The product is a ready-to-use capped and labeled 12mL tube that is filled with 2mL of Molecular Transport Medium (MTM) Solution.

- Minimum quantity 100 count boxes (8x8x8 corrugated box)
- Climate controlled room temperature storage
- Leak proof
- Shelf-life information- Coded for 30 months from Date of Manufacturing
- Case Dimensions are 16.75" x 9" x 42.625"
- Weight – 25.5 pounds per case