

USPTO PATENT FULL-TEXT AND IMAGE DATABASE



(11 of 11)

United States Patent
Fischer , et al.

8,080,645
December 20, 2011

****Please see images for: (Certificate of Correction) ****

Biological specimen collection/transport compositions and methods

Abstract

Disclosed are compositions for collecting, storing, and transporting populations of nucleic acids from biological specimens, and clinical, forensic, or environmental samples. Also disclosed are methods for using these compositions as one-step formulations for killing pathogens, inactivating nucleases, and releasing polynucleotides from other cellular components within the sample, and stabilizing the nucleic acids prior to further processing or assay. In particular embodiments, the invention provides a single, one-step, sample collection/transport/storage formulation containing a known quantity of a non-genomic, nucleic acid carrier molecule that serves as an internal reference control to monitor the fidelity of the collection/transportation medium, and measure the integrity of nucleic acids subsequently isolated and purified from the processed sample.

Inventors: Fischer; Gerald W. (Bethesda, MD), Daum; Luke T. (San Antonio, TX)

Assignee: Longhorn Vaccines & Diagnostics LLC (San Antonio, TX)

Family ID: 43011438

Appl. No.: 12/426,890

Filed: April 20, 2009

Prior Publication Data

Document Identifier

US 20090233309 A1

Publication Date

Sep 17, 2009

Related U.S. Patent Documents

<u>Application Number</u>	<u>Filing Date</u>	<u>Patent Number</u>	<u>Issue Date</u>
12243949	Oct 1, 2008		
PCT/US2008/078499	Oct 1, 2008		
60976728	Oct 1, 2007		

Current U.S. Class: **536/23.1; 549/1**

Current CPC Class: C12Q 1/6806 (20130101); C12Q 1/6848 (20130101); C12Q 1/6876 (20130101); C12N 15/1003 (20130101); C12Q 1/6848 (20130101); C12Q 1/6806 (20130101); C12N 15/1003 (20130101); C12Q 1/6851 (20130101); C12Q 2600/166 (20130101); C12Q 2523/113 (20130101); C12Q 2527/125 (20130101); C12Q 2531/113 (20130101); C12Q 2545/101 (20130101); C12Q 2523/113 (20130101); C12Q 2527/125 (20130101); C12Q 2531/113 (20130101); C12Q 2545/101 (20130101); C12Q 2527/125 (20130101); C12Q 2527/127 (20130101); C12Q 2545/101 (20130101)

Current International Class: C07H 21/02 (20060101); C07D 331/02 (20060101)

References Cited [\[Referenced By\]](#)

U.S. Patent Documents

4116777	September 1978	Takatsy et al.
4315073	February 1982	Brown et al.
4355102	October 1982	Quash
4588680	May 1986	Bucher et al.
4981782	January 1991	Judd et al.
5136019	August 1992	Judd et al.
5187060	February 1993	Cerutti et al.
5243030	September 1993	Judd et al.
5252458	October 1993	Liav et al.
5290686	March 1994	Kendal et al.
5316910	May 1994	Rota et al.
5589174	December 1996	Okuno et al.
5631350	May 1997	Okuno et al.
5663055	September 1997	Turner et al.
5719020	February 1998	Liav et al.
5766841	June 1998	Liav et al.
5945515	August 1999	Chomczynski
5958379	September 1999	Regenold et al.
6015664	January 2000	Henrickson et al.

6168915	January 2001	Scholl et al.
6242582	June 2001	Reece et al.
6280928	August 2001	Scholl et al.
6306582	October 2001	Scholl et al.
6376172	April 2002	Scholl et al.
6406842	June 2002	Scholl et al.
6458577	October 2002	Huang
6495316	December 2002	Scholl et al.
6503745	January 2003	Chand et al.
6573080	June 2003	Scholl et al.
6610474	August 2003	Huang
6627396	September 2003	Swanson et al.
6680308	January 2004	Hassan
6720409	April 2004	Okuno et al.
6734292	May 2004	Omura et al.
6811971	November 2004	Klepp et al.
6875600	April 2005	Scholl et al.
6881835	April 2005	Bai et al.
6893814	May 2005	Swanson et al.
6946291	September 2005	Scholl et al.
7494771	February 2009	Picard et al.
7718402	May 2010	Gayral et al.
2001/0021501	September 2001	Scholl et al.
2001/0034022	October 2001	Scholl et al.
2001/0036628	November 2001	Scholl et al.
2002/0054882	May 2002	Okuno et al.
2002/0055094	May 2002	Reece et al.
2002/0081567	June 2002	Henrickson et al.
2002/0169140	November 2002	Prendergast
2003/0119209	June 2003	Kaylor et al.
2003/0203357	October 2003	Huang
2003/0215796	November 2003	Scholl et al.
2003/0219442	November 2003	Mikayama et al.
2004/0009126	January 2004	Pilkiewicz et al.
2004/0071757	April 2004	Rolf
2004/0082549	April 2004	Jomaa
2004/0086849	May 2004	Shimasaki et al.
2004/0142319	July 2004	Yu et al.
2004/0170965	September 2004	Scholl et al.

2005/0009008	January 2005	Robinson et al.
2005/0042229	February 2005	Yang et al.
2005/0090009	April 2005	Cormier et al.
2005/0170334	August 2005	Mikayama et al.
2005/0187213	August 2005	Lang et al.
2006/0286557	December 2006	Basehore et al.
2007/0202497	August 2007	Renuart et al.
2007/0202511	August 2007	Chen et al.
2009/0098527	April 2009	Fischer et al.
2010/0055672	March 2010	Saghbini

Foreign Patent Documents

0 313 224	Apr 1989	EP
0 621 339	Oct 1994	EP
0 675 199	Oct 1995	EP
0 726 316	Aug 1996	EP
1 081 496	Mar 2001	EP
2150281	Oct 2000	RU
WO 92/03454	Mar 1992	WO
WO 92/16619	Oct 1992	WO
WO 94/09035	Apr 1994	WO
WO 94/17106	Aug 1994	WO
WO 97/05248	Feb 1997	WO
WO 03/026567	Apr 2003	WO
WO 2004/004658	Jan 2004	WO
WO 2004/055205	Jul 2004	WO
WO 2004/072270	Aug 2004	WO
WO 2004/084876	Oct 2004	WO
WO 2005/075642	Aug 2005	WO
WO 2005/085274	Aug 2005	WO
WO 2007/133682	Nov 2007	WO

Other References

Kutyavin et al. 3'-Minor groove binder-DNA probes increase sequence specificity at PCR extension temperatures. *Nucleic Acids Res.* (2000) vol. 28, No. 2, pp. 655-661. cited by examiner .

Morre, et al., "RNA Amplification by Nucleic Acid Sequence-Based Amplification with an Internal Standard Reliable Detection of *Chlamydia trachomatis* in Cervical Scrapings and Urine Samples," *J. of Clinical Microbiol.*, 34(12): 3108-3114 (1996). cited by other .

Rosenstraus, et al., "An Internal Control for Routine Diagnostic PCR: Design, Properties, and Effect on Clinical Performance," *J. of Clinical Microbiol.*, 36(1): 191-197 (1998). cited by other .

Daum, et al., Abstract and Poster--"Development of a Real Time Reverse-Transcription PCR (RRT-PCR) Assay for Detection of Influenza a H1N1 2009 From Clinical Respiratory Specimens," Pediatric Infectious Disease Conference ESPID, Nice, France, May 5-8, 2009. cited by other .

Daum, et al., Abstract --"A Molecular Transport Medium (MTM) for Pathogen Inactivation, Ambient Transport and Preservation of RNA from Clinical Samples," ICAAC, Boston, MA, Sep. 12-15, 2010. cited by other .

U.S. Appl. No. 12/243,949, filed Oct. 1, 2008, Fischer et al. cited by other .

U.S. Appl. No. 12/426,890, filed Apr. 20, 2009, Fischer et al. cited by other .

Chomczynski, P. and Sacchi, N., "Single-Step Method of RNA Isolation by Acid Guanidinium Thiocyanate-Phenol-Chloroform Extraction," *Anal. Biochem.*, 162:156-9 (1987). cited by other .

Tortora, et al., "Tools of Biochemistry 5A--Ways to Isolate and Purify Proteins and Other Macromolecules," *Microbiology--An Introduction*, pp. 152-155, 4.sup.th Ed., The Benjamin/Cummings Publishing Company, Inc., United States (1992). cited by other .

Ramanujam. R. et al., "Room-Temperature-Stable PCR Reagents." Cold Spring Harbor Laboratory Press. *PCR Methods and Appl.*, 3:75-76 (1993). cited by other .

Wolff, C. et al, "Single-Tube Nested PCR With Room-Temperature-Stable Reagents," Cold Spring Harbor Laboratory Press, *PCR Methods and Appl.*, 4:376-79 (1995). cited by other .

Matthews, et al., "Immunofluorescence and Fluorescent Antibody Techniques," *Biochemistry*, pp. 461-463, 2.sup.nd Ed., The Benjamin/Cummings Publishing Company, Inc., United States (1996). cited by other .

Schultz, C.L. et al.. "A Lysis, Storage, and Transportation Buffer for Long-Term, Room-Temperature Preservation of Human Clinical Lymphoid Tissue Samples Yielding High Molecular Weight Genomic DNA Suitable for Molecular Diagnosis," *Am. J. Clin. Pathol.*, 111(6):748-52 (1999). cited by other .

Daum, L.T. et al., "Genetic and Antigenetic Analysis of the First A/New Caledonia/20/99-Like H1N1 Influenza Isolates Reported in the Americas," *Emerg. Infect. Dis.*, 8(4):408-12 (Apr. 2002). cited by other .

De Moreau de Gerbehaye, A.I. et al., "Stable Hepatitis C Virus RNA Detection by RT-PCR During Four Days Storage," *BioMed Central, BMC Infectious Diseases*, 2:22 (2002). cited by other .

Daum, L.T. et al., "A Rapid, Single-Step Multiplex Reverse Transcription-PCR Assay for the Detection of Human H1N1, H3N2 and B Influenza Viruses," *J. of Clinic. Virol.*, 25(3): 345-50 (2002). cited by other .

Spackman, E. et al., "Development of a Real-Time Reverse Transcriptase PCR Assay for Type A Influenza Virus and the Avian H5 and H7 Hemagglutinin Subtypes," *J. of Clinic. Microbiol.*, 40(9): 3256-60 (2002). cited by other .

Blacksell, S.D. et al., "The Effect of Sample Degradation and RNA Stabilization on Classical Swine Fever Virus RT-PCR and ELISA methods," *J. Virol. Methods*, 118(1):33-7 (2004). cited by other .

Canas, L.C., "Clinical Laboratory: Selection, Collection, and Transport of Specimens for Viral Cultures," Department of the Air Force, Air Force Institute of Operational Health (AFIOH), Epidemiological Surveillance Division, SDE O1 44-5001, *Virol. Proc. Man.*, 1-8 (2005). cited by other .

Fouchier, R.A.M. et al., "Characterization of a Novel Influenza a Virus Hemagglutinin Subtype (H16) Obtained From Black-Headed Gulls," *J. of Virol.*, 79(5):2814-22 (Mar. 2005). cited by other .

- Krafft, A.E. et al., "Evaluation of PCR Testing of Ethanol-Fixed Nasal Swab Specimens as an Augmented Surveillance Strategy for Influenza Virus and Adenovirus Identification," *J. of Clinic. Microbiol.*, 43(4):1768-75 (Apr. 2005). cited by other .
- Pheng, O.C. et al., "Temperature Related Storage Evaluation of an RT-PCR Test Kit for the Detection of Dengue Infection in Mosquitoes," (Research Note), *Tropical Biomedicine*, 22(1):73-6 (2005). cited by other .
- "USB Taq PCR Master Mix in qPCR," USB Corporation, Tech Tips, 207 (2005). cited by other .
- Bright, R.A., et al., "Adamantane Resistance Among Influenza A Viruses Isolated Early During the 2005-2006 Influenza Season in the United States," *JAMA*, 295(8):891-4 (Feb. 22, 2006). cited by other .
- "Collecting, Preserving and Shipping Specimens for the Diagnosis of Avian Influenza A (H5N1) Virus Infection: Guide for Field Operations," WHO/CDS/EPR/ARO/2006.1 (2006). cited by other .
- Das, A. et al., "Development of an Internal Positive Control for Rapid Diagnosis of Avian Influenza Virus Infections by Real-Time Reverse Transcription-PCR with Lyophilized Reagents," *J. of Clinic. Microbiol.*, 44(9):3065-73 (Sep. 2006). cited by other .
- Daum, L.T., et al., "Molecular Analysis of Isolates From Influenza B Outbreaks in the U.S. and Nepal, 2005," *Arch. of Virol.*, 151:1863-1874 (2006). cited by other .
- Lin, B. et al., "Broad-Spectrum Respiratory Tract Pathogen Identification Using Resequencing DNA Microarrays," *Genome Res.*, 16(4): 527-35 (2006). cited by other .
- Mohany, J. et al., "Multiplex RT-PCR for Detecting Nineteen Respiratory Viruses," *J. of Clinic. Virol.*, 36: S9 (2006). cited by other .
- Daum, L.T. et al., "Real-Time RT-PCR Assays for Type and Subtype Detection of Influenza A and B Viruses," *Influenza & Other Resp. Viruses*, 1(4): 167-75 (2007). cited by other .
- "PCR-Ready Clear Supreme.TM.," Syntezza Bioscience Ltd., available at <http://www.syntezza.com/egt/PCR-Ready.sub.--Clear.sub.--Supreme.pdf> (2006). cited by other .
- Wang, Z., et al., "Identifying Influenza Viruses with Resequencing Microarrays," *Emerg. Infect. Dis.* 12(4):638-46 (2006). cited by other .
- "TechNotes Newsletter," Applied Biosystems, 14(4):1-37 (2007). cited by other .
- Blow, J.A. et al., "Viral Nucleic Acid Stabilization by RNA Extraction Agent," *J. of Virol. Meth.*, 150:41-4 (2008). cited by other .
- Daum, L.T., et al., "Poster--A Rapid, Simplified Collection-to-Detection System for Typing and Subtyping Influenza Viruses Using Real-Time RT-PCR and Culture," American Society for Microbiology (ASM) Conference on Emerging Technologies of Medical Importance for the Diagnosis of Infectious Diseases and the Detection of Pathogenic Microbes, Beijing, China, (2008). cited by other .
- Daum, L.T., et al., "Abstract--Quantitation of Influenza A Virus From Nasal and Lung Tissue of Cotton Rats Using Real-Time RT-PCR and Culture," 26.sup.th Annual Meeting of the European Society for Pediatric Infectious Diseases, Graz, Austria, (2008). cited by other .
- Daum, L.T., et al., "Abstract--Development and Clinical Evaluation of Rapid Real-Time RT-PCR Assays for Detection of Influenza A and B Viruses," 26.sup.th Annual Meeting of the European Society for Pediatric Infectious Diseases, Graz, Austria, (2008). cited by other .
- Daum, L.T., et al., "Poster--A Novel Specimen Collection Solution for Molecular Diagnostic Applications," The Pediatric Academic Societies (PAS) Annual Meeting, Honolulu, HI (2008). cited by other .
- Daum, L.T., et al., "Poster --Real-Time RT-PCR Detection of Influenza A Virus in Asymptomatic Culture-Negative Cotton Rats," The Pediatric Academic Societies (PAS) Annual Meeting,

Honolulu, HI (2008). cited by other .

Daum, L.T., et al., "Poster--Detection and Molecular Characterization of Clinical Influenza A and B Viruses from Original Nasal Wash Specimens Preserved in PrimeStore.TM.." The 3rd European Conference on Influenza, Vilamoura, Portugal (2008). cited by other .

Daum, L.T., et al., "Real-Time RT-PCR Detection of Influenza Virus Within Symptomatic and Asymptomatic Family Members," The 48th Annual IDSA/ICAAC, Washington D.C. (2008). cited by other .

European Patent Office, "PCT International Search Report, PCT Written Opinion of the International Searching Authority--Application No. PCT/US2007/078025," Oct. 28, 2008, 9 pages. cited by other .

"Abstracts--27th Annual Meeting of the European Society for Paediatric Infectious Disease, Brussels, Belgium, Jun. 9-13, 2009," The Ped. Infect. Dis. J., 28(6):e1, e75, e202, e229 (Jun. 2009). cited by other .

Borns, M. et al, "Most Accurate PCR Enzyme Improved With Hot Start Feature," Biocompare, available at <http://www.biocompare.com/technicalarticle/212/Most-Accurate-PCR-Enzyme-Improved-With-Hot-Start-Feature-from-Stratagene.html> (last visited Aug, 24, 2009). cited by other .

"KOD Hot Start DNA Polymerase," Novagen, available at <http://www.emdbiosciences.com/Products/ProductDisplay.asp?catno=71086> (last visited Aug. 24, 2009). cited by other .

"R.A.P.I.D.RTM. System," Idaho Technology Inc., available at <http://www.idahotech.com/RAPID/Rapid-Water.html> (last visited Aug. 24, 2009). cited by other .

"AgPath-iD.TM. One-Step RT-PCR Kit," Applied Biosystems, available at <http://www.abion.com/techlib/prot/bp.sub.--1005.pdf> (last visited Aug. 24, 2009). cited by other .

"PCR Optimization: Reaction Conditions and Components," Applied Biosystems, Part No. 4371091, Revision C, pp. 1-6 available at

<http://www3.appliedbiosystems.com/cms/groups/mcb.sub.--marketing/document-s/generaldocuments/cms.sub.--042520.pdf> (last visited Aug. 24, 2009). cited by other .

"Single Tube PCR Kit Manual," Takara Bio Inc., Cat. #RR021, V.02.09, pp. 1-6 available at <http://www.takara-bio.us/files/manuals/TAK.sub.--RR021.sub.--TM.pdf> (last visited Aug. 24, 2009). cited by other .

"Luminex Confirms Effectiveness of xTAG.RTM. Respiratory Viral Panel for Swine Flu Surveillance," Medical News Today, available at

<http://www.medicalnewstoday.com/printerfriendlynews.php?newsid=148498> (May 1, 2009). cited by other .

"Luminex Receives FDA Clearance for an Update to the xTAG.RTM. Respiratory Viral Panel Package Insert Package Insert to Include Data from Two New Publications on 2009 Influenza A/H1N1," available at <http://phx.corporate-ir.net/phoenix.zhtml?c=79403&p=irol-newsArticle&ID=1-307416&highlight=> (Jul. 14, 2009). cited by other .

European Patent Office, "PCT International Search Report, PCT Written Opinion of the International Searching Authority--Application No. PCT/US2008/078499," Aug. 4, 2009, 13 pages. cited by other.

Primary Examiner: Horlick; Kenneth R.

Assistant Examiner: Thomas; David

Attorney, Agent or Firm: Remenick PLLC

Parent Case Text

CROSS-REFERENCE TO RELATED APPLICATIONS

This application is a continuation-in-part of U.S. patent application Ser. No. 12/243,949, and PCT International Patent Application No. PCT/US2008/078499, both filed Oct. 1, 2008, now pending each of which claims the benefit of U.S. Provisional Patent Application No. 60/976,728, filed Oct. 1, 2007, now expired, the entire contents of each of which is specifically incorporated herein by express reference thereto.

Claims

What is claimed is:

1. An isolated, single-stranded ribonucleic acid molecule of at least 40 nucleotides in length of a sequence prepared by the in vitro transcription of a polynucleotide of the sequence of SEQ ID NO:1 comprising a nucleic acid segment that comprises: (a) a first sequence domain that specifically binds to a labeled probe of from about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length; and (d) a nucleotide sequence that is not present within a mammalian genome, or a genome of a bacterium, fungus, protozoan, virus that is pathogenic to a mammal; wherein the isolated ribonucleic acid molecule remains at least substantially non-degraded when placed in a solution comprising: (i) one or more chaotropes; (ii) one or more detergents; (iii) one or more reducing agents; (iv) one or more chelators; and (v) one or more buffers, collectively present in an amount sufficient to prevent the substantial degradation of a population of polynucleotides comprising the single-stranded ribonucleic acid molecule, at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days, and wherein the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion.
2. The isolated, single-stranded ribonucleic acid molecule of claim 1, of about 80 to about 100 nucleotides in length.
3. The isolated, single-stranded ribonucleic acid molecule of claim 1, of about 30 to about 100 nucleotides in length.
4. The isolated, single-stranded ribonucleic acid molecule of claim 1, wherein the first sequence domain shares at least a first common nucleotide sequence region of about 5 to 10 or more nucleotides with at least one of the second or third sequence domain.
5. The isolated, single-stranded ribonucleic acid molecule of claim 1, wherein the molecule remains at least substantially non-degraded under stringent hybridization conditions comprising: (i) pre-washing in a solution containing about 5.times.SSC, 0.5% SDS, and 1.0 mM EDTA (pH 8.0); (ii) hybridizing at a temperature of

about 60.degree. C. to about 70.degree. C. in 5.times.SSC overnight; and (iii) subsequently washing at about 65.degree. C. to 70.degree. C. for 20 min. with each of 2.times., 0.5.times. and 0.2.times.SSC containing 0.1% SDS), or equivalent hybridization conditions thereto.

6. The isolated, single-stranded ribonucleic acid molecule of claim 1, prepared by the in vitro transcription of a polynucleotide comprising at least 50 contiguous nucleotides from the sequence of SEQ ID NO 1.

7. An isolated, single-stranded ribonucleic acid molecule of at least 40 nucleotides in length of a sequence of SEQ ID NO:2 comprising a nucleic acid segment that comprises: (a) a first sequence domain that specifically binds to a labeled probe of from about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length; and (d) a nucleotide sequence that is not present within a mammalian genome, or a genome of a bacterium, fungus, protozoan, virus that is pathogenic to a mammal; wherein the isolated ribonucleic acid molecule remains at least substantially non-degraded when placed in a solution comprising: (i) one or more chaotropes; (ii) one or more detergents; (iii) one or more reducing agents; (iv) one or more chelators; and (v) one or more buffers, collectively present in an amount sufficient to prevent the substantial degradation of a population of polynucleotides comprising the single-stranded ribonucleic acid molecule, at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days, and wherein the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion.

8. The isolated, single-stranded ribonucleic acid molecule of claim 7, comprising at least 25 contiguous nucleotides from the sequence of SEQ ID NO 2.

9. The isolated, single-stranded ribonucleic acid molecule of claim 7, comprising at least 35 contiguous nucleotides from the sequence of SEQ ID NO 2.

10. The isolated, single-stranded ribonucleic acid molecule of claim 1, wherein the labeled probe comprises at least a first minor groove binder.

11. The isolated, single-stranded ribonucleic acid molecule of claim 1, wherein the labeled probe comprises a radioactive, luminescent, chemiluminescent, fluorescent, enzymatic, magnetic, or spin-resonance label.

12. The isolated, single-stranded ribonucleic acid molecule of claim 1, wherein the mammal is bovine, ovine, porcine, lupine, canine, equine, feline, ursine, murine, leonine, leporine, hircine, primate, or human.

13. An aqueous composition comprising: (a) the isolated single-stranded ribonucleic acid molecule of claim 1; and (b) (i) one or more chaotropes; (ii) one or more detergents; (iii) one or more reducing agents; (iv) one or more chelators; and (v) one or more buffers, collectively present in an amount sufficient to denature protein, inactivate nucleases, kill pathogens, and not degrade nucleic acids in a sample suspected of containing nucleic acids, when the sample is contacted with the composition.

14. The composition of claim 13, wherein (i) the one or more chaotropes are present in an amount from about 0.5 M to about 6 M; (ii) the one or more detergents are present in an amount from about 0.1% to about 1%

(wt./vol.); (iii) the one or more reducing agents are present in an amount from about 0.05 M to about 0.3 M; (iv) the one or more chelators are present in an amount from about 0.01 mM to about 1 mM; and (v) the one or more buffers are present in an amount from about 1 mM to about 1 M.

15. The composition of claim 14, wherein (i) the one or more chaotropes comprise guanidine thiocyanate, guanidine isocyanate, guanidine hydrochloride, or any combination thereof; (ii) the one or more detergents comprise sodium dodecyl sulfate, lithium dodecyl sulfate, sodium taurodeoxycholate, sodium taurocholate, sodium glycocholate, sodium deoxycholate, sodium cholate, sodium alkylbenzene sulfonate, N-lauroyl sarcosine, or any combination thereof; (iii) the one or more reducing agents comprise 2-mercaptoethanol, tris(2-carboxyethyl)phosphine, dithiothreitol, dimethylsulfoxide, or any combination thereof; (iv) the one or more chelators comprise ethylene glycol tetra acetic acid, hydroxyethylethylenediaminetriacetic acid, diethylene triamine penta acetic acid, N,N-bis(carboxymethyl)glycine, ethylenediaminetetraacetic, citrate anhydrous, sodium citrate, calcium citrate, ammonium citrate, ammonium bicitrate, citric acid, diammonium citrate, ferric ammonium citrate, lithium citrate, or any combination thereof; or (v) the one or more buffers comprise tris(hydroxymethyl)aminomethane, citrate, 2-(N-morpholino)ethanesulfonic acid, N,N-Bis(2-hydroxyethyl)-2-aminoethanesulfonic acid, 1,3-bis(tris(hydroxymethyl)methyl amino)propane, 4-(2-hydroxyethyl)-1-piperazine ethanesulfonic acid, 3-(N-morpholino) propanesulfonic acid, bicarbonate, phosphate, or any combination thereof.

16. The composition of claim 13, further comprising: (a) one or more surfactants selected from the group consisting of a silicone polymer, a polysorbate, and any combination thereof; (b) one or more short-chain alkanols selected from the group consisting of methanol, ethanol, propanol, butanol, pentanol, hexanol, and any combination thereof; or (c) one or more defoaming agents selected from the group consisting of a silicone polymer, a polysorbate, and any combination thereof.

17. The composition of claim 13, which (a) is buffered to a pH of about 6.0 to 7.0; (b) is substantially free of RNase and DNase activity; or (c) further comprises a population of isolated polynucleotides of bacterial, viral or fungal origin that comprises RNA, DNA, or any combination thereof.

18. The composition of claim 13, comprising: (a) (i) about 4 M guanidine thiocyanate; (ii) about 30 mM sodium citrate; (iii) about 0.25% (wt./vol.) sodium dodecyl sulfate; (iv) about 0.25% (wt./vol.) N-lauroyl sarcosine, sodium salt; (v) about 0.1 M 2-mercaptoethanol; and (vi) about 0.1% silicone polymer (wt./vol.); (b) (i) about 3 M guanidine thiocyanate; (ii) about 1 mM TCEP; (iii) about 10 mM sodium citrate; (iv) about 0.5% N-lauroyl sarcosine; (v) about 0.0002% silicone polymer; (vi) about 100 mM 2-amino-2-hydroxymethyl-propane-1,3-diol (TRIS); and (vii) about 0.1 mM EDTA; (c) (i) about 1 M to about 4 M guanidine thiocyanate; (ii) about 0.5 mM to 10 mM TCEP; (iii) about 1 mM to 100 mM sodium citrate; (iv) about 0.1% to about 1% SDS or NLS; (v) about 0.001% to about 0.0001% of a silicone polymer, (vi) about 10 mM to about 500 mM TRIS, (vii) about 0.1 mM to about 1 mM APCA, EDTA, EGTA, HEDTA, DTPA, NTA, or citrate; and (viii) about 10% to about 25% ethanol (vol./vol.); or (d) (i) about 3 M guanidine thiocyanate; (ii) 1 mM TCEP; about 10 mM sodium citrate; (iii) about 0.5% N-lauroyl sarcosine, sodium salt; (iv) about 0.0002% of a silicone polymer; (v) about 100 mM TRIS; (vi) about 0.1 mM EDTA; and (vii) about 10% to about 25% ethanol (vol./vol.).

19. A method for increasing the efficiency of obtaining a purified population of polynucleotides from a biological sample suspected of containing such polynucleotides, comprising contacting the sample with a composition comprising the isolated single-stranded ribonucleic acid molecule of claim 1, in an amount and for a time sufficient to increase the efficiency of obtaining the purified population of polynucleotides from the

biological sample.

20. The method of claim 19, wherein the biological sample is of clinical, veterinary, epidemiological, environmental, forensic, or pathological origin; or wherein it contains, or is suspected of containing one or more polynucleotides of viral, bacterial, fungal, or mammalian origin.

21. The method of claim 19, wherein the composition further comprises (a) one or more chaotropes; (b) one or more detergents; (c) one or more reducing agents; (d) one or more chelators; and (e) one or more buffers, collectively present in an amount sufficient to denature proteins, inactivate nucleases, kill pathogens, and not degrade nucleic acids in the sample.

22. The method of claim 19, wherein the integrity of a population of polynucleotides in the biological sample or the fidelity of at least a first sequence of one of the polynucleotides is (a) at least substantially maintained when the composition comprising the sample is stored at a temperature of about minus 20.degree. C. to about 40.degree. C. for a period of about 7 to about 14 days; or (b) at least substantially maintained when the composition comprising the sample is stored at a temperature of about minus 20.degree. C. to about 40.degree. C. for a period of about 14 to about 30 days.

23. The method of claim 19, wherein (a) one or more cells in the sample is lysed to release at least a first population of polynucleotides; (b) the sample further comprises one or more pathogens, at least a portion of which is killed or at least inactivated by the composition; or (c) the composition comprising the sample is stored substantially at ambient temperature from the time of collection to the time of analyzing a population of polynucleotides therein.

24. The method of claim 19, wherein (a) the composition comprising the sample is stored at a temperature of about 10.degree. C. to about 40.degree. C. substantially from the time of collection to the time of isolating, purifying, or characterizing a population of polynucleotides therein; or (b) less than about 5% of the population of polynucleotides contained in the sample is degraded after the composition comprising the sample has been stored at a temperature of about 10.degree. C. to about 40.degree. C. for a period of about 7 to about 30 days.

25. The method of claim 19, wherein the composition comprising the sample is stored at a temperature of (a) about minus 10.degree. C. to about 40.degree. C. for a period of about 24 to about 96 hrs; (b) about 0.degree. C. to about 35.degree. C. for a period of about 1 day to about 7 days; or (c) about 10.degree. C. to about 30.degree. C. for a period of about 7 days to about 30 days.

26. The method of claim 19, further comprising the step of analyzing or quantitating the purified population of polynucleotides obtained from the sample.

27. The method of claim 26, wherein the step of analyzing the purified population of polynucleotides obtained from the sample further comprises quantitating the amount of, or measuring the sequence fidelity of, the isolated single-stranded ribonucleic acid molecule of claim 1 contained within the composition.

28. The method of claim 27, wherein quantitating the total amount of, or measuring the sequence fidelity of the isolated single-stranded ribonucleic acid molecule of claim 1 involves determination of the PCR cycle threshold (C.sub.T) of the molecule.

29. An isolated, single-stranded ribonucleic acid molecule of at least 40 nucleotides in length comprising a nucleic acid segment that comprises: (a) a first sequence domain that specifically binds to a labeled probe of from about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length; and (d) the nucleotide sequence, prepared by the in vitro transcription of a polynucleotide comprising at least 30 contiguous nucleotides from the sequence of SEQ ID NO 1; wherein the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion.

30. An isolated, single-stranded ribonucleic acid molecule of at least 40 nucleotides in length comprising: a nucleic acid segment that comprises the sequence of SEQ ID NO 2 or the sequence of the in vitro transcription of SEQ ID NO 1 and a) does not substantially bind under stringent hybridization conditions to a mammalian genome, or to the genome of a bacterium, fungus, protozoan, virus that is pathogenic to a mammal, and b) comprises: (i) a first sequence domain that specifically binds to a labeled probe of from about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (ii) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; and (iii) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length, wherein the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion, and further wherein a population of the single-stranded ribonucleic acid molecules remains substantially non-degraded when stored at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days in a solution comprising a sufficient amount of each of (1) one or more chaotropes; (2) one or more detergents; (3) one or more reducing agents; (4) one or more chelators; and (5) one or more buffers.

31. An isolated, single-stranded ribonucleic acid molecule of claim 30, wherein at least 85% of a population of the single-stranded ribonucleic acid molecules remain substantially non-degraded when stored at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days in the solution.

32. The isolated, single-stranded ribonucleic acid molecule of claim 1, of about 60 to about 100 nucleotides in length.

33. The isolated, single-stranded ribonucleic acid molecule of claim 1, prepared by the in vitro transcription of a polynucleotide comprising at least 40 contiguous nucleotides from the sequence of SEQ ID NO 1.

34. The isolated, single-stranded ribonucleic acid molecule of claim 1, prepared by the in vitro transcription of a polynucleotide comprising at least 60 contiguous nucleotides from the sequence of SEQ ID NO 1.

35. The isolated, single-stranded ribonucleic acid molecule of claim 1, prepared by the in vitro transcription of a polynucleotide comprising at least 70 contiguous nucleotides from the sequence of SEQ ID NO 1.

36. The isolated, single-stranded ribonucleic acid molecule of claim 1, comprising at least 45 contiguous nucleotides from the sequence of SEQ ID NO 2.

37. The isolated, single-stranded ribonucleic acid molecule of claim 1, comprising at least 55 contiguous nucleotides from the sequence of SEQ ID NO 2.

38. The isolated, single-stranded ribonucleic acid molecule of claim 1, comprising at least 65 contiguous nucleotides from the sequence of SEQ ID NO 2.

39. An aqueous composition comprising: (1) about 1 M to about 4 M guanidine thiocyanate; (2) about 1 to about 100 mM sodium citrate; (3) about 0.25% to about 0.50 (wt./vol.) sodium dodecyl sulfate (SDS) or N-lauroyl sarcosine (NLS), sodium salt; (4) about 0.1 M 2-mercaptoethanol, or about 0.5 mM to about 10 mM TCEP; (5) about 0.0001% to about 0.1% silicone polymer (wt./vol.); and (6) an isolated, single-stranded ribonucleic acid molecule of at least 40 nucleotides in length comprising a nucleic acid segment that is not ordinarily found within a mammalian genome, or a genome of a bacterium, fungus, protozoan, virus that is pathogenic to a mammal, and that comprises: (a) a first sequence domain that specifically binds to a labeled probe of about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; and (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length, wherein: the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion; the ribonucleic acid molecule remains substantially full-length when placed in a solution comprising: (i) one or more chaotropes; (ii) one or more detergents; (iii) one or more reducing agents; (iv) one or more chelators; and (v) one or more buffers, each present in the solution in an amount that does not degrade a population of polynucleotides comprising the single-stranded ribonucleic acid molecule, when stored at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days; and the nucleotide sequence that is not ordinarily found within a mammalian genome comprises at least 30 contiguous nucleotides from SEQ ID NO:1, or at least 20 contiguous nucleotides from SEQ ID NO:2, or a combination thereof.

40. The aqueous composition of claim 39, further comprising about 0.1 mM to about 1 mM of a chelator selected from APCA, EDTA, EGTA, HEDTA, DTPA, NTA, and any combination thereof.

41. The aqueous composition of claim 39, further comprising about 10% to about 25% ethanol (vol./vol.).

42. The aqueous composition of claim 39, comprising: (1) about 1 M to about 4 M guanidine thiocyanate; (2) about 1 mM to 30 mM sodium citrate; (3) about 0.1% to about 1% (wt./vol.) sodium dodecyl sulfate (SDS) or N-lauroyl sarcosine (NLS), sodium salt; (4) about 0.5 mM to about 10 mM TCEP; (5) about 0.001% to about 0.0001% (wt./vol.) of a silicone polymer, (6) an isolated, single-stranded ribonucleic acid molecule of about 40 to 500 nucleotides in length comprising a nucleic acid segment that is not ordinarily found within a mammalian genome, or a genome of a bacterium, fungus, protozoan, virus that is pathogenic to a mammal, and that comprises: (a) a first sequence domain that specifically binds to a labeled probe of about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; and (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length, (7) about 10 mM to about 500 mM TRIS, (8) about 0.1 mM to about 1 mM APCA, EDTA, EGTA, HEDTA, DTPA, or NTA; and (9) about 10% to about 25% ethanol (vol./vol.).

43. An isolated deoxyribonucleic acid molecule of at least 40 nucleotides in length comprising a nucleic acid segment that comprises: (a) a first sequence domain that specifically binds to a labeled probe of from about 12 to about 35 nucleotides in length that is specific for the detection of the nucleic acid segment; (b) a second sequence domain that specifically binds to a forward PCR amplification primer of about 15 to about 35 nucleotides in length; (c) a third sequence domain that specifically binds to a reverse PCR amplification primer of about 15 to about 35 nucleotides in length; and (d) a nucleotide sequence that is not ordinarily found within a mammalian genome, or a genome of a bacterium, protozoan, fungus, virus that is pathogenic to a mammal; wherein: the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify the at least a first portion; the deoxyribonucleic acid remains substantially non-degraded when placed in a solution comprising: (i) one or more chaotropes; (ii) one or more detergents; (iii) one or more reducing agents; (iv) one or more chelators; and (v) one or more buffers, collectively present in an amount that does not degrade a population of polynucleotides comprising the single-stranded ribonucleic acid molecule, at a temperature of about 10.degree. C. to about 40.degree. C., for a period of about 1 day to about 90 days, and the nucleic acid sequence is at least about 80% identical to an at least 30 contiguous nucleotide sequence from the sequence of SEQ ID NO:1.

44. A method for determining the stability of a population of polynucleotides comprising: (a) placing a sample containing the population of polynucleotides into a one-step formulation including: (i) a known amount of an isolated nucleic acid molecule of at least 30 contiguous nucleotides from the sequence of SEQ ID NO 1; (ii) one or more chaotropes; one or more detergents; one or more reducing agents; one or more chelators; and one or more buffers, collectively, present in an amount sufficient to denature proteins, inactivate nucleases, kill pathogens, and not degrade nucleic acids in a sample suspected of containing nucleic acids, when the sample is contacted with the formulation; (b) storing the population of polynucleotides in the one-step formulation; (c) detecting the presence of the population of polynucleotides and the isolated nucleic acid molecule by one or more amplification-based strategies; (d) analyzing or quantitating the isolated nucleic acid molecule and the population of polynucleotides obtained from the sample; and (e) comparing the known amount of the isolated nucleic acid molecule to the amount of isolated nucleic acid present in the formulation for ascertaining the amount of nucleic acid degradation.

45. The method of claim 44, wherein at least 75% of the polynucleotides obtained from the sample are substantially non-degraded when the sample is stored for a period of about 7 to about 14 days in the one-step formulation substantially at a temperature of about minus 20.degree. C. to about 40.degree. C.

46. A method for determining the fidelity of the isolation of a population of polynucleotides comprising: (a) placing a sample containing the population of polynucleotides into a one-step formulation including: (i) an isolated nucleic acid molecule of at least 30 contiguous nucleotides from the sequence of SEQ ID NO 1 or at least 25 contiguous nucleotides from the sequence of SEQ ID NO 2; and (ii) one or more chaotropes; one or more detergents; one or more reducing agents; one or more chelators; and one or more buffers, collectively present in an amount sufficient to denature proteins, inactivate nucleases, kill pathogens, and not degrade nucleic acids in a sample suspected of containing nucleic acids, when the sample is contacted with the formulation; (b) storing the population of polynucleotides in the one-step formulation; (c) detecting the presence of the population of polynucleotides and the isolated nucleic acids by one or more amplification-based strategies; (d) measuring the amount of an isolated nucleic acid amplification product produced, wherein a reduction in the amount of the isolated nucleic acid amplification product produced relative to the amount of the isolated nucleic acid amplification product produced under similar conditions in a one-step formulation without the sample containing the population of polynucleotides indicates a lower amplification

efficiency; (e) measuring the amount of an amplification product from the population of polynucleotides and adjusting it to compensate for the reduction in the isolated nucleic acid amplified.

Description

FIELD OF THE INVENTION

The invention relates to compositions and methods for collection, transport, and storage of a biological sample containing a population of nucleic acids, which can subsequently be isolated, purified and/or characterized using one or more conventional methods or assays.

BACKGROUND OF THE INVENTION

In the field of molecular and diagnostic analysis, the ability to keep nucleic acids in a biological sample (and in particular, those contained in diagnostic samples obtained from human patients) stable, whether the specimen is taken in a remote field location, a doctor's office or in a laboratory, often determines whether the nucleic acids can be successfully analyzed. Nucleic acids in a biological sample quickly degrade and/or denature at room temperature and must generally be stored under freezing temperatures to remain stable; however, some degree of degradation still occurs over time. This problem is magnified when a specimen is collected at a remote field site, or a significant distance from a doctor's office or laboratory environment, and especially where there may be limited, or no, access to consistent and constant cooler/refrigerator/freezer conditions until the sample is analyzed, such as where access to electrical power, or refrigerator/freezer equipment is either unreliable, or non-existent. The problems associated with the collection and handling of biological specimens from which it is desirable to obtain nucleic acids are further exacerbated when the desired nucleic acids for downstream analysis include ribonucleic acid (RNA), which is particularly susceptible to degradation by endogenous or exogenous nuclease activity. Commercially-available specimen transport methods often use special transport media for biological samples for transport from point of collection to point of analysis, and in particular, packaging that imposes limitations on the time a sample may be stored, requires continual sample maintenance at low temperatures, even during transport or extended storage, and practically limits the time and distances possible between the collection site and the diagnostic laboratory.

In addition to concerns regarding specimen stability, there are often additional concerns regarding the handling and/or storage of reagents used in storing and transporting the collected samples. For example, the reagents themselves frequently require cold temperatures or other special care to maintain stability. Due to these stability issues, for example, transport of the reagents to a field site, storage at the field site before use, and transport of the biological specimens and reagents back to a testing site is a primary concern.

Another significant concern when working with biological specimens is the potential inoculation, release, or dissemination of live infectious pathogens or biological agents from the specimen into the environment. Specific protocols currently exist that are employed when handling samples that may be infectious or otherwise pose health or safety risks. If the sample is kept viable and/or biologically intact to preserve its integrity for testing, individuals involved in the collection, transfer, and testing process are potentially exposed to highly dangerous contagions. Additionally, innocent bystanders nearby a field or sample collection site (or nearby during transport) can be exposed if a release of the contagion occurs. As a result,

the required safety measures typically increase the expense and effort required to move such samples from one location to another.

Until recently, clinical laboratory methods for pathogen detection were labor-intensive, expensive processes that required highly knowledgeable and expert scientists with specific experience. The majority of clinical diagnostic laboratories employed the use of traditional culturing methods that typically require 3 to 7 days for a viral culture--and even longer for some other bacterial targets. Furthermore, traditional culturing requires collection, transport, and laboratory propagation and handling of potentially infectious biological agents such as Ebola, avian influenza, severe acute respiratory syndrome (SARS), etc.

The field of clinical molecular diagnostics changed drastically with the advent of polymerase chain reaction (PCR) in the mid eighties, however, and shortly thereafter with real-time PCR in the mid 90's. Nucleic-acid based detection platforms employing e.g., quantitative real-time PCR (qPCR) or reverse transcriptase PCR (RT-PCR) and quantitative, real-time, reverse transcriptase PCR (qRT-PCR) assays can deliver results in hours versus days required for traditional culture and isolation methods making molecular detection methods the mainstay of modern diagnostic laboratory analysis. Recent improvements in detection chemistries, such as new and improved reporting/quenching fluors, minor groove binders (MGB) (TaqMan MGB.TM., Applied Biosystems; for a general reference, see also e.g., Baraldi et al., *Pure Appl. Chem.*, 75 (2-3):187-194, 2003), and stabilized amplification reagents have paved the way for more sensitive and highly-specific nucleic acid detection assays, and have proved more timely, robust, and economical than antiquated cell culture-based methods. Advances in other nucleic acid detection strategies such as transcription-mediated amplification, ligase chain reaction (LCR), and so-called "laboratory-on-a-chip" multiplexed assays, have also contributed to the transition from culture vials to microarrays in the clinical laboratory.

Several commercial companies (e.g., Qiagen [Valencia, Calif., USA], Roche Applied Science [Indianapolis, Ind., USA], Gen-Probe [San Diego, Calif., USA], and bioMerieux [Durham, N.C., USA]) have developed instruments to automate the nucleic acid extraction process from sample isolation to molecular analysis. For example, the Tigris DTS.RTM. (Gen-Probe, San Diego, Calif., USA) automates the entire detection process, and in late 2004 was approved by the U.S. Food and Drug Administration (FDA) for simultaneously detecting *Chlamydia trachomatis* and *Neisseria gonorrhoeae* using Gen-Probe's APTIMA COMBO-2.RTM. amplified nucleic acid test (NAT) assay.

In view of the requirement for high-quality nucleic acid samples in contemporary detection and assay systems, there is now a need in the art for safe and facile collection, storage and transport systems that maintain the integrity and quality of nucleic acids contained within a variety of biological samples and specimens. Moreover, there is also now a need for collecting, preserving, and transporting samples (and particularly those containing harmful or pathogenic organisms) in remote or field locations under ambient environmental (i.e., non-ideal) conditions for extended periods of time without refrigerating, freezing, or otherwise special handling of the collection reagent(s), the biological sample itself, or the population of nucleic acids contained therein. Furthermore, there is now a need for less-expensive, and more-convenient collection/transport/storage media that minimize risk of pathogen exposure to workers or innocent bystanders, allow for the use of single-step formulations, and facilitate convenient ambient transportation of biological specimens over long distances or extended periods of time that contain high-fidelity, high-quality nucleic acid populations.

BRIEF SUMMARY OF THE INVENTION

The present invention encompasses new and useful compositions, as well as methods of making and employing them, that may advantageously improve conventional collection, lysis, transport and storage methods for the preparation of nucleic acids from one or more biological sources. Accordingly, the present invention advantageously provides a collection and preservation formulation to inactivate and lyse a biological specimen containing nucleic acids, and preserve nucleic acids (e.g., RNA or DNA) within the biological specimen, preferably all in a single reaction vessel, such that the integrity of the nucleic acids is at least substantially maintained, and preferably entirely maintained, so that a portion of the nucleic acids are readily available for molecular diagnostic analysis. An additional advantage of the present invention is that the formulation can enable the separated or released nucleic acids to remain at least substantially stable, without requiring consistent and constant cooler temperatures, such as refrigeration or freezing.

The one-step formulations disclosed herein accomplish the following principal functions: inactivation or killing of pathogens within the sample; lysis of cells and separation or release of nucleic acids from the cells; inactivation of endogenous or exogenous nucleases and other cellular enzymes to prevent degradation of the nucleic acids present in the sample; and facilitation of collection and handling of the sample at ambient temperatures, stabilization of the nucleic acids during subsequent transport and storage of the sample, preservation of sample integrity through the use of extra-genomic, nonsense RNA carrier molecules that also serve as an internal positive control (IPC) to monitor fidelity of the processed samples.

The ability to achieve all of these desirable functions in a single-step formulation, preferably in a single reaction zone or reaction vessel, is a particularly marked advantage over that presently available. Presently existing technologies do not include a single-step composition that provides for inactivation of biological components containing nucleic acids, release of nucleic acids through lysis of cells and separation or release of nucleic acids, maintenance of the integrity of the liberated population of nucleic acids, and convenient IPCs to quantitate specimen fidelity. The present invention both stabilizes and preserves the integrity of nucleic acids present in the specimen for diagnostic testing, while also providing convenient methods for monitoring the fidelity of the collection medium that employ a readily-detectable, quantifiable RNA carrier molecule.

The one-step formulations of the present invention allow for preferably simultaneous inactivation of biological components containing nucleic acids, lysis and separation or release of nucleic acids, stabilization, and preservation. In one embodiment, some or all of the inactivation, lysis and separation or release, stabilization, and preservation, are sequential. In a preferred embodiment, however, a majority or preferably all of these functions occur simultaneously. In all embodiments, the one-step formulation is combined with the sample to initiate these functions. This is in contrast to previous technology in which inactivation did not necessarily occur, and lysis, stabilization, and preservation occurred in a succession of separate steps, each step typically using one or more distinct reagents and protocols that were separately added.

Unlike prior nucleic acid isolation methods that required numerous sequential steps to minimize errors, avoid reagent incompatibility, and provide stepwise control of results, the present invention provides all these benefits in a single one-step formulation, and also adds the further benefits of maintaining the integrity of the nucleic acids, monitoring the fidelity of the storage/collection/transport process through the use of an internal nonsense RNA carrier molecule, facilitating both higher extraction and purification efficiencies of the sample nucleic acids, and improving their ultimate yield. The one-step formulations of the present invention preferably facilitate simultaneous inactivation of biological components containing nucleic acids, lysis and release of nucleic acids from cellular debris, stabilization, and preservation of nucleic acids reduces the chance for degradation of the population of nucleic acids contained within the biological sample that may

occur before, during, or after lysis.

The disclosed compositions were developed and optimized: 1) to facilitate preparation of high-quality nucleic acids from clinical or environmental specimens, 2) to inactivate, kill, or otherwise neutralize potentially infectious pathogens in a biological sample to facilitate safe handling and transport of the collected specimens, and 3) to stabilize released (i.e., `naked`) DNA/RNA for prolonged periods without hydrolysis or nuclease degradation of the released nucleic acids.

The compositions described herein are ideal for clinical, field and deployment use, or for high volume sample collection/extraction. Specimens collected in one or more of the disclosed compositions are biologically inactivated, and may be safely shipped, typically even without refrigeration or dry ice.

In certain embodiments, the addition of nucleic acids (e.g., RNA and/or DNA) is contemplated to be beneficial for a variety of purposes and applications of the disclosed methods: a) as a "carrier" (The addition of small amounts of supplemental RNA/DNA has been previously been shown to augment/increase the overall yield of samples/specimens, particularly original specimens that may contain low amounts of target, i.e., cells, viruses, bacteria); b) as an IPC for downstream molecular processes and to track or monitor the fidelity of the nucleic acid preparation from sample collection to detection; and c) for comparison to a `calibrator` for downstream quantitative analysis, e.g., qRT-PCR and the like. In such embodiments, one or more known or "control" nucleic acids could be added to the compositions in a final concentration of from about 1 ag to about 1 mg, more preferably from about 1 fg to about 1 .mu.g, and more preferably still, from about 1 pg to about 1 ng.

In an illustrative embodiment, the invention provides an isolated single-stranded (ss) or double-stranded (ds) RNA, DNA, PNA, or hybrid thereof that is useful: (a) as a carrier molecule for aiding in the recovery of polynucleotides from a biological sample suspected of containing nucleic acids, and/or (b) as an IPC (i.e., a "known," "reporter," "control," "standard," or "marker") sequence to monitor the integrity and fidelity of specimen collection and polynucleotide isolation/stabilization. In certain embodiments, the invention provides an isolated ds-RNA, ds-DNA, ds-PNA, or a hybrid thereof that is useful as a carrier molecule and/or an IPC. In other embodiments, the invention provides an isolated ssRNA, ssDNA, ssPNA, or a hybrid thereof that is useful as a carrier molecule and/or as an IPC sequence. In exemplary embodiments, the invention provides an isolated ssRNA molecule that is useful as both a carrier molecule and an IPC sequence.

Such molecules can be isolated from natural sources, prepared in the laboratory, or alternatively, a hybrid containing both native- and non-native sequences. As noted herein, because the compositions of the invention are particularly useful for the isolation and characterization of biological specimens obtained from mammalian (and in particular, human) sources that are suspected of containing polynucleotides of pathogen-origin, it is preferable that the sequence(s) employed as carrier and/or positive control compounds substantially contain a primary nucleotide sequence that is not ordinarily found within the genome of a mammal, or within the genome of an organism that is pathogenic to such a mammal. Exemplary mammals include, without limitation, bovines, ovines, porcines, lupines, canines, equines, felines, ursines, murines, leonines, leporines, hircines, and non-human primates.

Preferably, this non-mammalian, non-pathogen-specific carrier/reporter sequence is not cross-reactive, i.e., does not substantially, or preferably, do(es) not, hybridize to, mammalian or pathogen-specific sequences, and as such, non-coding, non-degenerate (i.e., nonsense) sequences are particularly preferred in the formulation of control/carrier sequences to minimize hybridization of the control/carrier sequence to a member of the

isolated population of polynucleotides obtained from the collected specimen. Exemplary carrier/control sequences therefore, do not substantially, or preferably, do(es) not, bind (e.g., hybridize under stringent hybridization conditions) to a population of polynucleotides isolated from a mammalian genome, or to a population of polynucleotides isolated from the genome of a bacterium, fungus, virus that is pathogenic to a mammal. Exemplary stringent hybridization conditions known to those of ordinary skill in the art include, without limitation, (a) pre-washing in a solution containing about 5.times.SSC, 0.5% SDS, and 1.0 mM EDTA (pH 8.0); (b) hybridizing at a temperature of from about 60.degree. C. to about 70.degree. C. in 5.times.SSC overnight; and (c) subsequently washing at about 65 to about 70.degree. C. for 20 min. with each of 2.times., 0.5.times. and 0.2.times.SSC containing 0.1% SDS), or equivalent hybridization conditions thereto.

While any practical length of carrier molecule may be utilized in the compositions of the invention, and in particular, in the formulations (e.g., PSS formulations) disclosed herein, preferably, the length of the carrier molecule is about 40 to about 900 nucleotides in length; more preferably, about 50 to about 800, or alternatively about 60 to about 700, nucleotides in length; still more preferably, about 70 to about 600, or alternatively, about 80 to about 500 nucleotides in length; and even more preferably still, about 90 to about 400 nucleotides, or alternatively, about 100 to about 300 or so nucleotides in length.

As such, all intermediate lengths of polynucleotides are contemplated to fall within the scope of the present disclosure, including, without limitation, carrier molecule sequences that are about 45 nucleotides in length, about 50 nucleotides in length, about 55 nucleotides in length, about 60 nucleotides in length, about 65 nucleotides in length, about 70 nucleotides in length, about 75 nucleotides in length, about 80 nucleotides in length, about 85 nucleotides in length, about 90 nucleotides in length, about 95 nucleotides in length, about 100 nucleotides in length, about 120 nucleotides in length, about 140 nucleotides in length, about 160 nucleotides in length, about 180 nucleotides in length, about 200 nucleotides in length, about 220 nucleotides in length, about 240 nucleotides in length, about 260 nucleotides in length, about 280 nucleotides in length, or even about 300 nucleotides or so in length.

While any practical, nonsense/minimally cross-hybridizing polynucleotide sequence is contemplated to be useful in the preparation of a carrier sequence to facilitate improved isolation of the target population of polynucleotides, in circumstances where a control sequence is desired, additional requirements of the primary nucleic acid sequence of the molecule are imposed, namely that the sequence must be detectable and/or quantifiable using one or more conventional molecular biology techniques.

To that end, the invention provides carrier/reporter molecules that contain at least a first sequence domain that specifically hybridizes (i.e., binds) to a suitably-detectable probe, including, without limitation, molecularly-labeled probes and derivatives thereof. Exemplary labeled probes are those that include radioactive, luminescent, chemiluminescent, fluorescent, enzymatic, magnetic, or spin-resonance labels known to those of ordinary skill in the molecular arts. In illustrative embodiments, the labeled probe contains at least a first minor groove binder.

In certain embodiments, to facilitate the binding of conventional detectable-label probes, the carrier/reporter molecules of the invention will contain at least a first sequence domain of from about 10 to about 60 nucleotides in length that specifically binds to at least a first detectable probe.

While the first sequence domain may be of any practical length within the entirety of the carrier sequence, preferably, the first sequence domain will be from about 12 to about 50 nucleotides in length; more

preferably, from about 14 to about 45 nucleotides in length; still more preferably, from about 16 to about 40 or so nucleotides in length, and more preferably still, from about 18 to about 30 or so nucleotides in length.

As such, all intermediate lengths of probe-hybridizing sequence domains are contemplated to fall within the scope of the present disclosure, including, without limitation, probe-binding domains that are about 13 nucleotides in length, about 14 nucleotides in length, about 15 nucleotides in length, about 16 nucleotides in length, about 17 nucleotides in length, about 18 nucleotides in length, about 19 nucleotides in length, about 20 nucleotides in length, about 21 nucleotides in length, about 22 nucleotides in length, about 23 nucleotides in length, about 24 nucleotides in length, about 25 nucleotides in length, about 26 nucleotides in length, about 27 nucleotides in length, about 28, or even about 29 or 30 or so nucleotides in length.

Where it is desirable that the carrier/reporter sequence be detectable by one or more amplification-based methodologies (including, without limitation, PCR-based methodologies), additional molecular characteristics for the carrier/reporter sequence are desirable, namely that the sequence, or at least a portion thereof, be amplifiable using one or more polymerase chain-based assays. For PCR or FRET-based methodologies, for example, it is desirable that the carrier/reporter molecule contains at least a second sequence domain that specifically binds to a forward PCR amplification primer (typically of from about 20 to about 40 nucleotides in length) and a third sequence domain that specifically binds to a reverse PCR amplification primer of (also typically of from about 20 to about 40 nucleotides in length). Preferably, the second and third sequence domains are operably positioned to facilitate a PCR-directed amplification of at least a first portion of the nucleic acid segment from the forward and reverse primers under conditions effective to amplify at least a first portion.

When it is desirable that the carrier/reporter sequence be detectable by one or more non-symmetric amplification-based methodologies (including, without limitation, asymmetric PCR, molecular hybridization, affinity-label-based methodologies, and the like), it is not essential that a pair of forward and reverse primer binding domains are required--only a single second sequence domain may be sufficient to detect that carrier sequence using appropriate diagnostic assays. In such embodiments, a single detection/replication primer (or polymerase-binding domain) may be employed. Typically, such sequence domains will be from about 15 to about 35 or so nucleotides in length, although any conventional targeting domain may be employed using conventional analytical methods. The design and primary sequence of such domains is considered routine in the molecular arts, and as such, is within the purview of the ordinarily-skilled artisan.

Depending upon the length and primary nucleotide of the IPC sequence, in certain embodiments, the first sequence domain shares at least a first common nucleotide sequence region with the second, third, or even both the second and third sequence domains. The overlap of such sequences may be only one or more nucleotides, or alternatively, larger regions of overlap consisting of 5 to 10 or more nucleotides.

In exemplary embodiments, the IPC sequence may be prepared by one or more suitable molecular biology techniques, including, e.g., by the *in vitro* transcription of a polynucleotide that comprises the sequence of the control, or alternatively, comprises a nucleic acid sequence that is complementary to the sequence of the control itself. As described in the following examples, one method of preparing the control sequence involves the preparation of a DNA duplex, from which a single-stranded RNA control sequence may be synthesized.

Exemplary single-stranded RNA control sequences include, but are not limited to, those that contain at least about 105, at least about 95, at least about 85, at least about 75, at least about 65, at least about 55, at least about 45, at least about 35, or at least about 25 or so contiguous nucleotides from the RNA sequence: 5'-

CCCUUAGCAGCACGUCAGUCAGGGAGCCAAUUUCAGAGCUCAGCGAGACAGUUUUA
UAGGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUCCAAAGUUUCAGUU-3' (SEQ ID
NO:2). In illustrative embodiments, the isolated single-stranded ribonucleic acid molecule is preparable by a process that includes the in vitro transcription of a polynucleotide that comprises, consists essentially of, or alternatively, consists of, the nucleic acid sequence of SEQ ID NO:2.

In one such embodiment, the single-stranded RNA control sequence may be synthesized in vitro from a DNA sequence that contains at least about 90, at least about 80, at least about 70, at least about 60, at least about 50, at least about 40, at least about 30, or at least about 20 or fewer contiguous nucleotides from the sequence: 5'-

CCCUUAGCAGCACGUCAGUCAGGGAGCCAAUUUCAGAGCUCAGCGAGACAGUUUUA
UAGGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUCCAAAGUUUCAGUU-3' (SEQ ID
NO:1). In exemplary embodiments, the isolated single-stranded ribonucleic acid molecule is preparable by a process that includes the in vitro transcription of a polynucleotide that comprises, consists essentially of, or alternatively, consists of the nucleic acid sequence of SEQ ID NO:1.

In one aspect of the invention, there is provided a method for detecting the presence or absence of the RNA carrier molecule within a plurality of polynucleotides obtained from a biological sample prepared using a disclosed sample collection/storage formulation, including, without limitation, one or more formulations as described herein, including, without limitation, PSS formulations. There is also a method for quantitating the amount of RNA carrier molecule within the sample and monitoring the efficacy of the formulations at stabilizing and protecting the molecular fidelity of the isolated polynucleotides. By comparing the amount of carrier molecule remaining in the sample after an elapsed period of time to that of the known quantity of carrier molecule present in the initial solution, one may determine the amount of degradation the carrier molecule has undergone while in the presence of the biological sample. One may use this correlation as an estimation of the extent of degradation that has taken place over time in the population of polynucleotides liberated from the original sample.

In another aspect, the present invention provides a method for rapidly detecting in a biological sample, a particular polynucleotide sequence, such as that of the IPC sequence. In an overall and general sense, this method comprises amplification of a population of nucleotides suspected of containing the particular sequence using conventional methods such as PCR and forward and reverse primers that are specific for the target sequence, hybridization of a specific probe set with the resulting single-stranded PCR product, performing melting curve analysis and analyzing the $T_{sub.m}$ change of the hybrid of the single-stranded PCR product with the hybridization probes.

In one embodiment, the present invention provides a method for rapidly detecting the presence of a polynucleotide (such as a carrier IPC sequence), using a PCR-based methodology, which generally comprises the steps of: (a) isolating polynucleotides from the sample to be analyzed; (b) amplifying the polynucleotides by PCR using a primer set that is specific to the target sequence; (c) hybridizing one or more labeled probes that are specific for the polynucleotide of interest with the single-stranded PCR product obtained from step (b); and (d) detecting the presence of the labeled probe in the sample, indicative of the presence of the specific target sequence within the population of isolated polynucleotides.

One such method for the detection of polynucleotides using a labeled "probe" sequence utilizes the process of fluorescence resonance energy transfer (FRET). Exemplary FRET detection methodologies often involve pairs of fluorophores comprising a donor fluorophore and acceptor fluorophore, wherein the donor

fluorophore is capable of transferring resonance energy to the acceptor fluorophore. In exemplary FRET assays, the absorption spectrum of the donor fluorophore does not substantially overlap the absorption spectrum of the acceptor fluorophore. As used herein, "a donor oligonucleotide probe" refers to an oligonucleotide that is labeled with a donor fluorophore of a fluorescent resonance energy transfer pair. As used herein, "an acceptor oligonucleotide probe" refers to an oligonucleotide that is labeled with an acceptor fluorophore of a fluorescent resonance energy transfer pair. As used herein, a "FRET oligonucleotide pair" will typically comprise an "anchor" or "donor" oligonucleotide probe and an "acceptor" or "sensor" oligonucleotide probe, and such pair forms a FRET relationship when the donor oligonucleotide probe and the acceptor oligonucleotide probe are both hybridized to their complementary target nucleic acid sequences. Acceptable fluorophore pairs for use as fluorescent resonance energy transfer pairs are well known to those skilled in the art and include, but are not limited to, fluorescein/rhodamine, phycoerythrin/Cy7, fluorescein/Cy5, fluorescein/Cy5.5, fluorescein/LC Red 640, and fluorescein/LC Red 705, and the like.

Primers useful in amplification of a particular IPC sequence may be designed using, for example, a computer program such as OLIGO.RTM. (Molecular Biology Insights Inc., Cascade, Colo., USA). Typically, oligonucleotide primers are from about 10 to about 60 or so nucleotides in length (including, without limitation, all intermediate integers, e.g., 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, or even 60 or more nucleotides in length), although primers of any practical length may be useful in the practice of the invention.

The invention also provides a method for increasing the efficiency of obtaining a purified population of polynucleotides from a biological sample suspected of containing such polynucleotides. In an overall and general sense, the method includes contacting the sample with a composition that comprises one or more of the formulations described herein (including, without limitation, PSS formulations), and a carrier and/or IPC nucleic acid segment, in an amount and for a time sufficient to increase the efficiency of obtaining the purified population of polynucleotides from the biological sample.

In exemplary embodiments, the integrity of a population of polynucleotides in the biological sample, and/or the fidelity of at least a first sequence of at least one of the polynucleotides obtained from the sample is at least substantially maintained (i.e., at least 75% of the nucleotides within the population are substantially full-length) when the composition comprising the sample is stored at a temperature of from about -20.degree. C. to about 40.degree. C. for a period of from about 7 to about 14 days or longer; alternatively at a temperature of from about -20.degree. C. to about 40.degree. C. for a period of from about 7 to about 14 days or longer; or alternatively at a temperature of from about 20.degree. C. to about 40.degree. C. for a period of from about 14 to about 30 days or more.

Alternatively, the integrity of a population of polynucleotides in the biological sample is at least substantially maintained such that at least about 80%, at least about 85%, at least about 90%, or at least about 95% or more of the nucleotides within the population are present in the solution when compared to the amount present in the solution when the sample was initially collected. In preferred embodiments, the integrity of the sample will be substantially maintained such that all, or almost all of the polynucleotides present in the initial sample will be maintained (i.e., not detectably degraded) over time.

In the practice of the disclosed methods, preferably from the time of collection to the time of isolating, purifying, or characterizing a population of polynucleotides therein, less than about 20% of the population of polynucleotides originally present in the collected sample will be degraded over time during subsequent

storage. Preferably, substantially less than about 15% of the population of polynucleotides originally present in the collected sample will be degraded over time during subsequent storage, more preferably, less than about 10% of the population of polynucleotides originally present in the collected sample will be degraded over time during subsequent storage, and more preferably still, less than about 5% of the population of polynucleotides originally present in the collected sample will be degraded over time during subsequent storage. In particularly preferred embodiments, not more than about 5%, about 4%, about 3%, about 2% or about 1% of the population of polynucleotides originally present in the collected sample will be degraded over time during subsequent storage. Such high-integrity preservation of sample quality is preferable, regardless of the conditions under which the sample is stored, and will be substantially maintained for a period of time of at least about 7 days, at least about 14 days, at least about 21 days, at least about 30 days, at least about 45 days, at least about 60 days, or even at least about 90 days or more.

The disclosed methods for collecting, storing, and/or transporting samples of biological origin from which one or more populations of nucleic acids may be obtained, may also further optionally include one or more steps, including without limitation, (a) analyzing and/or quantitating one or more polynucleotides obtained from the sample; (b) detecting the presence of the carrier molecule (i.e., IPC) in the collected sample; (c) quantitating the amount of the IPC in the collected sample; and/or (d) measuring the sequence fidelity of, or determining the percent degradation of the IPC in the collected sample.

While the presence of, integrity of, or sequence fidelity of, a particular polynucleotide sequence obtained from, or utilized in the practice of the present invention may be determined, using any conventional methodology known to those of ordinary skill in the molecular arts, in one embodiment, PCR amplification is utilized. Likewise, determination of the integrity of a polynucleotide of interest may include determination of the PCR cycle threshold (C.sub.T) under given conditions, and determination of the sequence fidelity, qualitative integrity of collected nucleic acids may be determined by conventional DNA or RNA sequencing methods, including, without limitation, the chemical-based methods of Maxam-Gilbert, the dideoxy chain termination method of Sanger et al., the dye fluorophore-based method of Mathies et al., or pyrosequencing techniques as described by Nyren and Ronaghi.

Exemplary formulations of the invention include a one-step collection solution that lyses, stabilizes, and preserves the integrity of nucleic acids prepared from a biological sample for subsequent RNA and/or DNA isolation, detection, quantitation, amplification, and/or analysis.

In certain embodiments, the nonsense, carrier/IPC compositions of the invention will be formulated in a one-step sample collection/storage/transport medium that includes: a) one or more chaotropes (each preferably present in the composition an amount from about 0.5 M to about 6 M); b) one or more detergents (each preferably present in the composition an amount from about 0.1% to about 1%); c) one or more chelators (each preferably present in the composition in an amount from about 0.01 mM to about 1 mM); d) one or more reducing agents (each preferably present in the composition in an amount from about 0.05 M to about 0.3 M); and e) one or more defoaming agents (each preferably present in the composition in an amount from about 0.0001% to about 0.3%).

Exemplary chaotropes include, without limitation, guanidine thiocyanate (GuSCN), guanidine hydrochloride (GuHCl), guanidine isothionate, potassium thiocyanate (KSCN), sodium iodide, sodium perchlorate, urea, or any combination thereof. Descriptions of additional exemplary chaotropes and chaotropic salts can be found in U.S. Pat. No. 5,234,809 (specifically incorporated herein in its entirety by express reference thereto).

Exemplary detergents include, without limitation, sodium dodecyl sulfate (SDS), lithium dodecyl sulfate (LDS), sodium taurodeoxycholate (NaTDC), sodium taurocholate (NaTC), sodium glycocholate (NaGC), sodium deoxycholate (NaDC), sodium cholate, sodium alkylbenzene sulfonate (NaABS), N-lauroyl sarcosine (NLS), salts of carboxylic acids (i.e., soaps), salts of sulfonic acids, salts of sulfuric acid, phosphoric and polyphosphoric acid esters, alkylphosphates, monoalkyl phosphate (MAP), and salts of perfluorocarboxylic acids, anionic detergents including those described in U.S. Pat. No. 5,691,299 (specifically incorporated herein in its entirety by express reference thereto), or any combination thereof.

Exemplary reducing agents include, without limitation, 2-mercaptoethanol (.beta.-ME), tris(2-carboxyethyl)phosphine (TCEP), dithiothreitol (DTT), formamide, dimethylsulfoxide (DMSO), or any combination thereof. In a preferred embodiment, the reducing agent includes or is TCEP.

Exemplary chelators include, without limitation, ethylene glycol tetraacetic acid (EGTA), hydroxyethylethylenediaminetriacetic acid (HEDTA), diethylene triamine pentaacetic acid (DTPA), N,N-bis(carboxymethyl)glycine (NTA), ethylenediaminetetraacetic (EDTA), citrate anhydrous, sodium citrate, calcium citrate, ammonium citrate, ammonium bicitrate, citric acid, diammonium citrate, potassium citrate, magnesium citrate, ferric ammonium citrate, lithium citrate, or any combination thereof. In preferred embodiments, the chelator includes EDTA, a citrate, or a combination thereof. In a more preferred embodiment, the chelator includes EDT.

The compositions of the invention can further include a defoaming agent to prevent the formation of bubbles that typically result from the presence of detergents in the formulation. Defoaming agents facilitate pipetting and handling of the disclosed compositions. Exemplary surfactants/defoaming agents include, without limitation, cocoamidopropyl hydroxysultaine, alkylaminopropionic acids, imidazoline carboxylates, betaines, sulfobetaines, sultaines, alkylphenol ethoxylates, alcohol ethoxylates, polyoxyethylenated polyoxypropylene glycols, polyoxyethylenated mercaptans, long-chain carboxylic acid esters, alkonolamides, tertiary acetylenic glycols, polyoxyethylenated silicones, N-alkylpyrrolidones, alkylpolyglycosidases, silicone polymers such as Antifoam A.RTM., or polysorbates such as Tween.RTM., or any combination thereof. In a preferred embodiment, a defoaming agent includes a silicone polymer.

Optionally, the compositions of the invention may further include one or more buffers (each preferably present in the final composition in an amount from about 1 mM to about 1 M). Exemplary buffers include, without limitation, tris(hydroxymethyl)aminomethane (Tris), citrate, 2-(N-morpholino)ethanesulfonic acid (MES), N,N-bis(2-hydroxyethyl)-2-aminoethanesulfonic acid (BES), 1,3-bis(tris(hydroxymethyl)methylamino)propane (Bis-Tris), 3-(cyclohexylamino)-1-propanesuhinic acid (CAPS), 3-(cyclohexylamino)-2-hydroxy-1-propanesuhic acid (CAPSO), 4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid (HEPES), 3-(N-morpholino)propanesulfonic acid (MOPS), 3-(N-morpholino)-2-hydroxypropanesulfonic acid (MOPSO), piperazine-N,N'-bis(2-hydroxypropanesulfonic acid (POPSO), N-[Tris(hydroxymethyl)methyl]-3-amino propanesulfonic acid (TAPS), N-[Tris(hydroxymethyl)methyl]-3-amino-2-hydroxypropansulfonic acid (TAPSO), N-[Tris(hydroxymethyl)methyl]-2-aminoethanesulfonic acid (TES), N,N-bis(2-hydroxyethyl)glycine (Bicine), N-[tris(hydroxymethyl)methyl]glycine (Tricine), N-2-acetamido-2-iminodiacetic acid (ADA), N-(2-acetamido)-2-aminoethanesulfonic acid (ACES), piperazine-1,4-bis(2-ethanesulfonic acid) (PIPES), bicarbonate, phosphate, or any combination thereof. In a preferred embodiment, the buffer includes a citrate.

The inclusion of one or more of such optional but preferred buffers is desirable to control the pH of the formulations, since it has been found that nucleic acid extraction is optimal in a pH range of about 5 to 7.

Preferably, the one or more buffers employed in the disclosed compositions are chosen to provide a significant buffering capacity in the range from a pH of about 6 to a pH of about 8, more preferably within a pH range of about 6 to about 7, and more preferably still, within a pH range of about 6.2 to about 6.8. In exemplary embodiments, the pH of a PrimeStore.TM. solution (also referred to herein as "PSS") is preferably about 6.9.+-.0.25.

The compositions of the invention may also further optionally include one or more short-chain (preferably from 1- to 6-carbon [i.e., C.sub.1-C.sub.6] alcohols) alkanols (each preferably present in the composition in an amount from about 1% to about 25%, although higher percentages of the alcohols may be employed if desired). Exemplary short-chain alkanols include linear and branched-chain alcohols, such as, without limitation, methanol, ethanol, propanol, butanol, pentanol, hexanol, or any combination thereof.

The compositions of the invention may also further optionally include one or more additional compounds or reagents including, without limitation, cationic functionalized zwitterionic compounds such as betaines (including, without limitation, N,N,N-trimethylglycine, cocamidopropyl betaine, and the like), albuminoids (including, without limitation, ovalbumin, and the serum albumins of bovine, equine, or human origin), and osmolytes (including, without limitation, trimethylamine N-oxide (TMAO), dimethylsulfoniopropionate, sarcosine, and saccharides or sugar alcohols including, without limitation, trehalose, maltose, rhamnose, sucrose, arabinose, fucose, mannitol, sorbitol, adonitol, and the like).

Preferably, the compositions of the invention provide sufficient buffering capacity to adequately stabilize the populations of polynucleotides obtained from a sample, and will, most preferably, be buffered to a pH of about 6.4 to 6.9 during formulation, and will maintain the isolated populations of polynucleotides in a similar pH range when the sample is contacted with the storage/collection formulations described herein.

The compositions of the present invention will typically at least substantially inactivate, and preferably entirely inactivate, any endogenous or exogenous RNAses or DNAses present in the sample, such that the nucleic acids of the sample are substantially free of any degradation, and preferably do not degrade, or lose integrity, during the collection, lysis, storage, and transport of the sample for subsequent in vitro or in vivo analyses.

Exemplary formulations of the storage/transport/collection compositions of the invention are described in the examples herein, and include, without limitation, a composition that includes about 4 M of a chaotrope (such as guanidine thiocyanate, guanidine hydrochloride, guanidine isocyanate, or any combination thereof), about 10 mM to 30 mM of a chelator (such as EGTA, HEDTA, DTPA, NTA, EDTA, citrate anhydrous, sodium citrate, calcium citrate, ammonium citrate, ammonium bicitrate, citric acid, diammonium citrate, ferric ammonium citrate, lithium citrate, or any combination thereof), about 0.25% of a detergent (such as SDS, LDS, NaTDC, NaTC, NaGC, NaDC, sodium cholate, NaABS, NLS, or any salt or combination thereof), about 0.1 M of a reducing agent (such as .beta.-ME, DTT, DMSO, formamide, TCEP, or any combination thereof), and about 0.1% of a surfactant/defoaming agent (such as a silicone polymer [e.g., Antifoam A.RTM.g] or a polysorbate [e.g., Tween.RTM.], or any combination thereof).

Additional exemplary formulations of the specimen collection compositions of the invention include, without limitation, a composition that includes about 3 M of a chaotrope (such as guanidine thiocyanate, guanidine hydrochloride, guanidine isocyanate, or any combination thereof), about 1 mM of 0.5 M reducing agent (such as, e.g., .beta.-ME, TCEP, formamide, DTT, DMSO, or any combination thereof), about 1 to about 10 mM of a chelator (such as, e.g., EGTA, HEDTA, DTPA, NTA, EDTA, citrate anhydrous, sodium citrate, calcium

citrate, ammonium citrate, ammonium bicitrate, citric acid, diammonium citrate, ferric ammonium citrate, lithium citrate, or any combination thereof), about 0.25% of a detergent (such as SDS, LDS, NaTDC, NaTC, NaGC, NaDC, sodium cholate, NaABS, NLS, or any salt or combination thereof), and optionally but preferably about 0.0002% of a defoaming agent (also referred to as an antifoaming agent) (such as a silicone polymer or a polysorbate, or any combination thereof) and about 100 mM of a buffer (such as Tris, MES, BES, Bis-Tris, HEPES, MOPS, bicarbonate, citrate, phosphate, or any combination thereof).

Another exemplary formulation of the disclosed polynucleotide isolation and stabilization compositions include, without limitation, a composition that includes about 1 to about 4 M of a chaotropic agent such as guanidine thiocyanate, guanidine hydrochloride, or guanidine isocyanate; about 0.5 to 100 mM of a chelating agent such as EDTA, or sodium citrate, or both; about 0.1 to about 1% of an anionic detergent such as SDS or N-lauroyl sarcosine, sodium salt; about 0.001 to about 0.0001% of a surfactant or wetting agent such as the silicone polymer, Antifoam A.RTM., e); about 10 to about 500 mM of a buffering agent such as Tris-HCl; and about 10 to about 25% of a short-chain alkanol such as ethanol.

In particular embodiments, the invention provides a composition that includes about 2.5 M guanidine thiocyanate; about 0.5 mM TCEP; about 10 mM sodium citrate; about 0.4% N-lauroyl sarcosine, sodium salt; about 0.0002% Antifoam A, about 100 mM Tris-HCl, about 0.1 mM EDTA; and about 23% ethanol.

The invention also provides a method for obtaining a population of polynucleotides from a sample suspected of containing nucleic acids. The method generally involves associating the sample with an amount of one of the disclosed compositions, under conditions effective to obtain a population of polynucleotides from the sample. The invention does not require separation of the population to "obtain" the sample, as later diagnosis may or may not need such separation.

The invention also provides a method of preparing a one-step aqueous formulation of the collection/lysis/transport/storage compositions described herein for the collection of nucleic acids such as RNA and/or DNA. In an overall sense, the method generally involves combining one or more chaotropes and nuclease-free water at a temperature of about 20.degree. C. to 90.degree. C. in a reaction zone; then combining the dissolved one or more chaotropes with one or more reducing agents, one or more chelators, and one or more detergents in the reaction zone to form an intermediate composition; optionally combining a silicone polymer with the intermediate composition in an amount sufficient to minimize foaming during further preparation of the one-step aqueous formulation; combining a sufficient amount of buffer to the intermediate composition to maintain a pH of about 6 to 6.9; optionally combining a second chelating agent to the reaction zone; then increasing the temperature of the second intermediate composition to about 60 to 95.degree. C. for about 1 to 30 minutes and lowering the temperature to ambient conditions; optionally then combining a C.sub.1-6alcohol with the contents of the reaction zone; and optionally adjusting the pH to be about 6.4 to 6.9.

In additional embodiments, the invention provides a method for preparing one-step aqueous formulations adapted to obtain a population of polynucleotides from a biological sample or specimen that is suspected of containing nucleic acids. This method generally involves at least the steps of: a) contacting the sample with an amount of the one-step aqueous formulation effective to:

- i) at least substantially kill or inactivate potentially-infectious pathogens in the sample;
- ii) at least lyse a portion of cells to release RNAs and/or DNAs from the sample; and

iii) at least substantially inhibit or prevent the released polynucleotides in the sample from further hydrolysis or enzymatic degradation, modification, or inactivation, so as to obtain the population of polynucleotides from the sample.

Preferably, the methods of the invention will include at least contacting the sample with an amount of one or more of the disclosed compositions at a temperature of from 0.degree. C. to about 40.degree. C. (more preferably at a temperature of 4.degree. C. to about 35.degree. C., and still more preferably at a temperature of 10.degree. C. to about 30.degree. C.) for a period of time of at least 24 hrs, more preferably, for a period of time of at least 48 hrs, at least 72 hrs, at least 96 hrs, or longer, without causing substantial deterioration, degradation, enzymatic cleavage, and/or nucleolytic digestion, modification, or processing of the nucleic acids contained within a sample contacted with such a composition.

In certain embodiments, the methods of the invention will include at least contacting the sample with an amount of one or more of the disclosed compositions at a temperature from about 0.degree. C. to about 40.degree. C. (more preferably at a temperature from about 4.degree. C. to about 35.degree. C., still more preferably at a temperature from about 10.degree. C. to about 30.degree. C., and more preferably still at a temperature from about 15.degree. C. to about 25.degree. C.) for a period of time of at least 7 days, more preferably, for a period of time of at least 14 days, at least 21 days, at least 28 days, or even longer without causing significant deterioration, degradation, enzymatic cleavage, and/or nucleolytic processing of the nucleic acids contained within a sample so processed. It should be understood that associating a sample with an inventive composition need only occur for a short time, but to avoid the need for immediate separation of the nucleic acids from the sample and the one-step composition of the invention all the materials may remain in contact for the time periods specified above without any substantial, or without any, degradation of the nucleic acids.

Preferably, the integrity of a population of polynucleotides released from the sample into the composition will be substantially maintained, even when the composition comprising the sample is stored at ambient temperatures, and even for prolonged periods of time, including, without limitation, storage for greater than about 10 days, greater than about 20 days, or even greater than about 30 days or more. Likewise, it is desirable that the integrity of a population of polynucleotides released from the sample into the composition will be substantially maintained, even when the composition comprising the sample is stored at subtropical and tropical temperatures--even for prolonged periods of time, including, without limitation, storage for greater than or equal to about 5 days, greater than or equal to about 10 days, greater than or equal to about 15 days, or even greater than or equal to about 20, about 25, about 30, about 35, about 40, about 60, or about 90 days or even greater.

In the practice of the present methods, it is preferable that at least one or more biological cells contained within the sample are substantially lysed to release at least a first population or first plurality of polynucleotides contained within such cells into the composition. Preferably, the components of the disclosed composition are sufficient to release such a population from all or substantially all of the remaining cellular/tissue and/or sample debris (including, without limitation, lipids, phospholipids, peptides, proteins, polysaccharides, lipopolysaccharides, polyols, cellular organelles, membrane components, and such like).

It is also desirable in the practice of the present methods that at least one or more exogenous or endogenous nucleases that may be present in, on, or about the sample itself, will be sufficiently inactivated by one or more components of the composition such that the resulting nucleic acids are not destroyed, damaged, or

nucleolytically cleaved when the biological cells contained within the sample are substantially lysed to release the population of polynucleotides from the cells. Preferably, one or more components of the disclosed composition are effective to kill, inactivate, or substantially inhibit the biological activity of a DNase or an RNase, when such a protein is present in the sample.

It is also desirable in the practice of the present methods that when one or more microbes, viruses, fungi, and/or other pathogens or organisms of interest are present in, on, or about the sample when collected, such microbes, viruses, fungi, and/or other pathogens or organisms of interest will be lysed, sufficiently inactivated, or substantially killed upon contact with the composition, which facilitates safe handling of the sample by the practitioner. Preferably, one or more components of the disclosed composition are effective to render a pathogenic sample substantially or preferably entirely, non-pathogenic without the need for adding additional components to the composition. However, in certain applications, it may also be desirable to include one or more additional anti-microbial, anti-viral, or anti-fungal agents to the compositions to render them substantially non-pathogenic, and thus, safe for handling by the practitioner.

Preferably, the composition containing the sample is at least sufficiently stable to permit storage of the sample in the composition at ambient, near-ambient, or even colder or warmer conditions at least substantially (or entirely) from the time of specimen or sample collection substantially until the time of analyzing or characterizing at least a first population of polynucleotides from within the sample. As used herein, "ambient temperature" can refer to temperatures of about 18.degree. C. to 25.degree. C., or in some embodiments, more preferably from about 20.degree. C. to about 22.degree. C.

In certain embodiments, the composition containing the sample may be stored at a temperature of about 0.degree. C. to about 40.degree. C., more preferably at a temperature of about 4.degree. C. to about 30.degree. C., more preferably, at a temperature of about 10.degree. C. to about 25.degree. C., at least substantially from the time of collection to the time that the polynucleotides obtained from the sample are further isolated, purified, or characterized using one or more conventional molecular biology methodologies.

In certain embodiments, the composition containing the sample suspected of containing nucleic acids will stabilize the nucleic acids to the extent that they either remain at least substantially non-degraded (i.e., at least substantially stable) even upon prolonged storage of the composition at ambient, refrigerator, or sub-zero temperatures. It will be desirable that this stability provides that at least about 70%, at least about 85%, more preferably at least about 90%, more preferably at least about 95%, or even more preferably, at least about 98% of the polynucleotides contained within the stored sample will not be degraded upon prolonged storage of the sample. In certain embodiments, substantially all of the polynucleotides contained within the sample will be stabilized such that the original integrity of the polynucleotides is preserved during the collection, lysis, storage, and transport of the processed sample.

In certain embodiments, the method will preferably provide a population of nucleic acids prepared from a biological sample in which less than about 15% of the polynucleotides contained in the sample will be degraded during the collection, lysis, storage, and transport of the sample after it has been stored in the composition at a temperature of from -20.degree. C. to about 40.degree. C. for a period of at least 24, 48, 72, or 96 hrs or longer after the sample was initially introduced into the composition.

In related embodiments, the method will preferably provide a population of nucleic acids prepared from a biological sample in which less than about 10% of the polynucleotides contained in the sample will be degraded during the collection, lysis, storage, and transport of the sample after it has been stored in the

composition at a temperature of from -20.degree. C. to about 40.degree. C. for a period of at least 24, 48, 72, or 96 hrs or longer after the sample was initially introduced into the composition.

Likewise, in some applications of the methodology disclosed herein, use of the disclosed compositions will preferably provide a population of nucleic acids that are prepared from a biological sample, wherein less than about 5% of the polynucleotides contained in the sample will be degraded during the collection, lysis, storage, and transport of the sample after it has been stored in the composition at a temperature from -20.degree. C. to about 40.degree. C. for a period of at least 24, 48, 72, or 96 hrs or longer after the sample was initially introduced into the composition.

In some instances, the population of nucleic acids prepared by the present methods may be maintained with sufficient integrity such that no more than about 1 or 2% of the sample will be degraded even when the composition is stored at a temperature from 0.degree. C. to about 40.degree. C. for periods of several days to several weeks. In fact, the inventors have shown that samples of nucleic acids isolated using the disclosed methods remain at least substantially stable, preferably stable, in their non-degraded form for periods of several weeks to even several months or more, even when the composition containing the nucleic acids is stored at a temperature from 10.degree. C. to about 40.degree. C. In one preferred embodiment, the upper limit on the above-noted temperature ranges is about 37.degree. C. Thus, the term "stable" as used herein may refer to the various embodiments noted above regarding the integrity of the population of nucleic acids after a particular time lapse at a given temperature.

The polynucleotides of the present invention, and particularly those useful in the preparation of IPCs, preferably contain at least a first polymerase transcription region operably linked to a positive control sequence (PCS). Such regions may be operatively positioned upstream of the PCS, such that the PCS is amplified when a suitable polymerase (e.g., T3, SP6, or T7 polymerase, and the like) is contacted with the PCS under suitable reaction conditions. The construct may further comprise a single primer, or alternatively, two or more primers (e.g., "forward" and "reverse" primers) may be used to facilitate expression of the PCS. Exemplary primers useful in the practice of the invention include, but are in no way limited to, those primer sequences that specifically bind to the PCS itself or to regions immediately upstream (5') and or downstream (3') of the PCS. In illustrative embodiments, the PCS will contain at least a first region to which a first detection probe (including, without limitation, luminescent, fluorescent, chemiluminescent, or FRET probes as described herein) specifically binds.

The polynucleotides useful in the preparation of IPCs may also further optionally comprise one or more native, synthetic, homologous, heterologous, or hybrid promoter(s), enhancer(s), regulatory element(s), linker(s), spacer(s), binding domain(s), or transcription activation site(s), etc.

In illustrative embodiments, the inventors have demonstrated that the disclosed sample collection/transport media may successfully be employed to permit the collection and storage of biological samples for periods of days, to weeks, to months, even when stored at ambient environmental conditions. As such, the disclosed formulations provide effective compositions for maintaining the fidelity and integrity of populations of isolated nucleic acids for about 1, about 2, about 3, about 4, about 5, about 6, about 7, about 8, about 9, or even about 10 weeks or more without significant deterioration in the quality of the polynucleotides contained within the sample. In fact, long-term storage of the samples is even possible under ambient and near-ambient conditions for periods of several months or more--a significant improvement over any conventional commercial formulations examined by the inventors during the present invention.

Data obtained using in vitro analyses of collected specimens have demonstrated that samples collected and stored in one of the disclosed collection/storage/transport formulations remain substantially non-degraded, and functionally-viable when stored in the solution for a period of a few days to a few weeks to a few months or longer. Comparison of results obtained after about 1, about 2, about 3, about 4, about 5, about 6, about 7, or even about 8 weeks or longer demonstrated the effectiveness of the solutions in preventing substantial nuclease activity against the isolated polynucleotides, and protecting the isolated molecules from contamination or molecular degradation.

Commercial Formulations and Kits

The present invention also provides kits and sample collection systems utilizing the disclosed compositions and collection/storage/transport solutions described herein. In particular embodiments, such sample collection systems may include a collection device, such as a swab, curette, or culture loop; and a collection vessel, such as a vial test tube, or specimen cup, that contains one or more of the compositions disclosed herein. The collection vessel is preferably releasably openable, such that it can be opened to insert the one-step compositions and closed and packaged, opened to insert the sample and optionally a portion of the collection device and closed for storage and transport, or both. The collection vessel may use any suitable releasable or openable mechanism, including without limitation a screw cap, snap top, press-and-turn top, or the like. Such systems may also further optionally include one or more additional reagents, storage devices, transport devices, and/or instructions for obtaining, collecting, lysing, storing, or transporting samples in such systems. In a preferred embodiment, the one-step compositions of the invention may already be disposed in the reaction zone into which the sample may be associated. In such embodiments, the invention requires only a collection device and the collection vessel. The kit may also include one or more isolation or extraction devices to help liberate and/or separate one or more populations or pluralities of nucleic acids contained within the sample from one or more other biomolecules or sample components to obtain at least partially, or substantially purified nucleic acids suitable for identification, detection, or further molecular analysis.

Kits may also be packaged for commercial distribution, and may further optionally include one or more collection, delivery, transportation, or storage devices for sample or specimen collection, handling, or processing. The container(s) for such kits may typically include at least one vial, test tube, flask, bottle, cup, or other suitable container or collection device, into which the composition(s) may be placed, and, preferably, suitably aliquoted for individual specimen collection, transport, and/or storage. The kit may also include a larger container, such as a case, that includes the smaller, individual containers noted above, along with other collection devices, equipment, reagents, instructions, and/or the like. The kit may also optionally include one or more additional buffers, compounds, or compositions, and may also further optionally include one or more instructions detailing use(s) of the kit in either the collection or storage of one or more biological, clinical, diagnostic, environmental, or forensic sample. Optionally, the kit may also further provide instructions for the transport of the sample once placed in one or more of the disclosed compositions, and may even include instructions or additional reagents detailing one or more subsequent analytical methods or assays employing the nucleic acids isolated from the sample or specimen. Such kits may also include multiples of the various collection devices and collection vessels and any other components to be included, so that the kits can be used to collect multiple samples from the same source or different sources. In one commercial application, the kits are packaged in sets of five or more for convenient sale and use.

BRIEF DESCRIPTION OF THE DRAWINGS

For promoting an understanding of the principles of the invention, reference will now be made to the

embodiments, or examples, illustrated in the drawings and specific language will be used to describe the same. It will, nevertheless be understood that no limitation of the scope of the invention is thereby intended. Any alterations and further modifications in the described embodiments, and any further applications of the principles of the invention as described herein are contemplated as would normally occur to one of ordinary skill in the art to which the invention relates.

The following drawings form part of the present specification and are included to demonstrate certain aspects of the present invention. The invention may be better understood by reference to the following description taken in conjunction with the accompanying drawings, in which like reference numerals identify like elements, and in which:

FIG. 1 shows the extraction efficiency of PSS (ver. 1). PSS (ver. 1 [depicted here as "One-step +"]) compared to the lysis solution provided in the RNAqueous.RTM.-Micro Kit (Ambion, Catalog #AM1931) using a standard amount of whole influenza A virus. For the comparison, either the one-step formulation or the lysis solution provided in the RNAqueous-Micro Kit was used for viral RNA lyses and then extracted according to manufacturer protocols. Replicate reactions were processed and analyzed by qRT-PCR using, for example, the Applied Biosystems ABI7500 sequence detection system (Life Technologies Corporation, Carlsbad, Calif., USA);

FIG. 2 shows the extraction efficiency of PSS (ver. 1) compared to commercial kits. Homogenized cotton rat nose(*) challenged with influenza A (H3N2) or a human clinical influenza A (H1N1) samples collected during the 2006-07 season were lysed in the PSS or lysed using the respective lyses solution, protocol, and extraction procedure from three commercially available kits: RNAqueous-Micro Kit (Ambion), QiaAmp.TM. Viral Mini Kit (Qiagen), and AI/NCD MaxMag.TM. (Ambion) Kit. Extraction efficiency was evaluated using the Applied Biosystems ABI7500 platform with the comparative C.sub.T method. The relative C.sub.T scores and viral copies detected were optimal when PSS (depicted as the "one-step formulation") was utilized in place of the respective lysis buffer provided in each standard commercially-available kit;

FIG. 3 shows the preservation of naked RNA in PSS vs. Ambion's RNA Storage Solution. Single-stranded Avian H5 RNA was stored in PSS, RNA storage solution (Ambion), or water at ambient temperature (22-24.degree. C.) for 96 hours. A total of 5 pg of RNA was extracted using the RNAqueous.RTM.-Micro Kit (Ambion, Cat#AM1931) according to the manufacturer's recommendations, and subsequently analyzed using qRT-PCR on the Applied Biosystems ABI7500 platform (Life Technologies). Values are given as C.sub.T using the absolute quantification method;

FIG. 4 shows an example of a PSS packaging format for clinical diagnostic collection. Directions of sample collection using a clinical collection swab (Copan Diagnostics) and a 5 mL collection tube containing 1.5 mL of PSS;

FIG. 5 shows an exemplary commercial formulation of the PSS. Three exemplary commercial collection solution formats are depicted: a 25 mL bottle, and the 5 mL and 1.5 mL test tubes;

FIG. 6 illustrates the ability of PSS to rapidly kill microorganisms. Shown is a comparison of cell growth of methicillin-resistant Staphylococcus aureus (MRSA) either in culture medium (TSB), or in a solution of PS. After 10 seconds in PSS, no viable bacterial pathogens were detected.

FIG. 7A shows the inactivation of chicken cloacal specimens in PSS (ver. 1). PSS (ver. 1) inactivated

microbial agents in .ltoreq.1 hr. Four original chicken cloacal samples were immersed in PSS (ver. 1) (top row) or water (bottom row) and subsequently plated on blood agar plates;

FIG. 7B demonstrates that PSS inhibits RNA base hydrolysis for 30 days at room temperature. RNA was incubated at room temperature (22-26.degree. C.) in PSS (gel lane 1 and 3) and water (gel lane 2 and 4), and subsequently RT-PCR amplified (1500 base pairs) at Day 0 and Day 30. PSS preserved collected RNA, and prevented RNA/DNA degradation at room temperature up to 30 days;

FIG. 8A depicts the qRT-PCR analysis of "naked" influenza A avian H5 RNA template preserved in PSS after incubation in RNA/DNA nucleases. H5 cRNA (2 ng) was incubated with ribonuclease A and T1, and DNaseI for 1 hour @ 37.degree. C. and extracted using the RNAaqueous.RTM.-Micro Kit (Ambion). Triplicate reactions were included for each reaction condition. qRT-PCR C.sub.T values of naked RNA preserved in PSS with added nucleases (average C.sub.T: 22.88) were similar to an equal quantity of template cRNA control (average C.sub.T: 23.70). Template cRNA reactions subjected to nuclease digestion without PSS were almost completely degraded (average C.sub.T 39.58);

FIG. 8B depicts the qRT-PCR analysis and gel electrophoresis of "naked" influenza A avian H5 RNA template preserved in PSS after incubation in RNA/DNA nucleases at 37.degree. C. for 7 days. Two nanograms of H5 cRNA was incubated with RNase A, RNase T1, and DNase I, then extracted using the RNAaqueous.RTM.-Micro Kit (Ambion) after 7 days. Duplicate reactions were included for each reaction. qRT-PCR C.sub.T values of naked RNA preserved in PSS with added nucleases (average C.sub.T: 33.51) were detected after 7 days. Template cRNA reactions subjected to nuclease digestion without PSS were completely degraded, and similar to NTC reactions.

FIG. 8C demonstrates that PSS formulations substantially inhibited digestion of the sample polynucleotides by endogenous nucleases. Gel electrophoresis of post-amplified product. Lane 3 is the PCR product from template RNA+PSS at 37.degree. C. for 7 days, and Lane 5 amplification of positive control RNA. Lane 5 (no amplification) contained RNA without PSS; Lanes 2 and 6 are 100-bp ladder, and NTC reactions, respectively.

FIG. 9 illustrates that the preservation of nucleic acids achieved in PSS is superior to that of other solutions. PSS (ver. 2 and ver. 2.2) Preservation of RNA from influenza A virus compared to Qiagen AVL buffer, ethanol, and viral transport medium (VTM) (Becton-Dickinson, Franklin Lakes, N.J., USA) at ambient temperature (22-25.degree. C.) for 30 days.

FIG. 10 shows the extraction efficiency of Influenza A virus preserved in PSS (ver. 2.2) for 30 days at various temperatures. Environmental (21-37.degree. C.); Freeze-thawed (-25.degree. C.; 32.times.); ambient temperature (22-26.degree. C.); and Lane 5: refrigerated (4.degree. C.).

FIG. 11 is a graph of C.sub.T vs. molar concentration using whole influenza A virus with TCEP as the reducing agent;

FIG. 12 is a graph of C.sub.T vs. molar concentration using H5N1 Avian influenza ssRNA with TCEP as the reducing agent;

FIG. 13A and FIG. 13B show the comparison between TCEP and .beta.-ME as reducing agent components in the PSS formulations, with water only serving as a control. In FIG. 13A, H5 avian influenza RNA was

employed, while in FIG. 13B, whole virus were used.

FIG. 14A shows the results of a study employing PSS in preserving nucleic acids from blood. PSS extraction efficiency of whole blood spiked with RNA compared to the lysis solution in the QIAamp.RTM. DNA Blood Mini Kit. 0.1 pg and 1 pg of influenza A RNA were spiked and extracted using PSS or AL Lysis buffer. At both RNA concentrations, PSS produced superior results as evident by qRT-PCR C.sub.T scores;

FIG. 14B tabulates data from the study shown in FIG. 14A involving the extraction of "naked" H5 avian influenza ssRNA from blood tubes. PSS extraction efficiency of whole blood spiked with RNA compared to the lysis solution in the QIAamp.RTM. DNA Blood Mini Kit using different blood anticoagulants. PSS was superior compared to the AL Lysis Buffer from Qiagen using Blood spiked with RNA in common anticoagulant blood-collection tubes; and

FIG. 14C tabulates data from the study shown in FIG. 14A involving the comparison of in PSS vs. a commercial extraction kit storage reagent (QIAamp.RTM., Qiagen, Inc.). Shown is the extraction efficiency of PSS for whole blood spiked with RNA compared to the lysis solution in the QIAamp.RTM. DNA Blood Mini Kit. 0.1 pg and 1 pg of Influenza A viral RNA were spiked and extracted using PSS or AVL Lysis buffer. At both RNA concentrations, PSS produced superior results as evident by qRT-PCR C.sub.T scores.

FIG. 15A and FIG. 15B show the qPCR amplification of 141 nucleotide single stranded DNA using an assay specific for detecting the IPC. FIG. 15A: Amplification was performed in duplicate reactions using 10.sup.-6 to 10.sup.-3 (i.e., 8-logs) of serially-diluted DNA template. FIG. 15B: IPC standard curve of DNA amplification. The assay slope was -3.39 and the calculated PCR efficiency was 97.23% using $E=10^{\text{sup.}(-1/\text{slope})-1}$. qPCR was performed using the Applied Biosystems ABI7500 platform.

FIG. 16A and FIG. 16B show the qRT-PCR amplification of a 130-nucleotide ssRNA carrier molecule capable of serving as an IPC. In FIG. 16A, shown: Amplification was performed in duplicate reactions using 5 picograms (pg) (2.75.times.10.sup.7 copies) to 5 attograms (ag) (.about.2-10 copies) of serially-diluted RNA template. In FIG. 16B, shown is the IPC standard curve for amplification of the carrier ssRNA molecule. The assay slope was -3.144 and the calculated PCR efficiency was 100% using $E=10^{\text{sup.}(-1/\text{slope})-1}$. qRT-PCR was performed using the Applied Biosystems ABI7500 platform.

FIG. 17 shows PSS containing synthetic ssRNA carrier/IPC molecules enhanced the preservation of naked RNA pathogen target. In this study, a naked influenza A RNA template (1.27 ng/.mu.l) was added to a human nasal wash specimen and preserved PSS with (closed squares) or without (closed circles) synthetic RNA augmentation at 37.degree. C. RNA was extracted, and then subsequently detected via qRT-PCR. The presence of the synthetic ssRNA carrier/IPC molecules in the formulation resulted in both an enhanced initial extraction (on day 0) of sample polynucleotides, and also a stable preservation of the isolated polynucleotides over extended periods of time (see e.g., day 31).

FIG. 18 illustrates that the ability of PSS to stabilize and preserve sample polynucleotides in solution is enhanced by the addition of synthetic ssRNA carrier/IPC nucleic acids to the PSS formulation. In this study, whole influenza A virus was preserved in PSS with (closed squares) or without (closed circles) synthetic ssRNA carrier/IPC augmentation at 37.degree. C. for 35 days. Whole flu virus was also preserved into VTM (closed triangles). RNA was extracted and detected using qRT-PCR. The presence of the synthetic carrier RNA molecule in the formulation resulted in both an enhanced initial extraction (on day 0) of sample polynucleotides, and also a stable preservation of the isolated polynucleotides over extended periods of time

(see e.g., day 35).

FIG. 19 illustrates the stability of the synthetic ssRNA carrier molecule at room temperature and 37.degree. C. In this study, synthetic RNA (0.1 pg/.mu.l) was added to PSS and stored at room temperature (RT) or 37.degree. C. for 14 days. Stability of the RNA was measured as a function of qRT-PCR using the IPC assay using the Applied Biosystems ABI7500 platform. Synthetic RNA was stable at RT and 37.degree. C. for 14 days with no significant loss in C.sub.T values. The average of duplicate reactions is shown.

DETAILED DESCRIPTION OF ILLUSTRATIVE EMBODIMENTS

Illustrative embodiments of the invention are described below. In the interest of clarity, not all features of an actual implementation are described in this specification. It will of course be appreciated that in the development of any such actual embodiment, numerous implementation-specific decisions must be made to achieve the developers' specific goals, such as compliance with system-related and business-related constraints, which will vary from one implementation to another. Likewise, it will also be appreciated that while such development efforts may be either complex or time-consuming, they would nevertheless be a routine undertaking for one of ordinary skill in the art having the benefit of the present disclosure.

The extended stabilization, collection, transport, and preservation imparted by the disclosed formulations are particularly advantageous when a sample or specimen is located in a geographical region that is remote from a testing facility. Remote locations, also referred to as field sites, encompass a variety of environments where diagnostic testing is typically not performed. These sites include doctors' offices, triage centers, airports, border crossings, outbreak areas, and a variety of outdoor locations. The disclosed compositions and methods for their use offer particular advantages in locations where there is no access to electricity and/or refrigeration, or where access is inconsistent. Because of the extended stability at room temperature, a sample can be taken from any remote location, for example without limitation at a malarial outbreak site in Africa, and the sample can be shipped to the United States or Europe for diagnostic analysis in a laboratory. Because the disclosed collection formulations are stable at room temperature or below, and preferably even at tropical or subtropical temperatures for a time, they can routinely be taken into the field without worry that the component reagents (such as RNA controls) themselves will degrade until a sample can be analyzed, typically at a remote location from the collection.

The compositions of the invention may be any suitable aqueous formulation as described herein, including but not limited to a solution, suspension (incl. colloidal suspension), slurry, emulsion, homogenate, or the like. A preferred aqueous formulation is a solution, and therefore the term "solution" has been used in the exemplary sense throughout the detailed description of the preferred embodiments to refer to any of the aqueous compositions of the invention.

Important aspects of the present invention concern isolated RNA and DNA segments and recombinant vectors and populations of polynucleotides comprising them, and the creation and use of recombinant methodologies to prepare RNA carrier molecules through the application of conventional molecular biology technologies, including, without limitation, qPCR, qRT-PCR and the like.

The present invention also concerns nucleic acid compositions, including, without limitation, DNA, RNA and PNA, isolatable from one or more biological samples or specimens using the specimen collection/storage/transport media described herein, or the isolated RNA or DNA segments utilized as carrier molecules and/or IPCs.

As used herein, the term "DNA segment" refers to a DNA molecule that has been isolated free of total genomic DNA of a particular species. Therefore, a DNA segment obtained from a biological sample using one of the compositions disclosed herein, refers to one or more DNA segments that have been isolated away from, or purified free from, total mammalian or human genomic DNA. Included within the term "DNA segment", are DNA segments and smaller fragments of such segments, as well as recombinant vectors, including, for example, plasmids, cosmids, phage, viruses, and the like.

Similarly, the term "RNA segment" refers to an RNA molecule that has been isolated free of total cellular RNA of a particular species. Therefore, RNA segments obtained from a biological sample using one of the compositions disclosed herein, refers to one or more RNA segments (either of native or synthetic origin) that have been isolated away from, or purified free from, other RNAs. Included within the term "RNA segment", are RNA segments and smaller fragments of such segments.

A nucleic acid segment comprising an isolated and/or purified polynucleotide refers to a nucleic acid segment isolated substantially away from other naturally occurring genes or protein-encoding sequences that may be present in the biological sample undergoing collection and subsequent analysis.

In particular embodiments, the invention concerns nucleic acid sample collection, storage, and transport media that comprise a first carrier molecule to aid in the isolation and stabilization of nucleic acids obtained from such samples. Preferably, the carrier molecules are single-stranded RNA molecules, although in some circumstances, RNA:RNA, or RNA:DNA duplex molecules may also be utilized.

While the particular molecular mechanism for how the incorporation of exogenous carrier molecules, and RNA molecules in particular, into the specimen collection/storage/transport formulations disclosed herein need not be elucidated for the routine practice of the invention, without being limited by any particular theory, the inventors believe that the presence of exogenously-supplied RNA carrier molecules in PSS and related formulations described herein aids in the precipitation and stabilization of polynucleotides that are liberated from the cellular material when contacted with PSS, and may also serve as a method for "coating" or "enrobing" the liberated sample nucleic acids in a protective molecular form that is favorable to polynucleotide agglomeration or precipitation.

Regardless of the mechanism of action, the inventors have demonstrated the surprising benefits afforded over conventional polynucleotide isolation solutions in the preparation of high-quality, high-fidelity populations of nucleic acids that are spared the deleterious effects routinely observed in conventional samples that are subjected to degradation and contamination during the collection, isolation, storage, and transport processes.

Carrier molecules and/or IPC-containing sequences may be synthesized in part or in whole, using conventional methodologies known in the art, or may be prepared from recombinant molecular biological methods (including, without limitation, in vitro transcription and/or translation methods, such as PCR amplification), or may be obtained from native biological sources (preferably from non-mammalian sources and sources which are not pathogenic to mammalian species). Alternatively, the carrier RNA molecules of the present invention may be hybrid molecules-containing both natural and synthetic portions operably linked using one or more conventional methods in the molecular genetic and recombinant DNA/RNA arts. In certain applications the RNA carrier molecules may be comprised of one or more transcription initiation sequences, primer or probe binding domains, functionalized binding sites, exo- or endo-nuclease cleavage sites, and/or functional motifs or active sites.

The term "a sequence essentially as set forth in SEQ ID NO:X" means that the sequence substantially corresponds to a portion of SEQ ID NO:X and has relatively few nucleotides (or amino acids in the case of polypeptide sequences) that are not identical to, or a biologically functional equivalent of, the nucleotides of SEQ ID NO:X. The term "biologically functional equivalent" is well understood in the art, and is further defined in detail herein. Accordingly, sequences that have about 85% to about 90%; or more preferably, about 91% to about 95%; or even more preferably, about 96% to about 99%; of nucleotides that are identical or functionally equivalent to one or more of the nucleotide sequences provided herein are particularly contemplated to be useful in the practice of the invention.

Suitable standard hybridization conditions for the present invention include, for example, hybridization in 50% formamide, 5.times.Denhardt's solution, 5.times.SSC, 25 mM sodium phosphate, 0.1% SDS and 100 .mu.g/ml of denatured salmon sperm DNA at 42.degree. C. for 16 h followed by 1 hr sequential washes with 0.1.times.SSC, 0.1% SDS solution at 60.degree. C. to remove the desired amount of background signal. Lower stringency hybridization conditions for the present invention include, for example, hybridization in 35% formamide, 5.times.Denhardt's solution, 5.times.SSC, 25 mM sodium phosphate, 0.1% SDS and 100 .mu.g/ml denatured salmon sperm DNA or E. coli DNA at 42.degree. C. for 16 h followed by sequential washes with 0.8.times.SSC, 0.1% SDS at 55.degree. C. Those of skill in the art will recognize that conditions can be readily adjusted to obtain the desired level of stringency.

Naturally, the present invention also encompasses DNA segments that are complementary, or essentially complementary, to one or more of the specific sequences set forth herein. Nucleic acid sequences that are "complementary" are those that are capable of base-pairing according to the standard Watson-Crick complementarity rules. As used herein, the term "complementary sequences" means nucleic acid sequences that are substantially complementary, as may be assessed by the same nucleotide comparison set forth above, or as defined as being capable of hybridizing to one or more of the specific nucleic acid segments disclosed herein under relatively stringent conditions such as those described immediately above.

As described above, the probes and primers of the present invention may be of any length. By assigning numeric values to a sequence, for example, the first residue is 1, the second residue is 2, etc., an algorithm defining all primers can be proposed: n to n+y

where n is an integer from 1 to the last number of the sequence and y is the length of the primer minus one, where n+y does not exceed the last number of the sequence. Thus, for a 25-mer, the probes correspond to bases 1 to 25, 2 to 26, 3 to 27 . . . and so on. For a 45-mer, the probes correspond to bases 1 to 45, 2 to 46, 3 to 47 . . . and so on. For a 60-mer, the probes correspond to bases 1 to 60, 2 to 61, 3 to 62 . . . and so on.

In certain embodiments, it will be advantageous to employ nucleic acid sequences of the present invention in combination with an appropriate detectable marker (i.e., a "label,") for determining hybridization. A wide variety of appropriate indicator compounds and compositions are known in the art, including, without limitation, fluorescent, radioactive, enzymatic or other ligands, such as avidin/biotin, etc., which are capable of being detected. In particular embodiments, one may desire to employ a fluorescent label or an enzyme tag such as urease, alkaline phosphatase or peroxidase, instead of radioactive or other environmentally less-desirable reagents. In the case of enzyme tags, colorimetric, chromogenic, or fluorogenic indicator substrates are known that can be employed to provide a method for detecting the sample that is visible to the human eye, or by analytical methods such as scintigraphy, fluorimetry, spectrophotometry, and the like, to identify specific hybridization with samples containing one or more complementary or substantially complementary

nucleic acid sequences.

In general, it is envisioned that the hybridization probes described herein will be useful both as reagents in solution hybridization, as in PCR, for detection of particular nucleotides, as well as in embodiments employing a solid phase.

Once hybridized, the nucleic acid:primer complex is contacted with one or more enzymes that facilitate template-dependent nucleic acid synthesis. Multiple rounds of amplification, also referred to as "cycles," are conducted until a sufficient amount of amplification product is produced.

Next, the amplification product is detected. In certain applications, the detection may be performed by visual inspection. Alternatively, the detection may involve indirect identification of the product via chemiluminescence, radioactive scintigraphy of incorporated radiolabel or fluorescent label or even via a system using electrical or thermal impulse signals.

A number of template dependent processes are available to amplify the marker sequences present in a given template sample. One of the best-known amplification methods is the polymerase chain reaction (referred to as PCR) which is described in detail in U.S. Pat. Nos. 4,683,195, 4,683,202 and 4,800,159, and each incorporated herein by reference in entirety.

Another method for amplification is the ligase chain reaction ("LCR"), disclosed, e.g., in EPA No. 320 308, and U.S. Pat. No. 4,883,750, each of which is incorporated herein in its entirety by express reference thereto.

An isothermal amplification method, in which restriction endonucleases and ligases are used to achieve the amplification of target molecules that contain nucleotide 5'-[.alpha.-thio]-triphosphates in one strand of a restriction site may also be useful in the amplification of nucleic acids in the present invention.

Strand Displacement Amplification (SDA) is another method of carrying out isothermal amplification of nucleic acids that involves multiple rounds of strand displacement and synthesis, i.e., nick translation. A similar method, called Repair Chain Reaction (RCR), involves annealing several probes throughout a region targeted for amplification, followed by a repair reaction in which only two of the four bases are present. The other two bases can be added as biotinylated derivatives for easy detection. A similar approach is used in SDA. Target specific sequences can also be detected using a cyclic probe reaction (CPR). In CPR, a probe having 3' and 5' sequences of non-specific DNA and a middle sequence of specific RNA is hybridized to DNA that is present in a sample. Upon hybridization, the reaction is treated with RNase H, and the products of the probe identified as distinctive products that are released after digestion. The original template is annealed to another cycling probe and the reaction is repeated.

Methods based on ligation of two (or more) oligonucleotides in the presence of nucleic acid having the sequence of the resulting "di-oligonucleotide," thereby amplifying the di-oligonucleotide, may also be used in the amplification step of the present invention.

Following any amplification, it may be desirable to separate the amplification product from the template and the excess primer for determining whether specific amplification has occurred. In one embodiment, amplification products are separated by agarose, agarose-acrylamide or polyacrylamide gel electrophoresis using conventional methods that are known to those of ordinary skill in the art.

Polynucleotide and Oligonucleotide Compositions

As used herein, "nucleic acid" or "polynucleotide" compositions include, but are not limited to, those that contain either single-stranded or double-stranded polynucleotides, such as for example, deoxyribonucleic acids (DNA), ribonucleic acids (RNA), peptide nucleic acids (PNA), or any combinations or derivatives thereof (including, e.g., genomic, extragenomic, plasmid, cosmid, recombinant, artificial, and/or synthetic). Such sequences may be coding or non-coding sequences, sense, non-sense, or anti-sense sequences, and may, but need not, be present within one or more populations or pluralities of polynucleotides (either of the present invention, and a polynucleotide may, but need not, be linked to other molecules and/or support materials).

Likewise, polynucleotides of the present invention, and particularly those functioning as carrier polynucleotides or IPC molecules in the disclosed specimen collection/transportation media, need not be identical, or even substantially homologous to the particular sequences employed in the various embodiments of the invention illustrated herein. While the inventors have illustrated the use of particular DNA and RNA control sequences as tools for monitoring the integrity and quality, and preferably, increasing the stability, of polynucleotides collected, transported, or stored in the compositions described herein, such control sequences need not contain the particular nucleotide sequences employed.

In fact, in certain circumstances, polynucleotides useful as internal control sequences in the disclosed collection/transport/storage compositions may comprise any suitable sequence that may be obtained, prepared, modified, or synthesized for such purpose. Moreover, it is preferable that such internal control sequences do not share significant homology or identity with the viral, bacterial, fungal, or other pathogenic species sought to be identified using the disclosed formulations.

Alternatively, polynucleotides suitable as IPCs may be homologous to, or even substantially identical to, one or more of the exemplary IPC sequences disclosed herein, and as such, may be capable of hybridizing to one of the specific IPC sequences disclosed herein (or a sequence that is complementary thereto) under moderately stringent, highly stringent, and/or even very highly stringent hybridization conditions.

The methods for nucleic acid hybridization are considered routine to those of ordinary skill in the molecular biological arts, and as such, a detailed discussion of analytical methods employing them need not be provided herein. However, as a guidance, "moderately stringent" hybridization conditions popularized by Southern et al. are generally considered in the art to include, e.g., pre-washing in a solution containing about 5.times. standard sodium citrate buffer (SSC), 0.5% sodium dodecyl sulfate (SDS), 1.0 mM ethylenediaminetetraacetic acid (EDTA) (e.g., pH 8.0); hybridizing at a temperature of from about 50.degree. C. to about 60.degree. C. in 5.times.SSC overnight; followed by washing twice at about 60 to 65.degree. C. for 20 min. with each of 2.times., 0.5.times. and 0.2.times.SSC containing 0.1% SDS). Likewise, "stringent" hybridization conditions typically include, e.g., pre-washing in a solution containing about 5.times.SSC, 0.5% SDS, 1.0 mM EDTA (pH 8.0); hybridizing at a temperature of from about 60.degree. C. to about 70.degree. C. in 5.times.SSC overnight; followed by washing twice at about 65 to 70.degree. C. for 20 min. with each of 2.times., 0.5.times. and 0.2.times.SSC containing 0.1% SDS). Similarly, representative examples of "highly-stringent" hybridization conditions include, but are not limited to, pre-washing in a solution containing about 5.times.SSC, 0.5% SDS, 1.0 mM EDTA (pH 8.0); hybridizing at a temperature of from about 70.degree. C. to about 75.degree. C. in 5.times.SSC overnight; followed by washing twice at about 70.degree. C. to about 75.degree. C. for 20 min. with each of 2.times., 0.5.times. and 0.2.times.SSC containing 0.1% SDS).

It will also be appreciated by those of ordinary skill in the art that, as a result of the degeneracy of the genetic

code, there are many nucleotide sequences that encode a given primary amino acid sequence. Some of these polynucleotides bear minimal homology to the nucleotide sequence of any native gene. Nonetheless, polynucleotides that vary due to differences in codon usage are specifically contemplated by the present invention.

IPCs may be prepared by conventional molecular biology recombination methodologies, or alternatively synthesized in whole or in part by conventional methods known in the art, including chemical synthesis (e.g., solid phase phosphoramidite chemical synthesis) and the like. Modifications in a polynucleotide sequence may also be introduced using standard mutagenesis techniques, such as oligonucleotide-directed site-specific mutagenesis. RNA molecules for IPCs may also be directly synthesized, or alternatively, be prepared by in vitro or in vivo transcription of DNA sequences using suitable systems (such as T3, T7, and SP6 polymerases and the like).

Polynucleotides of the present invention, and particularly those used to prepare IPCs, may be modified to increase stability either in vitro and/or in vivo. Such modifications include, without limitation, the addition of flanking sequences at the 5'-end, 3'-end, or both; the use of phosphorothioate or 2'-o-methyl rather than phosphodiesterase linkages in the backbone; and/or the inclusion of nontraditional bases such as inosine, queosine and wybutosine, as well as acetyl-, methyl-, thio-, or otherwise-modified forms of adenine, cytidine, guanine, thymine and uridine, or any combination thereof.

Nucleotide sequences as described herein may be joined or linked to a variety of other nucleotide sequences using established recombinant techniques. For example, a polynucleotide useful as an IPC may be produced by cloning into any of a variety of cloning vectors, including one or more of plasmids, phagemids, lambda phage derivatives and cosmids. Vectors of particular interest include expression vectors, replication vectors, probe generation vectors, and sequencing vectors. In general, a vector will contain an origin of replication functional in at least one organism, convenient restriction endonuclease sites and one or more selectable markers. Other elements will depend upon the desired use, and will be apparent to those of ordinary skill in the art. Alternatively, IPC sequences may be prepared through one or more template-dependent, or amplicon-directed recombinant production methodologies.

In particular embodiments, the present invention provides polynucleotide compositions that may be added to the disclosed collection/storage/transport media to provide one or more IPCs to monitor the fidelity and integrity of the sampling/storage/and transporting processes employed herein. Such polynucleotide compositions may contain one or more sequence domains to which specific polymerases, oligonucleotides primers, and/or recombinant probes may bind. In illustrative embodiments, exemplary IPCs may be employed, which contain a positive control sequence (PCS) that may be amplified using one or more specific amplification primer sequences, and detected and/or quantitated using one or more oligonucleotides probes that specifically hybridize to at least a first sequence domain of the PCS. In certain embodiments, detection probes may be utilized that employ fluorescent, chemiluminescent, or FRET-specific detection methodologies to detect and/or quantitate the amount of PCS in a given sample. Detection of the PCS in a sample can provide an internal reference standard for monitoring the effectiveness and fidelity of the one-step collection/storage/transport solutions employed herein.

Oligonucleotide primers and probes of the present invention may be designed for the selective amplification and detection of one or more PCS. Such primer sequences are suitable for use in hybridization methods, and in amplification methods such as PCR-based amplification methods (including, for example, real-time PCR analyses, RT-PCR and the like). Likewise, the disclosed oligonucleotide detection probes are suitable for

labeling with an appropriate label for detection and quantitation of the products resulting from the amplification of nucleic acids using one or more pairs of the amplification primers disclosed herein.

When labeled with appropriate markers, oligonucleotide detection probes are particularly suited for fluorescence-based detection of the IPC-specific nucleotide sequences contained in the collection/storage medium (including, without limitation, PSS formulations), including, for example, via FRET-based analyses. FRET-labeled detection probes are particularly useful in fluorimetric detection methodologies, including for example, the FRET-based microvolume fluorimetry devices. Use of one or more amplification and detection oligonucleotides is particularly contemplated in the combined real-time PCR/microvolume fluorimetry FRET-based methodologies (real-time PCR-FRET), and particularly in analyses facilitated by the "LightCycler.RTM." instrumentation as developed by Idaho Technology, Inc. (Salt Lake City, Utah, USA), and now manufactured and marketed by Roche Applied Science (Indianapolis, Ind., USA).

Internal Positive Control (IPC) Sequences

In illustrative embodiments, the invention also provides IPC sequences that comprise, consist essentially of, or consists of, nucleic acid sequences that are preferably at least about 80%, at least about 85%, at least about 90%, at least about 95%, at least about 96%, at least about 97%, at least about 98%, or at least about 99% or more identical to any one or more of the sequences disclosed above, and particularly those disclosed in any one of SEQ ID NO:2 through SEQ ID NO:13.

In related embodiments, the invention also provides IPC sequences that are synthesized in vitro using a template amplicon that comprises, consists essentially of, or consists of, nucleic acid sequences that are preferably at least about 80%, at least about 85%, at least about 90%, at least about 95%, at least about 96%, at least about 97%, at least about 98%, or at least about 99% or more identical to an at least 30 contiguous nucleotide sequence from any one or more of the sequences disclosed in SEQ ID NO:1, or SEQ ID NO:14 through SEQ ID NO:23.

Polynucleotide Amplification Kits

The present invention also provides kits for detecting and/or amplifying IPCs, and in particular RNA or DNA IPC sequences contained within one of the formulations disclosed herein, including without limitation, PSS formulations. Such kits typically comprise two or more components necessary for amplifying IPC-specific sequences, and such components may be compounds, reagents, containers and/or equipment. For example, one container within a kit may contain a first primer, while a second container within the kit may comprise a second primer. A third container within the kit may contain a set of hybridization probes, or one or more fluorescent probes for labeling the probes. In addition, the kits of the invention may also comprise instructions for use, e.g., instructions for using the primers in amplification and/or detection reactions as described herein, as well as one or more fluorescent molecules or other reagents as may be necessary, including for example, but not limited to, buffers enzymes, polymerases, RNases and such like.

The invention provides one or more carrier RNA or DNA compositions together with one or more collection/transport/storage solutions (including, without limitation, PSS formulations), or compatible excipients, buffers carriers, diluents, adjuvants, and/or other components, as may be employed in the isolation, detection, amplification, or quantitation of diagnostic tools.

Real-Time PCR-Based Fret Detection

Real-time PCR and FRET methodologies have been well described in the literature (see, for example, U.S. Pat. Nos. 4,996,143, 5,565,322, 5,849,489, 6,162,603, each of which is specifically incorporated herein in its entirety by express reference thereto). The LightCycler.RTM. platform represents a significant breakthrough in genetic mutation screening and analysis. This system incorporates a rapid, air-driven thermal cycling instrument that can perform 30 polymerase chain reaction (PCR) cycles in less than 20 minutes. It utilizes an in-line microvolume fluorimeter to detect and quantitate fluorescently-labeled hybridization probes, and provides the data necessary for determination of melting curve analyses. The LightCycler.RTM. platform provides innovative instrumentation to facilitate the development of genetic analysis tools, and to provide a rapid, qualitative method for the assay of specific nucleotide sequences, and genetic mutations. Detailed application of the instrumentation in amplification and detection methods may be found on the manufacturer's website, and in product application manuals. This technology has also been described, including for example PCT Intl. Appl. Publ. Nos. WO 97/46707, WO 97/46714 and WO 97/46712 (each of which is specifically incorporated in its entirety by express reference thereto).

Specimen Collection for Clinical Diagnostic Laboratories

Collection is first step in diagnostic platforms or molecular protocols requiring the detection of potentially minute amounts of nucleic acids from pathogens including viruses. To facilitate the dynamic advancements in nucleic acid based detection strategies and their integration into the mainstream diagnostic laboratories there is a colossal need for reliable, robust, and standardized collection systems developed specifically with the intent of being utilized for downstream nucleic acid based detection such as the aforementioned platforms. The invention can alternatively be adapted for transport of nucleic acids from a doctor's office or operating room, or alternatively transported to a regional center, such as a hospital.

A clinical or veterinary specimen or a forensic or environmental sample collection system may include one or more collection tools and one or more reagents for efficiently:

- 1) obtaining a high yield of suitable specimen beyond what is currently available in the art;
- 2) inactivating potentially infectious biological pathogens so that they are no longer viable and can be handled; shipped, or transported with minimal fear of pathogen release or contamination; or
- 3) effectively stabilizing and preserving lysed `naked` RNA/DNA polymers from hydrolysis or nuclease degradation for prolonged periods at ambient temperatures until samples can be processed at a diagnostic laboratory, and preferably for achieving two or more, or all three, of these goals.

As noted herein, the sample collection/transport/storage solutions of the present invention can provide a number of significant improvements over those available in the prior art. Exemplary benefits include, without limitation, one or more of the following: inactivation, killing, and/or lysis of microbes, viruses, or pathogens; destruction and/or inactivation of exogenous or endogenous nucleases, including, without limitation, RNase and/or DNase; compatibility with a variety of conventional nucleic acid extraction, purification, and amplification systems; preservation of RNA and/or DNA integrity within the sample; facilitation of transport and shipping at ambient temperatures, even over extended periods of time, or extreme temperature variations; and suitability for short- (several hours to several days), intermediate- (several days to several weeks), or long- (several weeks to several months) term storage of the isolated nucleic acids.

The disclosed compositions are particularly well suited for point-of-care, field studies, in-home health care or testing, triage/emergency and casualty assessment(s), mobile forensics, pathology, epidemiological sampling, crime-scene investigation, paternity testing, pre- and post-pregnancy genetic screening, rape/incest testing, family counseling, confidential testing of sexually transmitted diseases (including, without limitation, HIV, syphilis, Chlamydia, gonorrhea, and other venereal diseases and the like). Likewise, the disclosed compositions may also be of particular relevance for monitoring, etiology, and control of epidemic or pandemic diseases in human or animal populations, both domestically and abroad, including the collection and analysis of influenzavirus-containing specimens, particularly for predicting and managing viral "shift and drift," and identifying or managing emerging or ongoing epidemics, pandemics, and the like.

In certain embodiments, a nucleic acid isolated using the methods of the present invention may serve as a template in one or more subsequent molecular biological applications, assays, or techniques, including, without limitation, genetic fingerprinting; amplified fragment length polymorphism (AFLP) PCR; restriction fragment-length polymorphism analysis (RFLP); allele-specific oligonucleotide analysis (ASOA); microsatellite analysis; Southern hybridization; Northern hybridization; variable number of tandem repeats (VNTR) PCR; dot-blot hybridization; polymerase cycling assembly (PCA); nested PCR; quantitative PCR (qPCR); asymmetric PCR; DNA footprinting; single nucleotide polymorphism (SNP) genotyping; reverse-transcription PCR (RT-PCR); multiplex PCR (m-PCR); multiplex ligation-dependent probe amplification (MLPA); ligation-mediated PCR (LmPCR); methylation specific PCR (MPCR); helicase-dependent amplification (HDA); overlap-extension PCR (OE-PCR); whole-genome amplification (WGA); plasmid isolation; allelic amplification; site-directed mutagenesis; high-throughput genetic screening; or the like, or any combination thereof.

In one embodiment, the specimen collection/storage/transport formulations preserve, extend, facilitate, or maintain the stability of at least a first nucleic acid in a population of polynucleotides liberated from the specimen for an extended period of time, and particularly, e.g., for periods of time substantially longer than those afforded by conventional sample collection media.

Preferably, the composition will facilitate the preservation and stabilization of one or more nucleic acid segments from within such an isolated population of polynucleotides for periods of time including, without limitation, at least about 7 days, preferably at least about 14 days, more preferably, at least about 21 days, and more preferably still, at least about 28, 35, 42, 56, 63, or about 70 days or more, without need for refrigerating, or freezing, the collected sample—even when the sample is collected, or stored, under ambient environmental conditions. Such improved specimen preservation properties are particularly desirable when such samples are collected or stored in remote locations, or in temperate, sub-tropical or tropical climates.

In another embodiment the invention provides for the use of the disclosed compositions to prepare nucleic acids from a sample of biological origin in maintaining the integrity and fidelity of one or more nucleic acids contained within a population of polynucleotides released from the sample for extended periods including, without limitation, for at least about 5 days, at least about 10 days, preferably at least about 15 days, more preferably at least about 20 days, or more preferably still, at least about 25 days or more, without the requirement of storing the sample under refrigerating or freezing conditions.

Nucleic acids obtained in the practice of the disclosed methods are advantageously compatible with a number of conventional molecular and diagnostic isolation, purification, detection, and/or analytic methodologies. The disclosed compositions facilitate recovery storage, and transport of populations of stabilized, substantially non-degraded, polynucleotides for use in a variety of subsequent methodologies, including,

without limitation, nucleic acid isolation, purification, amplification, and molecular analytical and/or diagnostic testing, assay, analysis, or characterization, and the like.

Exemplary Commercial Kits of the Present Invention

The following outline provides exemplary commercial kits employing the PSS formulations of the present invention (FIG. 5).

Peel-Pouch Collection System

Five-mL tube containing 1.5 mL PrimeStore.TM. Solution; Collection swab (e.g., FlockedSWABS.RTM. [Copan Diagnostics, Inc., Murrieta, Calif., USA]); and instructions for collection and/or processing of samples. (packed, e.g., in 50 pouches/unit) (See FIG. 4A, FIG. 4B, FIG. 4C, FIG. 4D, and FIG. 4E for a schematic demonstration of such systems).

PrimeStore.TM. Stock Solutions

Once formulated, PrimeStore.TM. stock solution is stable at 4.degree. C. or below for periods of at least one year or more. Formulated PrimeStore.TM. solutions have also been shown to be stable at ambient temperature (e.g., about 20-30.degree. C.) for periods of many months.

Once a sample is contacted with a PrimeStore.TM. formulation as disclosed herein, it can reasonably expected to be stored indefinitely at temperatures of 0.degree. C. or below, at least one year or more under refrigeration (e.g., .apprxeq.4.degree. C. and at least 30 days or more at ambient temperature (e.g., about 20-30.degree. C.), without significant loss of nucleic acid composition, fidelity or integrity of the sample. For example, without limitation, the integrity of a population of polynucleotides obtained from the sample is at least substantially maintained, and preferably entirely maintained without detectable degradation, when the composition comprising the sample is stored at a temperature of from about -20.degree. C. to about 40.degree. C., for a period of from about 30 days to about 60 days.

Collection Kits for Environmental Samples

Exemplary commercial packages of the disclosed formulations for the collection of environmental samples include, without limitation, collection vessels containing effective amounts of the storage solution (e.g., 0.1-, 0.2-, 0.5-, 1-, 2-, or 3-mL collection vials each containing 0.1 mL, 0.2 mL, 0.25 mL, 0.5 mL, 0.75 mL, or 1 mL PrimeStore.TM. solution); and instructions for collection and/or processing of samples. The collection vessels may be sized or packaged in smaller or larger sizes as needed depending on the proposed specimen(s) or collection method(s). In certain applications, it may also be desirable to package individual collection devices into multiple containers (including, without limitation, commercial packs containing e.g., 10 or 20 collection devices/unit package).

Sequence Comparison, Identity, and Homology

The terms "identical" or percent "identity," in the context of two or more nucleic acid or polypeptide sequences, refer to two or more sequences or subsequences that are the same or have a specified percentage of amino acid residues or nucleotides that are the same, when compared and aligned for maximum correspondence, as measured using one of the sequence comparison algorithms described below (or other

algorithms available to persons of ordinary skill) or by visual inspection.

The phrase "substantially identical," in the context of two nucleic acids refers to two or more sequences or subsequences that have at least about 90%, preferably 91%, most preferably about 92%, 93%, 94%, 95%, 96%, 97%, 98%, 98.5%, 99%, 99.1%, 99.2%, 99.3%, 99.4%, 99.5%, 99.6%, 99.7%, 99.8%, 99.9% or more nucleotide residue identity, when compared and aligned for maximum correspondence, as measured using a sequence comparison algorithm or by visual inspection. Such "substantially identical" sequences are typically considered "homologous," without reference to actual ancestry.

For sequence comparison and homology determination, typically one sequence acts as a reference sequence to which test sequences are compared. When using a sequence comparison algorithm, test and reference sequences are input into a computer, subsequence coordinates are designated, if necessary, and sequence algorithm program parameters are designated. The sequence comparison algorithm can then be used to calculate the percent sequence identity for the test sequence(s) relative to the reference sequence, based on the designated program parameters.

Optimal alignment of sequences for comparison can be conducted, e.g., by the local homology algorithm (see e.g., Smith and Waterman, 1981), by the homology alignment algorithm (see e.g., Needleman and Wunsch, 1970), by the search similarity comparison method (see e.g., Pearson and Lipman, 1988), by computerized implementations of algorithms such as GAP, BESTFIT, FASTA, and TFASTA in the Wisconsin Genetics Software Package, Genetics Computer Group, Madison, Wis., USA, or by visual inspection. One example of an algorithm that is suitable for determining percent sequence identity and sequence similarity is the BLAST algorithm (Altschul et al., 1990) and BLOSUM62 scoring matrix (see, e.g., Henikoff and Henikoff, 1989). Software for performing BLAST analyses is publicly available through the National Center for Biotechnology Information (www.ncbi.nlm.nih.gov/).

In addition to calculating percent sequence identity, the BLAST algorithm also performs a statistical analysis of the similarity between two sequences (see, e.g., Karlin and Altschul, 1993). Another example of a useful sequence alignment algorithm is the PILEUP program, which creates a multiple sequence alignment from a group of related sequences using progressive, pairwise alignments. It can also plot a tree showing the clustering relationships used to create the alignment. PILEUP uses a simplification of the progressive alignment comparison method (see e.g., Feng and Doolittle, 1987), and employs a general alignment matrix similar to that described by Higgins and Sharp (1989).

EXEMPLARY DEFINITIONS

Unless defined otherwise, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art to which this invention belongs. Although any methods and compositions similar or equivalent to those described herein can be used in the practice or testing of the present invention, the preferred methods and compositions are described herein. For purposes of the present invention, the following terms are defined below:

The terms "about" and "approximately" as used herein, are interchangeable, and should generally be understood to refer to a range of numbers around a given number, as well as to all numbers in a recited range of numbers (e.g., "about 5 to 15" means "about 5 to about 15" unless otherwise stated). Moreover, all numerical ranges herein should be understood to include each whole integer within the range.

As used herein, the term "carrier" is intended to include any solvent(s), dispersion medium, coating(s), diluent(s), buffer(s), isotonic agent(s), solution(s), suspension(s), colloid(s), inert(s) or such like, or a combination thereof that is pharmaceutically acceptable for administration to the relevant animal or acceptable for a diagnostic purpose, as applicable. The use of one or more delivery vehicles for chemical compounds in general, and peptides and epitopes in particular, is well known to those of ordinary skill in the pharmaceutical arts. Except insofar as any conventional media or agent is incompatible with the active ingredient, its use in the diagnostic, prophylactic, and therapeutic compositions is contemplated. One or more supplementary active ingredient(s) may also be incorporated into, or administered in association with, one or more of the disclosed immunogenic compositions.

As used herein, the term "nucleic acid" includes one or more types of: polydeoxyribonucleotides (containing 2-deoxy-D-ribose), polyribonucleotides (containing D-ribose), and any other type of polynucleotide that is an N-glycoside of a purine or pyrimidine base, or modified purine or pyrimidine bases (including abasic sites). The term "nucleic acid," as used herein, also includes polymers of ribonucleosides or deoxyribonucleosides that are covalently bonded, typically by phosphodiester linkages between subunits, but in some cases by phosphorothioates, methylphosphonates, and the like. "Nucleic acids" include single- and double-stranded DNA, as well as single- and double-stranded RNA. Exemplary nucleic acids include, without limitation, gDNA; hnRNA; mRNA; rRNA, tRNA, micro RNA (miRNA), small interfering RNA (siRNA), small nucleolar RNA (snORNA), small nuclear RNA (snRNA), and small temporal RNA (stRNA), and the like, and any combination thereof.

As used herein, the terms "protein," "polypeptide," and "peptide" are used interchangeably, and include molecules that include at least one amide bond linking two or more amino acid residues together. Although used interchangeably, in general, a peptide is a relatively short (e.g., from 2 to about 100 amino acid residues in length) molecule, while a protein or a polypeptide is a relatively longer polymer (e.g., 100 or more residues in length). However, unless specifically defined by a chain length, the terms peptide, polypeptide, and protein are used interchangeably.

As used herein, "sample" includes anything containing or presumed to contain a substance of interest. It thus may be a composition of matter containing nucleic acid, protein, or another biomolecule of interest. The term "sample" can thus encompass a solution, cell, tissue, or population of one or more of the same that includes a population of nucleic acids (genomic DNA, cDNA, RNA, protein, other cellular molecules, etc.). The terms "nucleic acid source," "sample," and "specimen" are used interchangeably herein in a broad sense, and are intended to encompass a variety of biological sources that contain nucleic acids, protein, one or more other biomolecules of interest, or any combination thereof. Exemplary biological samples include, but are not limited to, whole blood, plasma, serum, sputum, urine, stool, white blood cells, red blood cells, buffy coat, swabs (including, without limitation, buccal swabs, throat swabs, vaginal swabs, urethral swabs, cervical swabs, rectal swabs, lesion swabs, abscess swabs, nasopharyngeal swabs, and the like), urine, stool, sputum, tears, mucus, saliva, semen, vaginal fluids, lymphatic fluid, amniotic fluid, spinal or cerebrospinal fluid, peritoneal effusions, pleural effusions, exudates, punctates, epithelial smears, biopsies, bone marrow samples, fluid from cysts or abscess contents, synovial fluid, vitreous or aqueous humor, eye washes or aspirates, pulmonary lavage or lung aspirates, and organs and tissues, including but not limited to, liver, spleen, kidney, lung, intestine, brain, heart, muscle, pancreas, and the like, and any combination thereof. In some embodiments, the sample may be, or be from, an organism that acts as a vector, such as a mosquito, or tick, or other insect(s).

Tissue culture cells, including explanted material, primary cells, secondary cell lines, and the like, as well as

lysates, homogenates, extracts, or materials obtained from any cells, are also within the meaning of the term "biological sample," as used herein. Microorganisms (including, without limitation, prokaryotes such as the archaeobacteria and eubacteria; cyanobacteria; fungi, yeasts, molds, actinomycetes; spirochetes, and mycoplasmas); viruses (including, without limitation the Orthohepadnaviruses [including, e.g., hepatitis A, B, and C viruses], human papillomavirus, Flaviviruses [including, e.g., Dengue virus], Lyssaviruses [including, e.g., rabies virus], Morbilliviruses [including, e.g., measles virus], Simplexviruses [including, e.g., herpes simplex virus], Polyomaviruses, Rubulaviruses [including, e.g., mumps virus], Rubiviruses [including, e.g., rubella virus], Varicellovirus [including, e.g., chickenpox virus], rotavirus, coronavirus, cytomegalovirus, adenovirus, adeno-associated virus, baculovirus, parvovirus, retrovirus, vaccinia, poxvirus, and the like), algae, protozoans, protists, plants, bryophytes, and the like, and any combination of any of the foregoing, that may be present on or in a biological sample are also within the scope of the invention, as are any materials obtained from clinical or forensic settings that contain one or more nucleic acids are also within the scope of the invention. The ordinary-skilled artisan will also appreciate that lysates, extracts, or materials obtained from any of the above exemplary biological samples are also within the scope of the invention.

As used herein, the term "buffer" includes one or more compositions, or aqueous solutions thereof, that resist fluctuation in the pH when an acid or an alkali is added to the solution or composition that includes the buffer. This resistance to pH change is due to the buffering properties of such solutions, and may be a function of one or more specific compounds included in the composition. Thus, solutions or other compositions exhibiting buffering activity are referred to as buffers or buffer solutions. Buffers generally do not have an unlimited ability to maintain the pH of a solution or composition; rather, they are typically able to maintain the pH within certain ranges, for example from a pH of about 5 to 7.

As used herein, the term "biological molecule" refers to any molecule found within a cell or produced by a living organism, including viruses. This may include, but is not limited to, nucleic acids, proteins, carbohydrates, and lipids. As used herein, a "cell" refers to the smallest structural unit of an organism that is capable of independent functioning and is comprised of cytoplasm and various organelles surrounded by a cell membrane. This may include, but is not limited to, cells that function independently such as bacteria and protists, or cells that live within a larger organism such as leukocytes and erythrocytes. As defined herein, a cell may not have a nucleus, such as a mature human red blood cell.

Samples in the practice of the invention can be used fresh, or can be used after being stored for a period of time, or for an extended period of time, including for example, cryopreserved, frozen, or refrigerated samples and the like, and may include material of clinical, veterinary, environmental or forensic origin, may be isolated from food, beverages, feedstocks, potable water sources, wastewater streams, industrial waste or effluents, natural water sources, soil, airborne sources, pandemic or epidemic populations, epidemiological samples, research materials, pathology specimens, suspected bioterrorism agents, crime scene evidence, and the like.

As used herein, the term "patient" (also interchangeably referred to as "host" or "subject") refers to any host that can serve as a source of one or more of the biological samples or specimens as discussed herein. In certain aspects, the donor will be a vertebrate animal, which is intended to denote any animal species (and preferably, a mammalian species such as a human being). In certain embodiments, a "patient" refers to any animal host, including but not limited to, human and non-human primates, avians, reptiles, amphibians, bovines, canines, caprines, cavines, corvines, epines, equines, felines, hircines, lapines, leporines, lupines, murines, ovines, porcines, racines, vulpines, and the like, including, without limitation, domesticated livestock, herding or migratory animals or birds, exotics or zoological specimens, as well as companion

animals, pets, and any animal under the care of a veterinary practitioner. The invention may also be used to monitor disease outbreak, progression, and epidemiological statistics for a variety of global populations, including, without limitation, wasting disease in ungulates, tuberculosis, ebola, SARS, and avian influenzas. In certain embodiments, the samples will preferably be of mammalian origin, and more preferably of human origin.

The term "chaotrope" or "chaotropic agent" as used herein, includes substances that cause disorder in a protein or nucleic acid by, for example, but not limited to, altering the secondary, tertiary, or quaternary structure of a protein or a nucleic acid while leaving the primary structure intact.

The term "e.g.," as used herein, is used merely by way of example, without limitation intended, and should not be construed as referring only those items explicitly enumerated in the specification.

The term "substantially free" or "essentially free," as used herein, typically means that a composition contains less than about 10 weight percent, preferably less than about 5 weight percent, and more preferably less than about 1 weight percent of a compound. In a preferred embodiment, these terms refer to less than about 0.5 weight percent, more preferably less than about 0.1 weight percent or even less than about 0.01 weight percent. The terms encompass a composition being entirely free of a compound or other stated property, as well. With respect to degradation or deterioration, the term "substantial" may also refer to the above-noted weight percentages, such that preventing substantial degradation would refer to less than about 15 weight percent, less than about 10 weight percent, preferably less than about 5 weight percent, etc., being lost to degradation. In other embodiments, these terms refer to mere percentages rather than weight percentages, such as with respect to the term "substantially non-pathogenic" where the term "substantially" refers to leaving less than about 10 percent, less than about 5 percent, etc., of the pathogenic activity.

As used herein, the term "heterologous" is defined in relation to a predetermined referenced nucleic acid sequence. For example, with respect to a structural gene sequence, a heterologous promoter is defined as a promoter that does not naturally occur adjacent to the referenced structural gene, but which is positioned by the hand of man in one or more laboratory manipulations that are routinely employed by those of ordinary skill in the molecular biological arts. Likewise, a heterologous gene or nucleic acid segment is defined as a gene or nucleic acid segment that does not naturally occur adjacent to the referenced sequence, promoter and/or enhancer element(s), etc.

As used herein, the term "homology" refers to a degree of complementarity between two or more polynucleotide or polypeptide sequences. The word "identity" may substitute for the word "homology" when a first nucleic acid or amino acid sequence has the exact same primary sequence as a second nucleic acid or amino acid sequence. Sequence homology and sequence identity can be determined by analyzing two or more sequences using algorithms and computer programs known in the art. Such methods may be used to assess whether a given sequence is identical or homologous to another selected sequence.

As used herein, "homologous" means, when referring to polynucleotides, sequences that have the same essential nucleotide sequence, despite arising from different origins. Typically, homologous nucleic acid sequences are derived from closely related genes or organisms possessing one or more substantially similar genomic sequences. By contrast, an "analogous" polynucleotide is one that shares the same function with a polynucleotide from a different species or organism, but may have a significantly different primary nucleotide sequence that encodes one or more proteins or polypeptides that accomplish similar functions or possess similar biological activity. Analogous polynucleotides may often be derived from two or more organisms that

are not closely related (e.g., either genetically or phylogenetically).

The terms "identical" or percent "identity," in the context of two or more nucleic acid or polynucleotide sequences, refer to two or more sequences or subsequences that are the same or have a specified percentage of nucleotides that are the same, when compared and aligned for maximum correspondence over a comparison window, as measured using a sequence comparison algorithm or by manual alignment and visual inspection.

A "primer" or "primer sequence" may include any nucleic acid sequence or segment that selectively hybridizes to a complementary template nucleic acid strand ("target sequence") and functions as an initiation point for the addition of nucleotides to replicate the template strand. Primer sequences of the present invention may be labeled or contain other modifications which allow the detection and/or analysis of amplification products. In addition to serving as initiators for polymerase-mediated duplication of target DNA sequences, primer sequences may also be used for the reverse transcription of template RNAs into corresponding DNAs.

A "target sequence" or "target nucleotide sequence" as used herein includes any nucleotide sequence to which one of said primer sequences hybridizes under conditions that allow an enzyme having polymerase activity to elongate the primer sequence, and thereby replicate the complementary strand.

As used herein, the term "operably linked" refers to a linkage of two or more polynucleotides or two or more nucleic acid sequences in a functional relationship. A nucleic acid is "operably linked" when it is placed into a functional relationship with another nucleic acid sequence. For instance, a promoter or enhancer is operably linked to a coding sequence if it affects the transcription of the coding sequence. "Operably linked" means that the nucleic acid sequences being linked are typically contiguous, or substantially contiguous, and, where necessary to join two protein coding regions, contiguous and in reading frame. Since enhancers generally function when separated from the promoter by several kilobases and intronic sequences may be of variable lengths; however, some polynucleotide elements may be operably linked but not contiguous.

"Transcriptional unit" refers to a polynucleotide sequence that comprises at least a first structural gene operably linked to at least a first cis-acting promoter sequence and optionally linked operably to one or more other cis-acting nucleic acid sequences necessary for efficient transcription of the structural gene sequences, and at least a first distal regulatory element as may be required for the appropriate tissue-specific and developmental transcription of the structural gene sequence operably positioned under the control of the promoter and/or enhancer elements, as well as any additional cis sequences that are necessary for efficient transcription and translation (e.g., polyadenylation site(s), mRNA stability controlling sequence(s), etc.

The phrases "isolated" or "biologically pure" refer to material that is substantially, or essentially, free from components that normally accompany the material as it is found in its native state. Thus, isolated polynucleotides in accordance with the invention preferably do not contain materials normally associated with those polynucleotides in their natural, or in situ, environment.

"Link" or "join" refers to any method known in the art for functionally connecting one or more proteins, peptides, nucleic acids, or polynucleotides, including, without limitation, recombinant fusion, covalent bonding, disulfide bonding, ionic bonding, hydrogen bonding, electrostatic bonding, and the like.

The term "pathogen" is defined herein as any sort of infectious agent, including e.g., viruses, prions,

protozoans, parasites, as well as microbes such as bacteria, yeast, molds, fungi, and the like.

As used herein, the term "plasmid" refers to a genetic construct that is composed of genetic material (i.e., nucleic acids). Typically, a plasmid contains an origin of replication that is functional in bacterial host cells, e.g., *Escherichia coli*, and selectable markers for detecting bacterial host cells comprising the plasmid. Plasmids of the present invention may include genetic elements as described herein arranged such that an inserted coding sequence can be transcribed and translated in eukaryotic cells. In addition, the plasmid may include one or more nucleic acids derived from natural or artificial sources.

As used herein, the term "polypeptide" is intended to encompass a singular "polypeptide" as well as plural "polypeptides," and includes any chain or chains of two or more amino acids. Thus, as used herein, terms including, but not limited to "peptide," "dipeptide," "tripeptide," "protein," "enzyme," "amino acid chain," and "contiguous amino acid sequence" are all encompassed within the definition of a "polypeptide," and the term "polypeptide" can be used instead of, or interchangeably with, any of these terms. The term further includes polypeptides that have undergone one or more post-translational modification(s), including for example, but not limited to, glycosylation, acetylation, phosphorylation, amidation, derivatization, proteolytic cleavage, post-translation processing, or modification by inclusion of one or more non-naturally occurring amino acids. Conventional nomenclature exists in the art for polynucleotide and polypeptide structures. For example, one-letter and three-letter abbreviations are widely employed to describe amino acids: Alanine (A; Ala), Arginine (R; Arg), Asparagine (N; Asn), Aspartic Acid (D; Asp), Cysteine (C; Cys), Glutamine (Q; Gln), Glutamic Acid (E; Glu), Glycine (G; Gly), Histidine (H; His), Isoleucine (I; Ile), Leucine (L; Leu), Methionine (M; Met), Phenylalanine (F; Phe), Proline (P; Pro), Serine (S; Ser), Threonine (T; Thr), Tryptophan (W; Trp), Tyrosine (Y; Tyr), Valine (V; Val), and Lysine (K; Lys). Amino acid residues described herein are preferred to be in the "L" isomeric form. However, residues in the "D" isomeric form may be substituted for any L-amino acid residue provided the desired properties of the polypeptide are retained.

"Protein" is used herein interchangeably with "peptide" and "polypeptide," and includes both peptides and polypeptides produced synthetically, recombinantly, or in vitro and peptides and polypeptides expressed in vivo after nucleic acid sequences are administered into a host animal or human subject. The term "polypeptide" is preferably intended to refer to all amino acid chain lengths, including those of short peptides of from about 2 to about 20 amino acid residues in length, oligopeptides of from about 10 to about 100 amino acid residues in length, and polypeptides including about 100 amino acid residues or more in length. The term "sequence," when referring to amino acids, relates to all or a portion of the linear N-terminal to C-terminal order of amino acids within a given amino acid chain, e.g., polypeptide or protein; "subsequence" means any consecutive stretch of amino acids within a sequence, e.g., at least 3 consecutive amino acids within a given protein or polypeptide sequence. With reference to nucleotide and polynucleotide chains, "sequence" and "subsequence" have similar meanings relating to the 5' to 3' order of nucleotides.

As used herein, the term "substantially homologous" encompasses two or more biomolecular sequences that are significantly similar to each other at the primary nucleotide sequence level. For example, in the context of two or more nucleic acid sequences, "substantially homologous" can refer to at least about 75%, preferably at least about 80%, and more preferably at least about 85% or at least about 90% identity, and even more preferably at least about 95%, more preferably at least about 97% identical, more preferably at least about 98% identical, more preferably at least about 99% identical, and even more preferably still, at least substantially or entirely 100% identical (i.e., "invariant").

Likewise, as used herein, the term "substantially identical" encompasses two or more biomolecular sequences

(and in particular polynucleotide sequences) that exhibit a high degree of identity to each other at the nucleotide level. For example, in the context of two or more nucleic acid sequences, "substantially identical" can refer to sequences that at least about 80%, and more preferably at least about 85% or at least about 90% identical to each other, and even more preferably at least about 95%, more preferably at least about 97% identical, more preferably at least about 98% identical, more preferably at least about 99% identical, and even more preferably still, at least substantially or entirely 100% identical (i.e., "non-degenerate").

The term "recombinant" indicates that the material (e.g., a polynucleotide or a polypeptide) has been artificially or synthetically (non-naturally) altered by human intervention. The alteration can be performed on the material within or removed from, its natural environment or state. Specifically, e.g., a promoter sequence is "recombinant" when it is produced by the expression of a nucleic acid segment engineered by the hand of man. For example, a "recombinant nucleic acid" is one that is made by recombining nucleic acids, e.g., during cloning, DNA shuffling or other procedures, or by chemical or other mutagenesis; a "recombinant polypeptide" or "recombinant protein" is a polypeptide or protein which is produced by expression of a recombinant nucleic acid; and a "recombinant virus," e.g., a recombinant influenza virus, is produced by the expression of a recombinant nucleic acid.

In accordance with the present invention, polynucleotides, nucleic acid segments, nucleic acid sequences, and the like, include, but are not limited to, DNAs (including and not limited to genomic or extragenomic DNAs), genes, peptide nucleic acids (PNAs) RNAs (including, but not limited to, rRNAs, mRNAs and tRNAs), nucleosides, and suitable nucleic acid segments either obtained from natural sources, chemically synthesized, modified, or otherwise prepared or synthesized in whole or in part by the hand of man.

Likewise, in accordance with long standing patent law convention, the words "a" and "an" when used in this application, including the claims, denotes "one or more."

EXAMPLES

The following examples are included to demonstrate illustrative embodiments of the invention. It should be appreciated by those of ordinary skill in the art that the techniques disclosed in the examples that follow represent techniques discovered to function well in the practice of the invention, and thus can be considered to constitute preferred modes for its practice. However, those of ordinary skill in the art should, in light of the present disclosure, appreciate that many changes can be made in the specific embodiments which are disclosed and still obtain a like or similar result without departing from the spirit and scope of the invention.

Example 1

Formulation of Exemplary Storage Solutions

The present example provides a general formulation of the PSS (PanFlu) compositions of the present invention. Additional formulations of the PSS compositions are also exemplified in Example 2 through Example 5.

Materials

Guanidine thiocyanate, sodium citrate, Antifoam A.RTM. Concentrate, and N-lauroylsarcosine, sodium salt, were all purchased from Sigma Chemical Co. (St. Louis, Mo., USA). Tris(2-carboxyethyl) phosphine

hydrochloride (TCEP) was obtained from Soltec Ventures Inc. (Beverly, Mass., USA). 2-amino-2-hydroxymethyl-propane-1,3-diol (TRIS) was obtained from Applied Biosystems/Ambion (Austin, Tex., USA). 2-[2-(Bis(carboxymethyl)amino)ethyl-(carboxymethyl)amino]acetic acid (EDTA) GIBCO.RTM. Ultra Pure was obtained from Invitrogen Corp. (Carlsbad, Calif., USA). All other reagents are available commercially from Sigma-Aldrich or USB Corporation.

TABLE-US-00001 TABLE 1 FORMULATION RANGES OF EXEMPLARY COMPONENTS FOR THE PREPARATION OF PRIMESTORE .TM. COMPOSITIONS

Compound	Final Concentration	Range
1. A chaotrope, e.g.: Guanidine thiocyanate	about 0.5 M to about 6 M	or Guanidine hydrochloride about 0.5 M to about 6 M
2. An anionic detergent, e.g.: N-lauroyl sarcosine (inter alia Na salt)	about 0.15% to about 1% (wt./vol.)	or Sodium dodecyl sulfate, same Lithium dodecyl sulfate, same Sodium glycocholate, same Sodium deoxycholate, same Sodium taurodeoxycholate, or same Sodium cholate
3. A reducing agent, e.g.: TCEP	about 0.5 mM to about 30 mM	or .beta.-ME, DTT, formamide, or DMSO about 0.05 M to about 0.3 M
4. A chelator, e.g.: Sodium citrate	about 0.5 mM to about 50 mM	or EGTA, HEDTA, DTPA, NTA, APCA, etc. about 0.01 mM to about 1 mM
5. A buffer (e.g., TRIS, HEPES, Bis-Tris, etc.)	about 1 mM to about 1 M	
6. An acid (e.g., HCl or citric acid)	q.s. to adjust to a pH of about 6 to 7, preferably 6.4 to 6.8	
7. Nuclease-free water	q.s. to desired final volume	
8. A surfactant/defoaming agent, e.g.: Antifoam A .RTM. or Tween .RTM.	about 0.0001% to about 0.3% (wt./vol.)	
9. An alkanol (e.g., methanol, ethanol, propanol, etc.)	about 1% to about 25% (vol./vol.)	
10. carrier/IPC RNA or DNA	about 1 pg to about 1 .mu.g/mL	

Example 2

Formulation of an Exemplary Storage Solution

The present example describes a first exemplary formulation of the compositions of the invention. This formulation has also been alternatively referred to by the inventors as "PrimeStore.TM. Solution" or "PSS" version 1.

TABLE-US-00002 TABLE 2 PREPARATION OF PRIMESTORE .TM. COMPOSITION (VER. 1)

Compound	Final Concentration
Guanidine thiocyanate	4 M
Sodium citrate	30 mM
Sodium dodecyl sulfate	0.25% (wt/vol.)
N-lauroyl sarcosine, sodium salt	0.25% (wt/vol.)
2-mercaptoethanol (.beta.-ME)	0.1 M
Antifoam A	0.1% (wt/vol.)
Citric acid	q.s. to adjust pH to 6.5
Nuclease-free water	11.82 mL

Example 3

Preparation of a Second Exemplary Storage Solution

The present example describes the preparation of another exemplary storage solution according to the present invention. This formulation has also been alternatively referred to by the inventors as PrimeStore.TM. version 2.

TABLE-US-00003 TABLE 3 PREPARATION OF PRIMESTORE .TM. COMPOSITION (VER. 2)

Compound	Quantity	Final Concentration
Guanidine thiocyanate	35.488 gm	3 M
TCEP	0.02867 gm	1 mM
Sodium citrate	0.2931 gm	10 mM
N-lauroyl sarcosine, sodium salt (NLS)	0.5 gm	0.5% salt (NLS)
Antifoam A (10% solution)	200 .mu.L	0.002%
TRIS (1 M)	10 mL	100 mM
EDTA (0.5 M)	20 .mu.L	0.1 mM
Hydrochloric acid (HCl)	q.s. to adjust pH to 6.7	--
Nuclease-free water	q.s. to 100 mL	--

Example 4

Preparation of a Third Exemplary Storage Solution

The present example describes the preparation of another exemplary storage solution according to the present invention. This formulation has also been alternatively referred to by the inventors as PSS version 2.2.

TABLE-US-00004 TABLE 4 PREPARATION OF PRIMESTORE .TM. COMPOSITION (VER. 2.2)

Compound	Quantity	Concentration
Guanidine thiocyanate	35.488 gm	2.5 M
TCEP	0.02867 gm	0.5 mM
Sodium citrate	0.2931 gm	10 mM
N-lauroyl sarcosine, sodium salt (NLS)	0.5 gm	0.4%
Antifoam A (10% solution)	200 .mu.L	0.002%
TRIS (1 M)	10 mL	100 mM
EDTA (0.5 M)	20 .mu.L	0.1 mM
Ethanol, molecular grade	23 mL	23% (vol./vol.) (96-100%)
Hydrochloric acid (HCl)	q.s.	to adjust pH to 6.7
Nuclease-free water	q.s.	to 100 mL

Exemplary Protocol for the Preparation of PSS (Ver. 2.2): 1. Add 40 mL of nuclease-free water to a clean beaker with a stir bar. 2. Place beaker on a hot plate/stirrer and adjust temperature to 60-65.degree. C. Set stirring speed to medium. 3. Add 35.488 gm of guanidine thiocyanate slowly to the water allowing it to dissolve as added. 4. Add 0.0287 gm of TCEP to beaker and increase stirrer speed to help dissolve crystals. 5. Add 0.2931 gm of sodium citrate to the beaker. 6. Add 0.5 gm of NLS to the solution. Increase stirrer speed to create a vortex in the beaker. This will pull the NLS into the solution and help dissolve the reagent. 7. Vortex a prepared 10% Antifoam A solution (1 mL Antifoam A Concentrate+9 mL nuclease-free water). Pipette 200 .mu.L of the 10% Antifoam A into the solution. 8. Pipette 10 mL of 1 M TRIS into the solution. 9. Pipette 20 .mu.L of 0.5 M EDTA into the solution. 10. Increase the temperature to bring the solution to 75-80.degree. C. and stir for 3-5 minutes. 11. Remove beaker from heat and allow solution to cool to room temperature (.apprxeq.22-25.degree. C.). 12. Add 23 mL of ethanol to the solution and mix thoroughly. 13. Adjust pH to 6.9 with HCl. 14. Pour solution into a clean 100 mL graduated cylinder. 15. Add nuclease-free water to bring total volume to 100 mL. 16. Transfer solution to a labeled sterile container. Store at room temperature (.apprxeq.22-25.degree. C.). *Note: Preferably, make sure each reagent is completely dissolved before adding the next.

Example 5

Comparison of PrimeStore.TM. Solution to Conventional Formulations

A sample of homogenized nasal tissue from a cotton rat (*Sigmodon hispidus*) challenged with influenza A (H3N2) or a human clinical influenza A (H1N1) sample collected as a human clinical nasal wash during the 2006-07 season were placed in PSS (ver. 1) and tested compared to the respective lysis formulation and protocol, and extraction procedure, from three commercially available kits: RNAqueous.RTM.-Micro (Ambion, Austin, Tex., USA), QIAamp Viral RNA Mini Kit (Catalogue #52904, Qiagen, Valencia, Calif., USA), and MagMax AI/ND Viral RNA Isolation Kit (Catalogue #AM1929, Ambion). Extraction efficiency was evaluated using the ABI7500 sequence detection system with the comparative C.sub.T method (see FIG. 2).

In FIG. 2, "delta Rn" represents the fluorescent reporter signal minus a baseline amount. As shown in FIG. 1 and FIG. 2, the relative C.sub.T scores and viral copies detected were optimal when the fixing formulation was used in place of the respective conventional lysis buffer for each commercial kit. In these two sample

types, the compositions of the invention worked better than the two conventional Kits for extraction purposes. The PrimeStore.TM. Solution (ver. 1) composition was also shown to be readily compatible with commercially available nucleic acid extraction kits. FIG. 1 illustrates RNA extraction results where the version 1 of PSS was used in conjunction with three commercially available kits: Qiagen Viral Mini, Ambion RNAqueous Mini, and Ambion AI/NCD MagMax. As illustrated by FIG. 1, when the lysis buffer of the extraction kit was replaced with the fixing formulation (denoted on the figure as "One-Step+), superior nucleic acid extraction was achieved when compared to extraction using kits according to standard protocol (denoted on the figure as "One-Step-". Extraction efficiency was measured by real time (r) reverse transcription (RT) polymerase chain reaction (PCR) [rRT-PCR].

FIG. 3 shows preservation of naked RNA in PrimeStore.TM. Solution compared to preservation in a prior solution, with water used as a control. As illustrated in FIG. 3, detection (by fluorescence) occurred at the earliest amplification cycle for RNA stored in PSS (ver. 1) at all time-points assayed.

Example 6

PrimeStore.TM. Solution for the Collection of Nasal Wash Specimens

A prospective clinical detection study was conducted using nasal wash specimens from: 1) symptomatic pediatric patients and 2) asymptomatic or symptomatic family members. Detection of influenza virus compared nasal wash specimens collected in PrimeStore.TM. Solution and VTM by qRT-PCR and traditional culture, respectively. Genetic characterization of influenza genes encoding hemagglutinin (HA), neuraminidase (NA), and matrix surface (MA) proteins were performed using select nasal wash specimens preserved in PrimeStore.TM. Solution to evaluate vaccine effectiveness and drug sensitivity within viral strains.

Influenza is a highly evolving, RNA-based respiratory virus responsible for more than 200,000 hospitalizations and about 36,000 fatalities each year in the United States. Widespread emergence of influenza drift variants among contemporary circulating human viruses prompted a change in all three-vaccine components for the upcoming 2008/09 season. Increased morbidity and mortality during the 2007/08 season included 72 influenza-associated pediatric deaths and continued drug resistance (oseltamivir [TamiFlu.RTM., Roche Laboratories, Inc., Nutley, N.J., USA] and adamantadine) within circulating strains.

Materials and Methods

A total of 100 pediatric (index) patients who met the clinical case criteria for influenza infection and 126 family contacts were enrolled in the study. Nasal washings were placed into PrimeStore.TM. Solution and Universal VTM and analyzed by rRT-PCR or culture analysis, respectively. rRT-PCR was performed using influenza type (A or B) and subtype (H3, H1, H5) specific primers/probes according to Daum et al. (2007). Further genetic characterization of selected clinical samples preserved in PSS was performed using standard RT-PCR and direct nucleotide sequencing of the hemagglutinin HA, NA, and MA viral proteins.

Results

Of the total samples evaluated (N=226; 100 index, 126 family contacts), 66 (29%) tested positive for influenza virus (45 H3N2, 2 H1N1 and 19 B) by rRT-PCR. rRT-PCR from nasal washings preserved in PSS detected influenza virus from 11 patients (9 Flu A and 2 Flu B) that were not detected by culture (Table 5 and

Table 6). Of these 11 specimens, five were from patients enrolled as family contacts.

Phylogenetic analysis of influenza A and B HA genes exhibited drifting compared to the 2007/08 vaccine strains and revealed a higher genetic homology to the 2008/09 Brisbane vaccine strains. Some genetic differences in viruses were noted among family members, particularly among influenza A (H3N2) strains. MA analysis revealed adamantane resistance in all influenza A H3N2 strains, but sensitivity in both H1N1 viruses. All influenza B strains (n=18) were sensitive to the neuraminidase inhibitor drugs zanamivir (Relenza.RTM. GlaxoSmithKline, Research Triangle Park, N.C., USA) and oseltamivir (Tamiflu.RTM. Roche) based on the presence of an aspartic acid (D) at amino acid 197 (influenza B numbering) in the NA gene.

qRT-PCR vs. Traditional Culture Methods

qRT-PCR is superior to traditional culture for the detection of influenza virus from original nasal wash specimens preserved in PSS: influenza was detected within 2 hours (c.f 2 to 7 days for conventional culture methods); and the analyses were more sensitive (11 specimens; 9 Flu A and 2 Flu B detected below culture limits). Moreover, the use of molecular diagnostic methods in lieu of conventional organism culture did not propagate potentially infectious viruses, and simultaneously provided the type and subtype of the influenza virus.

Genetic Analysis

Vaccine Relatedness

H3N2 Strains: Analysis of the HA1 gene of the influenza A (H3N2) hemagglutinin (HA) revealed genetic drift including five amino acid differences in all Texas strains compared to the 2007-08 A/Wisconsin/67/2005 vaccine strain. One HA1 mutation noted in all Texas strains (D122N) was a potential glycosylation site. All A/Texas (H3N2) strains exhibited a greater HA homology (99.0-99.7%) to the newly selected 2008-09 A/Brisbane/10/2007 strain.

H1N1 Strains: The hemagglutinin HA1 gene of the 2 influenza A (H1N1) exhibited 7 amino acid changes compared to A/Solomon Island/3/2006 vaccine strain. Four substitutions (R90K, T145V, K210T and E290K) were within known H1 antibody combining sites. Both Texas H1N1 strains exhibited greater HA homology (98.8% and 99.4%) to the newly selected 2008-09 A/Brisbane/59/2007 vaccine strain.

Influenza B strains: Analysis of the HA1 hemagglutinin and neuraminidase genes revealed all Texas strains were of the B/Yamagata lineage and genetically more homologous to the 2008-09 B/Brisbane/5/2007 vaccine strain than the 2007/08 B/Malaysia/2506/2004 vaccine strain.

Family Mutation

Amino acid changes were noted in the NA, HA1, M1 and M2e among family members. The HA1 Hemagglutinin showed the highest mutation of the influenza genes analyzed, with one family exhibiting five amino acid changes.

Analysis of the highly conserved 24 amino acid M2e viral surface proton pump showed some variation within families. One index patient strain contained 3 unique amino acid M2e substitutions that were `wild-

type` within family member strains.

Antiviral Susceptibility

Adamantane: Matrix gene (MA) genetic analysis, specifically a serine-to-asparagine substitution at position 31 (S31N), revealed adamantane resistance in all influenza A (H3N2) strains but sensitivity in both influenza A (H1N1) viruses.

Neuraminidase Inhibitors: All Texas influenza A (H3N2) isolates were shown to be sensitive to oseltamivir through genetic analysis of E119V, R292K, and N294S substitutions in the NA gene. Genetic analysis of the influenza B NA gene revealed that all Texas strains contained an aspartic acid (D) residue at position 197, and are thus likely sensitive to oseltamivir.

The protocols and tests herein can be adapted for other microbes like tuberculosis, malaria, staphylococcus, and the like and other pathogens where there is a need to know antimicrobial susceptibility quickly.

Example 7

Influenza Sample Collection Using PrimeStore.TM. Solutions

The compositions of the present invention provide a single sample collection, transport, and storage reagent that facilitate: 1) procuring high quality nucleic acids from clinical or environmental specimens, 2) inactivation of potentially infectious biological pathogens for safe handling and transport of specimens, and 3) stabilization and preservation of released `naked` RNA/DNA preventing hydrolysis/nuclease degradation for prolonged periods at ambient temperatures. The results of one such study are presented in the following example. This example illustrates the effectiveness of PSS (ver. 2.2) in killing pathogenic microbe(s).

TABLE-US-00005 TABLE 5 Influenza Subtype Detection: rRT-PCR vs. Culture Total Flu A Influenza A Influenza B Samples (N = 47) rRT-PCR (N = 47) Culture (N = 40) rRT-PCR (N = 19) Culture (N = 17) Total Flu B Samples (N = 19) Index Patients (28) 28/28 (100%) 23/28 (82%) 16/16 (100%) 14/16 (88%) Index Patients (16) Family Contacts (19) 19/19 (100%) 15/19 (79%) 3/3 (100%) 3/3 (100%) Family Contacts (3)

TABLE-US-00006 TABLE 6 Positive Influenza Detection: rRT-PCR vs. Culture Total Samples Flu Positive rRT-PCR Culture (N = 226) (N = 66) (N = 66) (N = 66) Index Patients (100) 44/100 39/39 (100%) 37/39 (94%) Family Contacts (126) 22/126 27/27 (100%) 18/27 (67%)

Methods

qRT-PCR was used to assay influenza A (H5N1) virus nucleic acid preserved in PrimeStore.TM. Solution. A time-course study at room temperature was carried out to evaluate the integrity of clinical specimens, cloacal samples, and cloned template avian influenza A virus (H5) RNA stored and extracted from PSS, VTM, RNA storage solution, or nuclease-free water. PrimeStore.TM. Solution extraction efficiency was compared to three commercially available nucleic acid extraction kits. Furthermore, the ability of RNA contained in PSS to resist nuclease degradation was evaluated.

Results

PSS (ver. 2, but lacking ethanol) inactivated microbial agents while preserving released RNA/DNA from clinical material, i.e., nasal washes, throat swabs, or environmental samples. Clinical specimens or environmental samples placed in this solution were stabilized at room temperature for up to 30 days while degradation of nucleic acids occurred in other transport media. PSS is compatible with commercially available RNA isolation kits and produced an increased nucleic acid yield.

Example 8

Killing of MRSA (ATCC33592) in PrimeStore.TM. Solution

This example illustrates the effectiveness of PSS (ver. 2.2) in killing a bacterial pathogen. Methicillin-resistant *Staphylococcus aureus* (MRSA) strain ATCC33592 was diluted 10-fold and 1000-fold into PSS (ver. 2.2) and quantitated (see FIG. 6).

Experimental Protocol

Day Procedure

0 Transfer MRSA (ATCC33592) from a Culti-loop.RTM. (Remel) to 1.5 mL of TSB in a 15-ml conical test tube. Incubate at 37.degree. C. for approximately 15 min. Gently vortex suspension and transfer 100 .mu.L to a blood-agar plate. Incubate the plate overnight at 37.degree. C.

1 Observe heavy and uniform colony growth after 12 hr incubation. Transfer .about.10% of colonies to 300 mL of tryptic soy broth (TSB) in a sterile, 1-liter flask. Place flask on shaker at 37.degree. C. and 200 rpm.

After approximately 4-6 hrs' incubation, transfer .about.50 mL of bacterial suspension to new 1-liter flask containing 300 mL fresh TSB.

After approximately 4-6 hrs' incubation, transfer .about.100 .mu.L of culture into 900 .mu.L of TSB (1:10 control dilution). From this suspension 10 .mu.L were transferred to 990 .mu.L of TSB (1:1000 control dilution).

Transfer 100 .mu.L of culture into 900 .mu.L of PrimeStore.TM. Solution (1:10 PSS dilution). From this suspension 10 .mu.L were transferred to 990 .mu.L of TSB (1:1000 PSS dilution).

Immediately after transfer into TSB or PrimeStore.TM. Solution (ver. 2.2), the suspensions were gently vortexed and 100 .mu.L were plated from both dilutions of control and PrimeStore.TM. suspensions onto blood agar plates. The time-zero time point was actually about two minutes following addition of the bacteria to the TSB or PSS.

The suspensions in TSB and PSS (ver. 2.2) were kept at room temperature.

An additional 100 .mu.L was plated onto blood agar plates at 5, 15, 30, 60, 120 and 240 minutes after the preparation of the suspensions in TSB and PSS.

Bacterial suspensions on the plates were allowed to dry, the plated inverted and the plates incubated overnight at 37.degree. C.

A titration of the shaker culture was also performed by mixing 100 μ L of the suspension from the shaker culture with 900 μ L of TSB (10⁻¹ dilution). Serial 10-fold dilutions were prepared through 10⁻⁹. 100 μ L samples were plated onto blood agar from all dilutions except 10⁻¹.

2 Plates were observed for bacterial colonies. All plates were stored at 4.degree. C. for later observation, if necessary.

Results

The results are presented in Table 7 and Table 8. Briefly, the bacterial suspension contained approximately 4.7.times.10⁹ cfu/ml. Thus, the 1:10 dilution contained approximately 4.7.times.10⁸ cfu/ml and the 1:1000 dilution contained 4.7.times.10⁶ cfu/ml. At all time points, the bacteria suspended in TSB were too numerous to count. At all time points the bacteria suspended in PrimeStore.TM. Solution and plated onto blood agar plates had no detectable colonies.

TABLE-US-00007 TABLE 7 KILLING OF MRSA (ATCC33592) BY PRIMESTORE .TM. SOLUTION (VER. 2.2) Incubation Time In TSB In Primestore (Minutes) 1:10 1:1000 1:10 1:1000 0 TNTC TNTC 0 0 5 TNTC TNTC 0 0 15 TNTC TNTC 0 0 30 TNTC TNTC 0 0 60 TNTC TNTC 0 0 120 TNTC TNTC 0 0 240 TNTC TNTC 0 0 TNTC = too numerous to count.

TABLE-US-00008 TABLE 8 TITRATION OF MRSA ATCC33592 FROM SUSPENSION CULTURE Dilution CFU/plate CFU/ml 1.E+01 TNTC 1.E+02 TNTC 1.E+03 TNTC 1.E+04 TNTC 1.E+05 TNTC 1.E+06 TNTC 1.E+07 35 3.5 .times. 10⁹ 1.E+08 6 6 .times. 10⁹ 1.E+09 0 NOTE: CFU/ml calculations are corrected to include the plating volume of 0.1 mls Final Conc: 4.7 .times. 10⁹/ml TNTC = too numerous to count. CFU = colony forming units.

An additional study was performed to determine the time of exposure necessary for killing MRSA ATCC33592 when diluted 10-fold into PSS (Ver. 2.2), and to determine the effect of dilution of the bacteria after exposure to PSS, but before plating.

Experimental Protocol

Day Procedure

0 Transfer MRSA (ATCC33592) from TNTC plate from the study described above into 4 mL of TSB. These plates had been stored at 4.degree. C. for approximately 48 hrs. Bacteria were vortexed gently and placed at room temperature for approximately 10 min. before use. 0.1 mL of bacterial suspension was transferred to 0.9 mL PSS and vortexed gently. After approximately 60 sec, the bacteria in PrimeStore.TM. were again vortexed gently and 0.1 mL of bacterial suspension was transferred into 0.3 mL of TSB (1:4 dilution). 100 μ L of bacteria in PrimeStore.TM. Solution (designated "neat") and from the 1:4 dilution into TSB were plated onto blood agar plates (5% sheep RBCs in TSA). This process was repeated at 5 and 15 min., and then again with dilutions made into TSB instead of PrimeStore.TM. Solution. The liquid bacterial suspensions on the blood agar plates were allowed to dry at room temperature and then incubated overnight at 37.degree. C. 1 After approximately 16 hrs' incubation, the plates were removed from the incubator and colonies counted.

Results

The bacterial suspension contained an unknown number of colony forming units (cfu) per mL. At all time points the bacteria suspended in tryptic soy broth (TSB) were too numerous to count (TNTC). At all time points, the bacteria suspended in PrimeStore.TM. compositions and plated onto blood agar plates produced no colonies (Table 9).

TABLE-US-00009 TABLE 9 KILLING OF MRSA (ATCC33592) BY PRIMESTORE .TM. SOLUTION
Incubation Time In TSB In TSB (minutes) neat 1:4 neat 1:4 1 TNTC TNTC 0 0 5 TNTC TNTC 0 0 15 TNTC
TNTC 0 0 TNTC = too numerous to count.

Example 9

Additional Studies Evaluating PrimeStore.TM. Solutions

The data in FIG. 7B illustrate the ability of PSS to inactivate microbes. Shown is a study in which chicken cloacal specimens were collected in PSS (Ver. 1). PrimeStore.TM. Solution inactivated the microbial agents in .ltoreq.1 hr. Four original chicken cloacal samples were immersed in PSS or water and subsequently plated on blood agar plates. These results showed that the disclosed composition could quickly kill or inactive microorganisms in the sample.

The data in FIG. 7C illustrate the ability of PSS to inhibit RNA base hydrolysis for 30 days at room temperature. RNA was incubated at room temperature (22-26.degree. C.) in PrimeStore.TM. (gel lane 1 and 3) and water (gel lane 2 and 4), and subsequently RT-PCR amplified (1500 base pairs) at Day 0 and Day 30. PSS preserved collected RNA, and prevented RNA/DNA degradation at room temperature up to 30 days (see also, e.g., Table 11).

Flu Inhibition Assay

The reagents for this assay include:

Trypsin medium containing: (a) 45 mL Sterile N/C EMEM; (b) 3 mL stock 7.5% Na Bicarbonate (2%); (c) 1.5 mL SPG (1%); (d) 75 .mu.L trypsin (0.05%); (e) 1.5 mL fungizone (1%); and (f) 150 .mu.L gentamicin.

Crystal violet containing: (a) 150 ml glutaric dialdehyde; (b) 2 gm crystal violet; and (c) 2850 mL deionized water

Preparation of Serum Samples for Assay

Thaw and vortex serum samples: For each sample, label the lid of a corresponding Spin-X.TM. tube. Combine 450 .mu.L non-complete EMEM with 50 .mu.L serum into a Spin-X.TM. tube. Warm tubes containing the sera and EMEM in a 56.degree. C. water bath for 30 min. Centrifuge tubes at 8000 RPM for 2 min. at room temperature. Label and place samples into a -20.degree. C. freezer until assayed.

Dilution Plates

Load 160 .mu.L of each neat compound or serum sample into wells A1 through A12. Load the remaining wells with 120 .mu.L trypsin medium. Using a multi-channel pipette, draw 40 .mu.L of neat sample from row A and dilute into the corresponding wells in row B. Repeat dilution for each row, mixing well after each

transfer. At row H, after mixing the transfer from row G, draw up 40 μL from each well and discard. Obtain virus stock (10^{sup.6}) from -80^{degree}. C. freezer and thaw. Dilute virus stock in trypsin media to a 10^{sup.3} dilution. After serial dilutions are completed, transfer 120 μL of influenza virus (10^{sup.4} TCID₅₀ per ml) to all wells in the dilution plate. This results in a total of 240 μL in all wells. Incubate dilution plate(s) at room temperature for 1 hour.

MDCK Cell Plates

Sterilize and place glass reservoir, comb dispenser and tubing inside the fume hood. Inside the hood, connect the tubing to the reservoir and fill nozzle of the comb. Connect the aspirator tube to the vacuum nozzle on the comb. Place the reservoir on an elevated surface and turn on the aspirator. Put PBS into the reservoir (1 L or more may be needed depending on the number of plates). Wash the cell plates 3 times with the PBS comb (aspirate the medium, then press the button for roughly 1 second to wash the wells, repeat twice). Using a multi-channel pipette, transfer 50 μL from each well in column 1 of the dilution plate to columns 1 through 4 of the cell plate. Transfer 50 μL from each well in column 2 of the dilution plate to columns 5 through 8 of the cell plate. Transfer 50 μL from each well in column 3 to columns 9 through 12 of the cell plate. Repeat transfer to additional cell plates for remaining samples. Incubate cell plates for 1 hour at 37^{degree}. C. After incubation period, add 50 μL trypsin medium to all wells of the cell plates. Return plates to incubation chamber, and incubate for 4 days post-infection.

Staining

Add 100 μL of crystal violet to all wells. Let sit for 1 hour. Rinse plates in cold running water and air dry.

TABLE-US-00010 TABLE 10 5 mM 10 mM 25 mM 35 mM 50 mM TITRATION OF TCEP USING WHOLE INFLUENZA A VIRUS 30.353 24.58 24.52 24.14 25.9582 30.2261 22.74 24.26 22.74 26.0337 30.28955 23.66 24.39 23.44 25.99595 0.089732 1.301076 0.183848 0.989949 0.053387 Titration of TCEP Using H5 Avian ssRNA 27.2 25.25 25.63 27.3 28.3039 26.73 24.89 25.36 27.62 26.6854 26.965 25.07 25.495 27.46 27.49465 0.33234 0.254558 0.190919 0.226274 1.144452

Time-Course Study of the Long-Term Stability of PrimeStore Compositions

The following data demonstrate the effectiveness of various PrimeStore compositions in preserving nucleic acid integrity over a thirty-day period with samples stored at room temperature. PrimeStore compositions have been compared to water alone, ethanol alone, commercial buffers such as VTM and AVL.

TABLE-US-00011 TABLE 11 30-DAY TIME-COURSE COMPARISON STUDY OF VARIOUS COMPOSITIONS PS-V1 PS-V1 PS-V2.2 VTM Water EtOH AVL (Year old) (new lot) PS-V2 (w/EtOH) DAY 1 27.0225 26.1403 18.4463 24.2698 24.2607 23.9524 23.4426 20.2102 24.42 25.6044 18.3206 24.4789 24.3716 23.9615 23.7387 20.063 25.72125 25.87235 18.4463 24.37435 24.31615 23.95695 23.59065 20.1366 AVG 1.840245 0.378939 0.088883 0.147856 0.0784181 0.006435 0.209374 0.10408612- STDEV DAY 6 29.1988 29.3053 27.4058 37.9226 27.2379 27.165 24.53 22.4887 28.6799 28.7916 27.0781 40 26.4857 26.7658 24.4418 22.4676 28.93935 29.04845 27.24195 38.9613 26.8618 26.9654 24.4859 22.47815 AVG 0.366918 0.363241 0.231719 1.468944 0.5318857 0.282277 0.062367 0.01491995- STDEV DAY 12 27.997 28.151 26.9011 40 30.8352 31.0478 25.8926 22.2074 28.0062 28.2211 26.2139 38.0439 30.4502 30.1935 25.3037 22.0025 28.0016 28.18605 26.5575 39.02195 30.6427 30.62065 25.59815 22.10495 AVG 0.006505 0.049568 0.485924 1.383172 0.2722361 0.604081 0.416415 0.14488618- STDEV DAY 20

27.9851 28.7713 27.1105 40 30.1844 27.193 25.7407 20.8364 28.4067 27.7929 27.0105 40 30.2465 27.2274
 25.6213 20.2843 28.1959 28.2821 27.0605 40 30.21545 27.2102 25.681 20.56035 AVG 0.298116 0.691833
 0.070711 0 0.0439113 0.024324 0.084429 0.39039365 STDEV DAY 30 29.23 31.9168 33.012 40 29.1993
 30.2386 23.0589 20.9348 29.9067 31.3252 32.3001 40 28.827 29.6081 22.9662 20.4973 29.56835 31.621
 32.65605 40 29.01315 29.92335 23.01255 20.71605 AVG 0.478499 0.418324 0.503389 0 0.2632559
 0.445831 0.065549 0.30935922 STDEV PS-N1 (year old) = One-year old PSS (ver. 1). PS-V1 (new lot) =
 Fresh PSS (ver. 2). PS-V2 = Fresh PSS (ver 2) (without ethanol). PS-V2.2 (w/EtOH) = Fresh PSS (ver. 2.2)
 (i.e., with ethanol).

Example 10

Storage Solutions Containing IPCs

As noted herein, in certain embodiments it is desirable to include a nucleic acid carrier molecule and/or an IPC sequence to aid in preparation, stabilization, and quantitation of the isolated polynucleotides.

In one embodiment, the inventors have employed a single-stranded RNA molecule that comprises the sequence of SEQ ID NO:2 which contains the following sequence:

TABLE-US-00012 (SEQ ID NO: 2) 5'-
 CCCUUAGCAGCACGUCAGUCAGGGAGCCAAUUUCAGAGCUCAGCGAG
 ACAGUUUUAUAGGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUC
 CAAAGUUUCAGUU-3'

Data presented in FIG. 15A to FIG. 19 demonstrate the effectiveness of such IPC molecules in stabilizing the isolated population of polynucleotides, and acting as an assayable control sequence. In particular, FIG. 15A and FIG. 15B illustrate the real-time PCR amplification of a single stranded DNA using an assay specific for the detection of the IPC. Likewise, FIG. 16A and FIG. 16B illustrate a similar qRT-PCR amplification of the 130-nt ssRNA.

The data in FIG. 17 demonstrated that a PSS containing a synthetic RNA enhanced the preservation of a naked RNA pathogen target. In that study, a naked influenza A RNA template was added to a human nasal wash specimen, and then stored in PSS that either contained, or lacked the synthetic RNA. The results demonstrated that the presence of the synthetic RNA facilitated enhanced initial extraction of, and long-term preservation of the target sample.

Similarly, the data in FIG. 18 demonstrated that the stability of PrimeStore.TM. solution was enhanced by the addition of a synthetic carrier RNA. In this study, whole influenza A virus was preserved in PSS that either contained, or lacked, the synthetic RNA. The results demonstrated that the synthetic RNA enhanced initial extraction of, and long-term preservation of the target polynucleotide.

Finally, the data in FIG. 19 demonstrated the stability of synthetic RNA at room temperature and human body temperature (37.degree. C.). In this study, the synthetic RNA carrier molecule was added to PSS, and stored for two weeks at either temperature. Stability of the RNA was measured by qRT-PCR. The synthetic carrier RNA was stable at both room temperature and at 37.degree. C. for at least two weeks without significant decline in C.sub.T values.

It is important to note that IPCs useful in the practice of the present invention need not include one of the illustrative sequences described herein, nor do the IPCs even need be substantially homologous to any of the IPC sequences enclosed herein. To illustrate this point, the following sequences represent variants of SEQ ID NO:2 that are also functional as carrier RNA/IPC sequences, despite having sequence degeneracy:

TABLE-US-00013 (SEQ ID NO: 3) 5'-

CUUAGCAGCACGUCAGUCAGGGAGCCAAUUCAGAGCUCAGCGAGACA
 GUUUAUAGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUCCAAG UUUCAGU-3'.
 (SEQ ID NO: 4) 5'-UUAGCAGCACGUCAGUCAGGGAGCCAAUCAGAGCUCAGCGAGACAGU
 UAUAGGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUCCAAGUU UCA-3'. (SEQ ID
 NO: 5) 5'-CCUAGCAGCACGUCAGUCAGGGAGCCAAUUCAGAGCUCAGCGAGACAG
 UUAUAGGCAUGGCAUCAGCUACGCUCGCUCAGGCCUAGUCAGGUCCAAGUU CAG-3'. (SEQ ID
 NO: 6) 5'-AGCAGCACGUCAGUCAGGGAGCCAAUUUCAGAGCUCAGCGAGACAGU
 AUUUAUGGCAUGCAUCAGCUACGCUCGCUCAGCCUAGUCACGUCCAUGUA UCAGU-3'. (SEQ
 ID NO: 7) 5'-CUAGCAGCAGUCAGCAGGGAGCCAAUUUCAGAUUCAGCGAGACAGUUU
 AUAGCCGCAUGCCAUCAGCUACGCUCGCUCAGCCUAGUCAGUCAAGUUUU CAGU-3'. (SEQ ID
 NO: 8) 5'-UAGCAGCACGUCAGUCAGGGAGCCAAUUCAGAGCUCAGCGAGACAGU
 UAGGCAUGCCAUCAACUACGCUCGCUCAGCCUAGUCACGUCCAAGUUCAG AA-3'. (SEQ ID NO:
 9) 5'-CCUAGCAGCACGUCGUCAGGAGCCAAUUCAGAGCUCAGGAGACAGUU
 UAUAGCAUGCAUCAGCUAGCUCGCUCAGCUAGUCAGUCAAGUUCAGU U-3'. (SEQ ID NO: 10)
 5'-CUUACGCAGCACCGGUCAGUAUUCGCGGAGCCU AUUCAGAGCUCAGC
 GAGACAGUU AUAGCAUGCAUCAGCUACCCUCGCUCAGGCCUGUCAGGUCAG UUCGAUU-3'.
 (SEQ ID NO: 11) 5'-UAGCAGCACGUCAGUCAGAGCAUCAGAGCUCAGCGAGACAGAAUUA
 AGCCAUGCAUCAGCUGCUCAGCUAGUCAGUCAAGUCCAGCU-3'. (SEQ ID NO: 12) 5'-
 CCCUAGCAGCACGUCAGUCAGGGAGCCAAUUUCAGAGCUCAGCAGU
 UAUAGGCAUGGCAUCAGCUACGCUCGCAGUAGUCAGGUCCAAGUUCAG U-3'. (SEQ ID NO:
 13) 5'-CAGCACGUCAGUCAGAGCAUCAGAGCUCAGCGAGACAGUAUAGGCAU
 GCAUCAGCACGCUCUCAGGCCUAGUCAGCUCGAAAGUCAGAAU-3'.

Example 11

In Vitro Transcription of ssRNA IPC Sequences

As noted above, the IPCs of the invention may be directly chemically synthesized using conventional methods, or alternatively, prepared using recombinant DNA technology. The present example provides an example of the latter method, and exemplifies preparation of a dsDNA amplicon using PCR and in vitro RNA transcription assays. Using such recombinant techniques, an oligonucleotide is produced for adding to PSS as a stable, IPC that acts as a normalizer for qPCR of pathogen-specific sequences. The 141-nucleotide (nt) sequence is depicted 5' to 3' and shows: 1) T7 RNA transcription start (single underline) site, 2) internal forward and reverse primers (double underline), and 3) real-time internal probe sequence (bold). RNA is synthesis by PCR amplification using T7 forward (single underline) and IPC reverse (double underline) primers and transcribed in vitro using the Green T7 transcription start site. The resultant single stranded RNA polymer is 130-nt in length, and the IPC sequence is located (bound by the forward and reverse primers) within the RNA transcript.

TABLE-US-00014 (SEQ ID NO: 1) 5'-

ATCGTATTAATACGACTCACTATAGGGAATCGTCGTGCAGTCAGTCC

CTCGGTAAAGTCTCGAGTCGCTCTGTCAAAAATATCCGTACCGTAGTCGA
TGCGAGCGAGTCCGATCAGTCCAG GTTTCAAAGTCAAATGACTA-3'.

The specifications for the illustrative DNA sequence employed in this study are as follows:

The total number of DNA nucleotides is 141, while the total number of RNA nucleotides is 130.

TABLE-US-00015 % A = 26.95 (38 nucleotides) % G = 22.70 (32 nucleotides) % T = 26.24 (37 nucleotides)
% C = 24.11 (34 nucleotides) G + C content: % A + T = 53.19 (75 nucleotides) % C + G = 46.81 (66
nucleotides) Base Count: 38 (A), 34 (C), 32 (G), and 37 (T).

The primers and probe sequences for the illustrative IPC assay is located within the synthetic DNA/RNA polymer (shown below). The resulting amplification is 100-bp in length. The primers and probe assay were designed to operate using a two-step thermocycling protocol consisting of one cycle at 50.degree. C. for 20 min and 95.degree. C. for 5 min for reverse transcription and hot start activation, respectively. This is then followed by 40 cycles of 95.degree. C. for 15 sec and 60.degree. C. for 30 sec.

TABLE-US-00016 Forward primer: (SEQ ID NO: 24) 5'-GTGCAGTCAGTCCCTCGGTAA-3'. Reverse primer: (SEQ ID NO: 25) 5'-TTGACTTTGAAACCTGGACTGATC-3'. Probe: (SEQ ID NO: 26) 5'-(FAM)-AAATATCCGTACCGTAGTCG-(MGB)-3'.

TABLE-US-00017 TABLE 12 PANEL OF NUCLEIC ACIDS USED TO EVALUATE CROSS-REACTIVITY OF IPC QPCR ASSAY Viruses Bacteria Other Influenza A (H3) S. pyogenes Total human DNA Influenza A (H1) S. aureus Total human RNA Influenza B B. pertussis Parainfluenza 1, 2, 3, 4 S. pneumoniae CMV (HHV5) N. gonorrhoeae Herpes (HHV1) H. influenzae B Adenovirus 4, 7 E. coli Echovirus RSV Enterovirus

As noted above, it is desirable to formulate an IPC sequence that is both non-genomic, and that does not significantly hybridize to a mammalian genome, or to the genome of pathogenic species of interest.

To that end, the IPC produced in vitro by PCR amplification from SEQ ID NO:1, was screened against a panel of total genomic sequences from a number of viral and bacterial species, as well as total human RNA and DNA. Table 12 lists the species screened; none cross-hybridized to the IPC sequence indicating its desirability as a "nonsense" control sequence that lacked significant homology to any of the tested organisms.

As noted above, the IPCs of the present invention need not be prepared from the precise illustrative DNA amplicon disclosed herein as SEQ ID NO:1. Additional examples of DNA sequences useful in the in vitro preparation of suitable carrier RNA molecules include, without limitation, one or more of the following sequences. In each instance, the polymerase transcription site (here, the T7 transcription site) is shown in single underline, while the sequences of exemplary forward and reverse PCR primer binding domains are shown in double underline. Exemplary sequence domains to which suitable labeled molecular probes are bound are shown in bold.

TABLE-US-00018 (SEQ ID NO: 14) 5'-
X.sub.nTATTAATACGACTCACTATAGGGX.sub.nGTGCAGTCAGTCCCTCGGTT
AAAGTCTCGAGTCGCTCTGTCAAAAATATCCGTACCGTAGTCGATGCGAGC
GAGTCCGATCAGTCCAGGTTTCAAAGTCAAXn-3', wherein X is any nucleotide and n is any integer

from 0 to about 500. (SEQ ID NO: 15) 5'-
 ATCGTATTAATACGACTCACTATAGGGAATCGTCGTGCAGTCAGTCC
 CTCGGTTAAAGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGA
 TGCGAGCGAGTCCGATCAGTCCAGGTTTCAAAGTCAAATGACTA-3'. (SEQ ID NO: 16) 5'-
 ATCGTATTAATACGACTCACTATAGGGAATCGTCGTGCAGTCAGTCC
 CTCGGTTAAAGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGA
 TGCGAGCGAGTCCGATCAGTCCAGGTTTCAAAGTCAAATGACTA-3'. (SEQ ID NO: 17) 5'-
 ATCGTATTAATACGACTCACTATAGGGAATCGTCGTGCAGTCAGTCC
 CTCGGTTAAAGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGA
 TGCGAGCGAGTCCGATCAGTCCAGGTTTCAAAGTCAAATGACTA-3'. (SEQ ID NO: 18) 5'-
 ATATTAATACGACTCACTATAGGGAGTGCAGTCAGTCCCTCGGTTAA
 AGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGA
 GTCCGATCAGTCCAGGTTTCAAAGTCAAAT-3'. (SEQ ID NO: 19) 5'-
 ATATTAATACGACTCACTATAGGGAGTGCAGTCAGTCCCTCGGTTAA
 AGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGA
 GTCCGATCAGTCCAGGTTTCAAAGTCAAAT-3'. (SEQ ID NO: 20) 5'-
 ATATTAATACGACTCACTATAGGGAGTGCAGTCAGTCCCTCGGTTAA
 AGTCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGA
 GTCCGATCAGTCCAGGTTTCAAAGTCAAAT-3'. (SEQ ID NO: 21) 5'-
 TATTAATACGACTCACTATAGGGGTGCAGTCAGTCCCTCGGTTAAAG
 TCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGAGT
 CCGATCAGTCCAGGTTTCAAAGTCAA-3'. (SEQ ID NO: 22) 5'-
 TATTAATACGACTCACTATAGGGGTGCAGTCAGTCCCTCGGTTAAAG
 TCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGAGT
 CCGATCAGTCCAGGTTTCAAAGTCAA-3'. (SEQ ID NO: 23) 5'-
 TATTAATACGACTCACTATAGGGGTGCAGTCAGTCCCTCGGTTAAAG
 TCTCGAGTCGCTCTGTCAAATATCCGTACCGTAGTCGATGCGAGCGAGT
 CCGATCAGTCCAGGTTTCAAAGTCAA-3'.

All of the compositions and methods disclosed and claimed herein can be made and executed without undue experimentation in light of the present disclosure. While the compositions and methods of this invention have been described in terms of exemplary embodiments, it will be apparent to those of ordinary skill in the art that variations may be applied to the composition, methods and in the steps or in the sequence of steps of the method described herein without departing from the concept, spirit and scope of the invention. More specifically, it will be apparent that certain agents that are both chemically- and physiologically-related may be substituted for the agents described herein while the same or similar results would be achieved. All such similar substitutes and modifications apparent to those of ordinary skill in the art are deemed to be within the spirit, scope and concept of the invention as defined by the appended claims. Accordingly, the exclusive rights sought to be patented are as described in the claims below.

SEQUENCE LISTINGS

1

261100DNAArtificial SequenceSynthetic Oligonucleotide 1gtgcagtcag tcctcggtt aaagtctcga gtcgctctgt
 caaaatatcc gtaccgtagt 60cgatgcgagc gaggccgagc agtcagggtt tcaaatgcaa 1002110RNAArtificial
 SequenceSynthetic Oligonucleotide 2cccuuagcag cagcugaguc agggagccaa uuucagagcu cagcgagaca

guuuuauagg 60cauggcauca gcuacgcucg cucagguag ucagguccaa aguuucaguu 1103104RNAArtificial
 SequenceSynthetic Oligonucleotide 3cuuagcagca cgucagucag ggagccaauu cagagcucag cgagacaguu
 uauagcaugg 60caucagcuac gcucgcucag gcuagucagg uccaaguuuu cagu 1044100RNAArtificial
 SequenceSynthetic Oligonucleotide 4uuagcagcac gucagucagg gagccaauca gagcucagcg agacaguuau
 aggcauggca 60ucagcuacgc ucgucagggc uagucagguc caaguuuca 1005100RNAArtificial SequenceSynthetic
 Oligonucleotide 5ccuagcagca cgucagucag ggagccauuc agagcucagc gagacaguuu uaggcauggc 60aucagcuacg
 cucgcucagg cuagucaggu ccaaguucag 1006102RNAArtificial SequenceSynthetic Oligonucleotide
 6agcagcacgu cagucagggg gccaauuuca gagcucagcg agacaguauu uaggcaugca 60ucagcuacgc ucgucagcc
 uagucagcgc cauguauca gu 1027101RNAArtificial SequenceSynthetic Oligonucleotide 7cuagcagcag
 ucagcagggg gccaauuuca gaucagcgcag acaguuuuau gccgcaugcc 60aucagcuacg cucgcuucag cuagucaguc
 aaguuuucag u 101899RNAArtificial SequenceSynthetic Oligonucleotide 8uagcagcagc ucagucaggg
 agccaauca gagcucagcg agacaguuag gcaugccauc 60aacuacgcuc gcucagccua gucacgucca aguucagaa
 99995RNAArtificial SequenceSynthetic Oligonucleotide 9ccuagcagca cgucgucagg agccaauuca gagcucagga
 gagaguuuau agcaugcauc 60agcuagcucg cucagcuagu cagucaaguu caguu 9510104RNAArtificial
 SequenceSynthetic Oligonucleotide 10cuuacgcagc accggucagu auucgcggag ccuauucaga gcucagcgcag
 acaguuuuag 60caugcaucag cuaccucgc ucaggcuguc aggcaguuuc gauu 1041193RNAArtificial
 SequenceSynthetic Oligonucleotide 11uagcagcagc ucagucagag caucagagcu cagcgagaca gaauuaaagc
 ccaugcauca 60gcucgucguc agcuagucag uccaagucca gcu 931297RNAArtificial SequenceSynthetic
 Oligonucleotide 12ccuuagcag cagcucaguc agggagccaa uuucagagcu cagcaguuau aggcauggca 60ucagcuacg
 ucgucaguuu cagguccaaa guucaguu 971389RNAArtificial SequenceSynthetic Oligonucleotide 13cagcagcuca
 gucagagcau cagagcucag cgagacaguu uaggcaugca ucagcagcgc 60cucagguag ucagcucgaa agucagaa
 8914126DNAArtificial SequenceSynthetic Oligonucleotide N = any nucleotide 14ntattaatac gactactat
 agggngtgca gtcagtcct cggtaaagt ctcgagtcgc 60tctgtcaaaa taccgtacc gtagtcgatg cgagcgagtc ctagctcc
 aggttcaaa 120gtcaan 12615141DNAArtificial SequenceSynthetic Oligonucleotide 15atcgtattaa tacgactcac
 tatagggaaat cgtcgtgcag tcagtcctc ggtaaagtc 60tcgagtcgct ctgtcaaaa atccgtacc gtagtcgatg gagcgagtc
 gatcagtc 120ggttcaaag tcaatgact a 14116141DNAArtificial SequenceSynthetic Oligonucleotide
 16atcgtattaa tacgactcac tatagggaaat cgtcgtgcag tcagtcctc ggtaaagtc 60tcgagtcgct ctgtcaaaa atccgtacc
 tagtcgatg gagcgagtc gatcagtc 120ggttcaaag tcaatgact a 14117141DNAArtificial SequenceSynthetic
 Oligonucleotide 17atcgtattaa tacgactcac tatagggaaat cgtcgtgcag tcagtcctc ggtaaagtc 60tcgagtcgct ctgtcaaaa
 atccgtacc gtagtcgatg gagcgagtc gatcagtc 120ggttcaaag tcaatgact a 14118127DNAArtificial
 SequenceSynthetic Oligonucleotide 18atattaatac gactactat agggagtgca gtcagtcct cggtaaagt ctcgagtcgc
 60tctgtcaaaa taccgtacc gtagtcgatg cgagcgagtc ctagctcc aggttcaaa 120gtcaat 12719127DNAArtificial
 SequenceSynthetic Oligonucleotide 19atattaatac gactactat agggagtgca gtcagtcct cggtaaagt ctcgagtcgc
 60tctgtcaaaa taccgtacc gtagtcgatg cgagcgagtc ctagctcc aggttcaaa 120gtcaat 12720127DNAArtificial
 SequenceSynthetic Oligonucleotide 20atattaatac gactactat agggagtgca gtcagtcct cggtaaagt ctcgagtcgc
 60tctgtcaaaa taccgtacc gtagtcgatg cgagcgagtc ctagctcc aggttcaaa 120gtcaat 12721123DNAArtificial
 SequenceSynthetic Oligonucleotide 21tattaatac actactata ggggtgcag cagtcctcg gtaaagtct cgagtcgctc
 60tgcaaaaata tccgtaccg agtcgatgc agcgagtc atcagtcag gttcaaagt 120caa 12322123DNAArtificial
 SequenceSynthetic Oligonucleotide 22tattaatac actactata ggggtgcag cagtcctcg gtaaagtct cgagtcgctc
 60tgcaaaaata tccgtaccg agtcgatgc agcgagtc atcagtcag gttcaaagt 120caa 12323123DNAArtificial
 SequenceSynthetic Oligonucleotide 23tattaatac actactata ggggtgcag cagtcctcg gtaaagtct cgagtcgctc
 60tgcaaaaata tccgtaccg agtcgatgc agcgagtc atcagtcag gttcaaagt 120caa 1232421DNAArtificial
 SequenceSynthetic Oligonucleotide Forward Primer 24gtgcagtcag tcctcggtt a 212524DNAArtificial
 SequenceSynthetic Oligonucleotide Reverse Primer 25ttgacttga aacctggat gatc 242620DNAArtificial
 SequenceSynthetic Oligonucleotide Probe 26aatatcct accgtagtcg 20

* * * * *

